

**PCT**WORLD INTELLECTUAL PROPERTY ORGANIZATION  
International Bureau

## INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

<b>(51) International Patent Classification <sup>6</sup> :</b> <b>C07K 14/00</b>	<b>A2</b>	<b>(11) International Publication Number:</b> <b>WO 99/03882</b> <b>(43) International Publication Date:</b> 28 January 1999 (28.01.99)
<b>(21) International Application Number:</b> PCT/GB98/02058 <b>(22) International Filing Date:</b> 13 July 1998 (13.07.98) <b>(30) Priority Data:</b> 9714760.7 15 July 1997 (15.07.97) GB 9721366.4 9 October 1997 (09.10.97) GB <b>(71) Applicant:</b> ZENECA LIMITED [GB/GB]; 15 Stanhope Gate, London W1Y 6LN (GB). <b>(72) Inventors:</b> AIYAR, Jayashree; 1800 Concord Pike, P.O. Box 15437, Wilmington, DE 19850-8437 (US). LOGSDON, Naomi, Jean; 1800 Concord Pike, P.O. Box 15437, Wilmington, DE 19850-8437 (US). <b>(74) Agent:</b> PHILLIPS, Neil, Godfrey, Alasdair; Zeneca Pharmaceuticals, Intellectual Property Dept., Mereside, Alderley Park, Macclesfield, Cheshire SK10 4TG (GB).		<b>(81) Designated States:</b> AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB, GE, GH, GM, HU, ID, IL, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, UZ, VN, YU, ZW, ARIPO patent (GH, GM, KE, LS, MW, SD, SZ, UG, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG).  <b>Published</b> <i>Without international search report and to be republished upon receipt of that report.</i>
<b>(54) Title:</b> METHODS		
<b>(57) Abstract</b>  A novel human leukocyte calcium activated potassium channel polypeptide is described which is expressed at a high level in activated T-cells. A full length cDNA which encodes the novel calcium activated potassium channel polypeptide is disclosed as well as the interior structural region and the amino acid residue sequence of the native biological molecule. Methods are provided to identify compounds that modulate the biological activity of a human leukocyte calcium-activated potassium channel.		

**FOR THE PURPOSES OF INFORMATION ONLY**

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AL	Albania	ES	Spain	LS	Lesotho	SI	Slovenia
AM	Armenia	FI	Finland	LT	Lithuania	SK	Slovakia
AT	Austria	FR	France	LU	Luxembourg	SN	Senegal
AU	Australia	GA	Gabon	LV	Latvia	SZ	Swaziland
AZ	Azerbaijan	GB	United Kingdom	MC	Monaco	TD	Chad
BA	Bosnia and Herzegovina	GE	Georgia	MD	Republic of Moldova	TG	Togo
BB	Barbados	GH	Ghana	MG	Madagascar	TJ	Tajikistan
BE	Belgium	GN	Guinea	MK	The former Yugoslav Republic of Macedonia	TM	Turkmenistan
BF	Burkina Faso	GR	Greece	ML	Mali	TR	Turkey
BG	Bulgaria	HU	Hungary	MN	Mongolia	TT	Trinidad and Tobago
BJ	Benin	IE	Ireland	MR	Mauritania	UA	Ukraine
BR	Brazil	IL	Israel	MW	Malawi	UG	Uganda
BY	Belarus	IS	Iceland	MX	Mexico	US	United States of America
CA	Canada	IT	Italy	NE	Niger	UZ	Uzbekistan
CF	Central African Republic	JP	Japan	NL	Netherlands	VN	Viet Nam
CG	Congo	KE	Kenya	NO	Norway	YU	Yugoslavia
CH	Switzerland	KG	Kyrgyzstan	NZ	New Zealand	ZW	Zimbabwe
CI	Côte d'Ivoire	KP	Democratic People's Republic of Korea	PL	Poland		
CM	Cameroon	KR	Republic of Korea	PT	Portugal		
CN	China	KZ	Kazakstan	RO	Romania		
CU	Cuba	LC	Saint Lucia	RU	Russian Federation		
CZ	Czech Republic	LJ	Liechtenstein	SD	Sudan		
DE	Germany	LK	Sri Lanka	SE	Sweden		
DK	Denmark	LR	Liberia	SG	Singapore		
EE	Estonia						

## **METHODS**

### **ABSTRACT OF THE DISCLOSURE**

5 A novel human leukocyte calcium activated potassium channel polypeptide is described which is expressed at a high level in activated T-cells. A full length cDNA which encodes the novel calcium activated potassium channel polypeptide is disclosed as well as the interior structural region and the amino acid residue sequence of the native biological molecule. Methods are provided to identify compounds that modulate the biological activity of a human leukocyte  
10 calcium-activated potassium channel.

### **FIELD OF THE INVENTION**

The present invention relates to nucleic acid and amino acid sequences of a novel human lymph node derived calcium activated potassium channel polypeptide and to the use of these  
15 sequences to identify compounds that modulate the activity of human leukocytes. The invention also related to the diagnosis, study, prevention, and treatment of disease related to dysfunctional leukocytes.

### **BACKGROUND OF THE INVENTION**

20 Hematopoietic stem cells give rise to all formed elements of the blood including the leukocytes which are generally implicated in various disease conditions including acute and chronic inflammation, asthma, allergies, graft rejection, proliferative disorders, anemias, septic shock and related disorders, psoriasis, neurodegenerative diseases with immunological components, as well as autoimmune diseases including rheumatoid arthritis, type-1 diabetes  
25 mellitus, multiple sclerosis, myasthenia gravis, systemic lupus erythematosus, Sjogren's syndrome, mixed connective tissue disease, and experimental allergic encephalomyelitis (EAE). Leukocytes comprise neutrophils, eosinophils, basophils, lymphocytes (which include T-cells and B-cells), and monocytes. Leukocytes can traverse every epithelium in the body regardless of its shape, thickness, or permeability. Gallin, E. K., *et al.*, Inflammation: Basic  
30 Principles and Clinical Correlates, Second Edition, 441-458 (1992).

The central lymphoid organs comprise bone marrow and thymus, where stem cells give rise to a diverse progeny of blood cell components including lymphocytes. Phagocytes, neutrophils, and monocytes develop in the bone marrow. After leaving the marrow they enter the bloodstream, where they form a circulating and marginating pool. Monocytes migrate into the  
5 tissues and body cavities, where they differentiate into macrophages. Lymphatic tissues in mammals are classified as either central or peripheral where they generate functional precursor lymphocytes or provide microenvironments for interaction of lymphoid cells with antigen and accessory cells. Lymphocytes are any of the mononuclear, non-phagocytic leukocytes, that are the body's immunologically competent cells and their precursors.

10

Movement of T and B lymphocytes and mononuclear cells through and within lymphatic tissue is vitally important for the intricate inductive and regulatory cellular interactions which take place during the initiation and course of *in vivo* immune responses including inflammation. At the top of the inflammatory response cascade are T-lymphocytes. T-cells  
15 respond to noxious substances by activation and sending out more than a dozen different chemical signals which attract inflammatory cells. Vogel, G., Science, 276:1643 (1997); *Immunophysiology*, Edited by Oppenheim, J.J., et al., Oxford (1990).

The T-cell directs the participation of leukocytes including B-cells and macrophages in overall  
20 physiological humoral and cellular immune responses. T-cells are responsible for directing cells to proliferate, migrate, and differentiate into appropriate effector cells via inter- and intracellular signal transduction. The importance of the bioactive integrity of a normal T-lymphocyte is readily perceived in clinical medicine today a result of the striking demonstration of the consequences of T-cell dysfunction. *Textbook of Internal Medicine*,  
25 Chp. 146, Lippincott, Phil. (1989).

The emigration of inflammatory cells from the blood into the tissues represents one of the most important components of inflammatory response. However, it is not the actual accumulation of cells that is so critical but, rather, what they do in the tissue upon their  
30 arrival. Certain leukocytes (e.g., neutrophils and eosinophils), when appropriately stimulated, move from blood to tissues, and within seconds, release their contents into an endocytic



vacuole or, by fusion with the plasma membrane, to the exterior of the cell. Thus, a major function in inflammatory processes are accomplished. The category of phagocytic inflammatory cells comprises neutrophils, eosinophils, and the mononuclear phagocytic series. Gallin, E. K., *et al.*, *Inflammation: Basic Principles and Clinical Correlates*, Second Edition, 441-458 (1992).

The chronic inflammatory disease, rheumatoid arthritis, for instance, is believed to be mediated by activated T-cells that infiltrate the synovial membrane and initiate a series of inflammatory processes. Panayi, G.S., *et al.*, *The Importance of the T-Cell in Initiating and Maintaining the Chronic Synovitis of Rheumatoid Arthritis*, *Arthritis Rheum*, 35:729 (1992). Accumulating evidence also indicates that the autoimmune disease multiple sclerosis (MS) is mediated by autoreactive T-lymphocytes. Stinissen, P., *et al.*, *Crit. Rev. Immunol.*, 17(1):33 (1997). Autoreactive T-cells have been demonstrated to undergo *in vivo* activation and clonal expansion in patients with MS. Zhang, J., *et al.*, *J. Mol. Med.*, 74(11):653 (1996). In diabetes mellitus, abnormal T-cells systematically destroy pancreatic islet cells such that they prove incapable of producing insulin. Moreover, another propelling recent development in the implication of overactive T-Cells is the recognition that a particular subset of T-lymphocytes appear to be a major culprit in asthma and other allergic diseases, by responding with undue vigor to apparently harmless invaders (rates of asthma per capita in the developing world have increased dramatically in the last several decades; doubling in the U.S. since 1980). *New Clues to Asthma Therapies*: Vogel, G., *Science*, 276:1643 (1997).

Leukocytes receive and respond to a variety of stimuli as part of their role in host defenses and immune function. Physiologic T-cell activation, for example, is achieved when T-cells encounter antigen presenting cells (APCs) bearing their cognate ligand (i.e., a particular peptide bound to a specific class I or class II MHC molecule). *Development and Function of Lymphocytes* (Paul, W. E., *Inflammation: Basic Principles and Clinical Correlates*, Second Edition. Edited by J. I. Gallin, I. M. Godlstein, and R. Snyderman. Raven Press, Ltd., New York (1992)).

Ion channels and fluxes have long been suspected to play a role in lymphocyte signal transduction. By controlling ion fluxes across the plasma membrane, potassium channels mediate changes in intracellular ion concentrations and membrane potential in response to a variety of stimuli. Calcium-activated potassium (K(Ca)) channels cause lymphocytes to  
5 hyperpolarize in response to the elevation of intracellular  $\text{Ca}^{2+}$  triggered by antigen-receptor engagement. The expression of both voltage-gated K(V) and K(Ca) channels is up-regulated as cells progress towards division following mitogenic stimulation. The regulatory aspects of lymphocyte potassium channels including why they are required for T-cell activation, as well as causal relationships among ion channels, membrane potential, mitogen stimulation and  
10 lymphokine gene expression has been reviewed. Lewis, R.S. and Cahalan, M.D., *Potassium and Calcium Channels in Lymphocytes*, Annual Review of Immunol., 13:623 (1995).

Recent advances in cloning of calcium activated potassium channels have greatly expanded the understanding of the structure of these macromolecules. Full-length cDNAs as well as  
15 peptide components of a previously unidentified potassium channel from the brain have been identified. Kohler, M., *et al.*, Science, 273:1709 (1996).

Thus far, the consensus is that at least two types of calcium-activated potassium K(Ca) channels, readily distinguished by their different conductances and pharmacological profiles,  
20 are expressed in a lineage-specific pattern in lymphocytes. The genes for these channels have not yet been cloned. Lewis, R.S. and Cahalan, M.D., *Potassium and Calcium Channels in Lymphocytes*, Annual Review of Immunol., 13:623 (1995).

T- and B-cell mitogens stimulate an increase in both the voltage-gated and calcium activated  
25 potassium channel density, and there is a positive correlation between  $\text{K}^+$ -channel density and proliferative activity in thymocytes. Moreover, pharmacological studies clearly suggest a requirement for functional  $\text{K}^+$  channels in the activation of T and B cells. *Id.*

A number of marketed drugs function as potassium channel antagonists. The most important  
30 of these include the compounds Glyburide, Glipizide and Tolbutamide. These potassium channel antagonists are useful as antidiabetic agents. Potassium channel antagonists are also

utilized as Class III antiarrhythmic agents and to treat acute infarctions in humans. A number of naturally occurring toxins are known to block potassium channels including apamin, iberiotoxin, charybdotoxin (CTX), margatoxin, noxiustoxin, kaliotoxin, dendrotoxin(s), mast cell degranulating (MCD) peptide, and beta -bungarotoxin (beta -BTX).

5

A variety of pharmacological agents have been demonstrated to block the voltage-gated potassium channel and, in parallel, inhibit activation and proliferation of T and B-lymphocytes. Lewis, R.S. and Cahalan, M.D., *Potassium and Calcium Channels in Lymphocytes*, Annual Review of Immunol., 13:623 (1995).

10

Unfortunately, no compound has yet been found to block *calcium activated potassium channels* exclusively, especially in human leukocytes. Accordingly, there remains a need to identify a calcium activated potassium channel of leukocyte origin as a pharmacological target for screening candidate compounds for the modulation of cell activity. Such modulators are potentially useful in treating disorders manifested by dysfunctional leukocytes. There continues to be a need for an ability to screen for such compounds. Selective channel blockers are moreover needed to address the contribution of calcium activated potassium channels to the membrane potential and to signaling during T-cell activation. Grissmer, S., *et al.*, *Calcium-activated Potassium Channels in Resting and Activated Human T Lymphocytes*, J. Gen. Physiol, 102:601 (1993). The central role of potassium channels in regulating numerous cell functions makes them particularly important targets for therapeutic development. Specific immune therapies designed to control T-cell activation and proliferation are likely to improve the clinical course of various autoimmune diseases.

## 25 SUMMARY OF THE INVENTION

The present invention is directed to an isolated and purified polynucleotide molecule, which encodes a polypeptide of a potassium channel, or a biologically active derivative thereof comprising a nucleic acid sequence encoding the polypeptide having the sequence substantially as depicted in SEQ ID NO:3 or a biologically active fragment thereof. Isolated and purified polynucleotides for the present invention include but are not limited to SEQ ID

30

NO:1, SEQ ID NO:7 (calcium-activated potassium channel cDNA) and SEQ ID NO:2 (calcium-activated potassium channel structural coding region).

In addition, the current invention is directed to a purified polypeptide comprising the amino acid sequence substantially as depicted in SEQ ID NO:3 which functions as a polypeptide subunit of a human calcium-activated potassium channel.

The invention is further directed to an expression vector for expression of a polypeptide of a calcium-activated potassium channel in a recombinant host, wherein said vector contains a polynucleotide comprising a nucleic acid sequence encoding the polypeptide of a calcium-activated potassium channel having the sequence substantially as depicted in SEQ ID NO:3 or a biologically active derivative thereof.

Further the invention is directed to a host cell containing an expression vector for expression of a polypeptide of a calcium-activated potassium channel, wherein said vector contains a polynucleotide comprising a nucleic acid sequence encoding the polypeptide of a calcium-activated potassium channel having the sequence substantially as depicted in SEQ ID NO:3 or a biologically active derivative thereof.

The instant invention is further directed to a method of identifying compounds that modulate a potassium channel activity, comprising:

- (a) combining a candidate compound modulator of the activity of a potassium channel with the potassium channel having the sequence substantially as depicted in SEQ ID NO:3, and
- (b) measuring an effect of the modulator on the channel.

The invention is also directed to a method of identifying compounds that modulate the activity of a potassium channel, comprising:

(a) combining a candidate compound modulator of a potassium channel activity with a host-cell expressing the polypeptide of a potassium channel having the sequence substantially as depicted in SEQ ID NO:3, and

5 (b) measuring an effect of the modulator on the channel.

The present invention is also directed to active compound identified by means of the aforementioned methods, wherein said compounds modulate the biological activity of a potassium channel.

10

Further, the invention is directed to a pharmaceutical composition comprising a compound active in the aforementioned methods, wherein said compound is a modulator of a potassium channel.

15 Additionally, the invention is directed to a novel treatment of a patient in need of such treatment for a condition which is mediated by a potassium channel, comprising administration of a potassium channel modulating compound active in the aforementioned method.

20 The invention is further directed to an antisense polynucleotide molecule comprising substantially the complement of SEQ ID NO:2 or a biologically-effective portion thereof, as well as a method for inhibiting the expression of a potassium channel in a cell comprising administering an effective amount of the antisense molecule.

25 The current invention is also drawn toward an antibody specific for a purified polypeptide comprising the amino acid sequence substantially as depicted in SEQ ID NO:3, and a diagnostic composition for the identification of a polypeptide sequence comprising the amino acid sequence substantially as depicted in SEQ ID NO:3.

30

**BRIEF DESCRIPTION OF THE FIGURES**

Figure 1 displays SEQ ID NO:1 which is the 2262 base full-length cDNA nucleic acid sequence which encodes the novel human lymph node derived calcium activated potassium channel polypeptide (hKCa4) described herein.

5

Figure 2 displays SEQ ID NO:2 which is the translated structural region, **ATG** to **TAG**, of the cDNA nucleic acid sequence which encodes the novel human lymph node derived calcium activated potassium channel polypeptide described herein.

10 Figure 3 displays SEQ ID NO:3 which is the 427 amino acid residue sequence of the novel human lymph node derived calcium activated potassium channel polypeptide described herein.

Figure 4 shows SEQ ID NO:4 which is the amino acid residue sequence of the HSK1 brain  
15 derived small-conductance calcium activated potassium channel polypeptide (Kohler, M., *et al.*, Science, 273:1709 (1996)).

Figure 5 shows SEQ ID NO:5 which is the amino acid residue sequence of the RSK2 brain  
20 derived small-conductance calcium activated potassium channel polypeptide (Kohler, M., *et al.*, Science, 273:1709 (1996)).

Figure 6 shows SEQ ID NO:6 which is the amino acid residue sequence of the RSK3 brain  
derived small-conductance calcium activated potassium channel polypeptide (Kohler, M., *et al.*, Science, 273:1709 (1996)).

25

Figure 7 shows a comparison between the amino acid residue sequence of the novel human lymph node derived calcium activated potassium channel polypeptide described herein (SEQ ID NO:3) (designated HKCA4), and the amino acid residue sequences of the brain derived small-conductance calcium activated potassium channel polypeptides HSK1 (SEQ ID NO:4),  
30 RSK2 (SEQ ID NO:5), and RSK3 (SEQ ID NO:6). Conserved amino acid residues are boxed. Dashes represent gaps introduced to optimize the alignment. Sequences shown in this figure

were produced using the multisequence alignment program of DNASTAR software (DNASTAR Inc, Madison WI).

Figure 8 displays a hydrophobicity plot of the 427 amino acid residue sequence of the novel human lymph node derived calcium activated potassium channel polypeptide (SEQ ID NO:3). Six (6) transmembrane regions are indicated, S1 through S6. Kyte, J., Doolittle, R.F., J. Mol. Biol., 157:105 (1982). (DNASTAR Inc, Madison WI).

Figure 9 shows a strong mRNA signal of the the novel human lymph node derived calcium activated potassium channel described herein upon T-cell activation in Northern blot analysis.

Figure 10 shows SEQ ID NO:8 which is a partial cDNA clone of the novel human lymph node derived calcium activated potassium channel polypeptide described herein.

Figure 11 shows SEQ ID NO:9 which is a partial cDNA clone of the novel human lymph node derived calcium activated potassium channel polypeptide described herein.

Figure 12 shows SEQ ID NO:7 a 2238 base full-length cDNA nucleic acid sequence which comprises SEQ ID NO:2 and encodes the novel human lymph node derived calcium activated potassium channel polypeptide described herein (a 5' UTR variation of SEQ ID NO:1, the initiator codon and stop codon are underlined).

Figure 13 demonstrates KCa currents (biological activity) in HEK 293 cells transfected with SEQ ID NO:2 (a) Dependence of whole cell current on 1  $\mu$ M ionomycin perfusion and block by CTX (100 nM). Superimposed currents were evoked by 200 ms voltage ramps from -100 to +40 mV ( $E_{\text{hold}} = -50$  mV). The pipette solution contained 100 nM  $[\text{Ca}^{2+}]_{\text{free}}$ , and the bath solution contained 160 mM  $\text{K}^+$ . (b) Concentration-dependent block of KCa current by CTX ( $\bullet$ ), TEA ( $\blacksquare$ ), and Clotrimazole ( $\blacktriangle$ ). Protocol is the same as in (a), except that  $[\text{Ca}^{2+}]_{\text{free}}$  was buffered to 1  $\mu$ M in the pipette and ionomycin was omitted from the bath. Currents were measured at -95 mV on the ramp and normalized between control amplitude and that obtained during perfusion with 100 nM CTX in the same experiment. Each point is the mean  $\pm$  SEM

of 3-5 experiments. Solid lines are fits to a Hill equation of the following form:  $100 / [1 + (K_d/x)^n]$ , where  $x$  is concentration,  $K_d$  the concentration producing 50 percent inhibition, and  $n$  is the slope factor of the line. See text for fitted parameters (c)  $K^+$  selectivity of KCa current. Each point represents the mean  $\pm$  SEM voltage ( $n = 3-6$ ) where control current  
 5 converged with current in CTX during voltage ramps at the designated  $K^+_{(out)}$  concentration. Solid line is a linear regression to the data (slope 57 mV).

Figure 14 displays single channel records and calcium-sensitivity of the novel channel (SEQ ID NO:2): (a) Single-channel current from a cell-attached patch with 160  $K^+$  in the pipette.  
 10 The bath contained Ringer with 1  $\mu M$  ionomycin. The patch was ramped over 100 ms from -120 to +60 mV. A blank trace (i.e., no channel openings) was digitally subtracted. The dashed line represents a slope conductance of 31 pS. (b) Representative traces obtained from an inside-out patch excised from an hKCa4-transfected cell ( $E_{hold} = -80$  mV). The pipette contained 160 mM  $K^+$  and patch was perfused with the  $[Ca^{2+}]_{free}$  indicated. (c) Relationship  
 15 between open probability and cytoplasmic  $[Ca^{2+}]_{free}$  determined using the experimental protocol in (b). Single channel open probabilities were determined by the relative areas of Gaussian fits to amplitude-frequency histograms. Points are the mean  $\pm$  SEM from designated number of experiments. Solid line is a fit to the Hill equation described in FIG.13b. (d) Summary table of data obtained from experiments shown in FIG.13 & 14 for SEQ ID NO:2-  
 20 expressed in HEK cells compared to published data for the human T cell KCa channel.

Figure 15 displays the nucleic acid sequence of SEQ ID NO:7 as well as the encoded amino acid sequence (SEQ ID NO:3). Predicted transmembrane regions are boxed and shaded. Consensus sites for glycosylation are indicated by a triangle, serine-threonine phosphorylation  
 25 sites are circled, the polyadenylation signal is underlined, and the stop codon is indicated by an asterisk. (b) Hydropathy plot of SEQ ID NO:3 showing six predicted transmembrane segments and a pore between S5 and S6. (c) Amino acid dendrogram of SEQ ID NO:3 with the recently cloned small-conductance KCa channels, hSK1, rSK2 and rSK3 (SEQ ID NO:4, SEQ ID NO:5 and SEQ ID NO:6 respectively).



Figure 16 shows Northern blot analysis of SEQ ID NO:7 message: (a) in multiple tissues, (b) in resting and activated human T lymphocytes, showing up-regulation of message, (c) patch clamp analysis of representative T cells that were subject to Northernblots in (b) showing up-regulation of K<sub>Ca</sub> channel currents upon activation. Each superimposed current trace was elicited by a 200 ms voltage ramp from -100 to +40 mV ( $E_{\text{hold}} = -50$  mV). Trace 1: Resting T cell, 2: Resting T cell plus 100 nM CTX, 3: Activated T cell, 4: Activated T cell plus 100 nM CTX. Components of K<sub>v</sub> and K<sub>Ca</sub> currents are shown.

Figure 17 shows patch clamp current recording of stable HEK293/SEQ ID NO:2 transfectants.

Figure 18 shows a dose-response curve for CTX-binding hKCa4 (SEQ ID NO:3) channels.

Figure 19 shows Western Blot analysis of anti-SEQ ID NO:3 antibodies.

Figure 20 demonstrates down-regulation of the calcium activated potassium channel (SEQ ID NO:3) currents in T cell-lines induced to anergy.

Figure 21 illustrates deletion analysis of the SEQ ID NO:3 COOH-tail to identify the calmodulin interaction site.

Figure 22 shows the full-length cDNA (SEQ ID NO:13) of the murine Kca4 ortholog which includes the structural coding region (underlined ATG → TAG).

Figure 23 shows the native amino acid sequence (SEQ ID NO:14) of the murine Kca4 homolog.

Figure 24 shows amino acid alignment of mouse and human KCa4 homolog sequences.

**DETAILED DESCRIPTION OF THE INVENTION**

Unless defined otherwise, all technical and scientific terms used herein have the same meaning as is commonly understood by one of skill in the art to which this invention belongs. All publications and patents referred to herein are incorporated by reference.

5

Biological activity as used herein refers to the ability of a calcium activated potassium channel to allow transmembrane potassium ion flow and/or transport or regulate transmembrane potassium ion flow and/or transport. Biological activity as used herein is also intended to encompass pharmacological activity, *per se*, as in the ability of a subunit to bind or associate  
10 with at least one other channel peptide subunit, ligand, or cofactor (e.g., calmodulin and/or calcium) and/or otherwise modulate the biological activity of a potassium channel.

Biologically active fragment as used herein includes peptides which have been truncated with respect to either the N- or C-termini, or both; or the corresponding 5' or 3' end of the  
15 corresponding nucleic acid coding regions, or both, which fragments perform substantially the same function or encode peptides which perform substantially the same function in substantially the same way as the precursor. Biologically active fragment as used herein is also intended to encompass pharmacologically active dominant negative mutants contemplated herein.

20

Nucleic acid sequence as used herein refers to an oligonucleotide, nucleotide or polynucleotide sequence, and fragments or portions thereof, and to DNA or RNA of genomic or synthetic origin which may be double-stranded or single-stranded whether representing the sense or antisense strand. Similarly, amino acid sequence as used herein refers to peptide or  
25 protein sequences or portions thereof. The term 'purified' refers to molecules, either nucleic or amino acid sequences, that are removed from their natural environment and isolated or separated from at least one other component with which they are naturally associated.

'Substantially as depicted' as used herein refers to functional derivative proteins, variant  
30 peptides and DNA sequences that may have changes but perform substantially the same

biochemical function in substantially the same way; however, 'substantially as depicted' as used herein also refers to dominant negative mutant versions.

Direct administration as used herein refers to the direct administration of nucleic acid  
5 constructs which encode embodiments (e.g., SEQ ID NO:3, dominant negative mutant, modulator compound molecule, antisense molecule, antibody molecule) of the present invention or fragments thereof; and the direct administration of embodiments of the present invention or fragments thereof, and the *in vivo* introduction of molecules of the present invention preferably via an effective eukaryotic expression vector in a suitable pharmaceutical  
10 carrier. Polynucleotides and therapeutic molecules of the present invention may also be delivered in the form of nucleic acid transcripts.

The term 'modulation' is used herein to refer to the capacity to either enhance or inhibit a functional property of a subunit or potassium channel. The term 'modulation' is also used  
15 herein to refer to the capacity to affect the biophysical activity of a cell. Modulate physiology as used herein refers to the biophysiological regulation of cells and/or tissue and the treatment of pathophysiological disorders related thereto.

Modulation or regulation of biological activity and/or pharmacological activity as used herein  
20 refers to binding, blocking, antagonization, repression, neutralization, or sequestration, of a potassium channel biomolecular structure including but not limited to the novel calcium activated potassium channel polypeptide described herein; as well as *up regulation* or agonization or activation of a potassium channel by a compound identified by means described herein.

25

The nucleic acid sequence also provides for the design of antisense molecules useful in downregulating, diminishing, or eliminating expression of the genomic nucleic acid sequence in cells including, but not limited to, leukocytes, endothelial cells, microglia, and tumor or cancer cells.

30

Expression vector as used herein refers to nucleic acid vector constructions which have components to direct the expression of heterologous protein coding regions including coding regions of the present invention through accurate transcription and translation in host cells.

Expression vectors usually contain a promoter to direct polymerases to transcribe the  
5 heterologous coding region, a cloning site at which to introduce the heterologous coding region, and usually polyadenylation signals. Expression vectors include but are not limited to plasmids, retroviral vectors, viral and synthetic vectors.

Transformed host cells as used herein refer to cells which have coding regions of the present  
10 invention stably integrated into their genome, or episomally present as replicating or nonreplicating entities in the form of linear nucleic acid or transcript or circular plasmid or vector.

15

#### **Calcium activated potassium channels**

Calcium-activated potassium (K(Ca)) channels cause lymphocytes to hyperpolarize in response to the elevation of intracellular  $\text{Ca}^{2+}$  triggered by antigen-receptor engagement. The expression of both voltage-gated K(V) and calcium-activated potassium K(Ca) channels is up-  
20 regulated as cells progress towards division following mitogenic stimulation. Calcium signal is required to enable a number of essential activation events closely linked to the opening and closing of calcium and potassium channels as demonstrated by  $\text{Ca}^{2+}$  imaging in single cells. Sustained signaling and oscillations depend absolutely on plasma-membrane  $\text{Ca}^{2+}$  channels that are activated indirectly by the depletion of intracellular calcium stores. Lewis, R.S. and  
25 Cahalan, M.D., *Potassium and Calcium Channels in Lymphocytes*, Annual Review of Immunol., 13:623 (1995).

The surface density of potassium channels in the plasma membrane of T-lymphocytes is regulated during development and during activation of mature cells by mitogens. For reviews,  
30 see Cahalan, M.D., *et al.*, Curr. Top. Membr., 39:357-394; Lewis, R.S., *et al.*, Trends Neurosci., 11:214-218. Expression of calcium activated potassium K(Ca) channels in human

T-cells is profoundly up-regulated by mitogens, increasing from ~20 per resting cell to >500 channels per T-cell blasts treated with phytohemmagglutinin (PHA). Thus, the relative contribution of K(Ca) channels to the membrane potential appears to be greatly enhanced in activated cells and perhaps in memory cells. Lewis, R.S. and Cahalan, M.D., *Potassium and Calcium Channels in Lymphocytes*, Annual Review of Immunol., 13:623 (1995).

Grissmer, S., *et al.*, have further demonstrated that mitogenic stimulation of human T-cells results in an increase in the cytosolic free  $\text{Ca}^{2+}$  concentration, that is vital for subsequent events leading to cell proliferation and secretion of lymphokines. Calcium activated potassium channels (K(Ca)) were characterized in resting and activated human peripheral blood T-lymphocytes using simultaneous patch-clamp recording and fura-2 monitoring of cytosolic  $\text{Ca}^{2+}$  concentration. Whole-cell experiments, using EGTA-buffered pipette solutions to raise  $[\text{Ca}^{2+}]_i$  to 1  $\mu\text{M}$ , revealed a 25-fold increase in the number of conducting calcium-activated potassium channels per cell, from an average of 20 in resting T cells to >500 channels per cell in T cell blasts after mitogenic activation. *Calcium-activated Potassium Channels in Resting and Activated Human T Lymphocytes*, J. Gen. Physiol, 102:601 (1993) which is incorporated herein by reference. The high flow rate of calcium activated potassium channels, and discrete opening and closing, permits measurement of currents through single channels using standard patch-clamp techniques. Reported results indicate that expression of the charybdotoxin (CTX) -sensitive calcium activated potassium channel increases dramatically after activation. Activated T-cell blasts express many more calcium activated potassium channels than resting HPB T-lymphocytes. The increased number of K(Ca) channels reflects an increase primarily in the surface density of channels. J. Gen. Physiol, 102:601 (1993).

Charybdotoxin, (CTX), a 37 amino-acid peptide purified from the venom of scorpion *Leiurus quinquestratus* (Gimenez-Gallego et al., PNAS, 85: 3329-3333, 1988) has been shown to block mitogen-induced proliferation of human peripheral T-lymphocytes (Lin C et al., FASEB J. 6: 1693, 1992). This effect of CTX has been demonstrated to act via blockade of potassium channels which maintain the cell's resting membrane potential (Leonard R et al., PNAS 89: 10094-98, 1992). CTX is a blocker of both the calcium-activated potassium channel (KCa)

- and the voltage-gated potassium channel (Kv1.3) in T-cells (Grissmer et al., J. Gen. Physiol. 102: 601-630, 1993). CTX affects only the calcium-dependent activation pathways in lymphocytes (Garcia ML, J Bioenerg. Biomembr 23: 615-646, 1991). Radiolabeled CTX has been used as a ligand for high-throughput screening (HTS) of drugs for Kv1.3 (Hill RJ, Mol. Pharm. 48: 98-104, 1995; Deutsch C et al., J. Biol. Chem. 266: 3668-3674, 1991). Although other peptide toxins that selectively block Kv1.3 have been identified, there is no pharmacological agent that selectively blocks the KCa channel in lymphocytes with high potency. With the cloning of the novel calcium activated potassium channel described herein, it is now possible to generate a cell-line that over-expresses the channel which can then be used in high throughput screening (HTS) to screen for selective blockers of the calcium-activated potassium channel in T-lymphocytes. These channel-blocking agents could also be used as research tools to understand the physiological role of these channels in a variety of cells of hematopoietic and other origin.
- 15 Selective channel blockers are needed to address the contribution of calcium activated potassium channels to the membrane potential and to signaling during T-cell activation. Grissmer, S., et al., *Calcium-activated Potassium Channels in Resting and Activated Human T Lymphocytes*, J. Gen. Physiol, 102:601 (1993).
- 20 The calcium activated potassium channels described by Grissmer, S., et al., were moreover reported to resemble those found in red blood cells, the cell type in which K(Ca) channels were first discovered through tracer flux measurements. Gardos, G., *The Function of Calcium in the Potassium Permeability of Human Erythrocytes*, Biochim. Biophys. Acta., 30:653 (1958); reviewed by Schwarz and Passow, *Calcium Activated Potassium Channels in Erythrocytes*, Annual Review of Physiology, 45:359 (1983); see also, Brugnara, C., et al., *Therapy with Oral Clotrimazole Induces Inhibition of the Gardos Channel and Reduction or Erythrocyte Dehydration in Patients with Sickel Cell Disease*, J. Clin. Investigation, 97(5):1227 (1996).
- 30 U.S. Patent No. 5,397,702, *Assay For and Treatment of Autoimmune Diseases*, issued March 14, 1995; U.S. Patent No. 5,637,470, *Screening array using cells expressing recombinant*

*alpha and beta subunits of the mammalian large-conductance (maxi-K) potassium channel*, issued June 10, 1997; U.S. Patent No. 5,607,843, *Nucleic Acid Sequence Encoding Apamin Binding Protein*, issued March 4, 1997; and U.S. Patent No. 5,602,169, *3-substituted oxindole derivatives as potassium channel modulators*, issued Feb. 11, 1997 are, each of which, as  
5 well as methods described, incorporated herein by reference.

### **IK (intermediate conductance) calcium activated potassium channels**

The IK channel is expressed in hematopoietic and endothelial cells, certain cells of epithelial origin and in red blood cells. In the red blood cells, where the channel has been denominated  
10 the Gardos channel, a rise in the concentration of intracellular  $\text{Ca}^{2+}$  opens the channel and causes potassium loss and cell dehydration, a condition which is exacerbated in sickle cell anemia. Promising therapeutic approaches for sickle cell anemia involve specifically blocking the IK channel. IK channels have also been implicated in the microvasculature of the kidney, where they may be responsible for the vasodilatory effects of bradykinin. The decrease in  
15 blood pressure during septic shock is caused by an increased NO production by the endothelial cells, and the IK channels in these cells are responsible for maintaining the  $\text{Ca}^{2+}$  influx activating the  $\text{Ca}^{2+}$ -sensitive NO-synthase. In brain capillary endothelial cells, IK channels, activated by endothelin that is produced by neurons and glia, shunt excess  $\text{K}^{+}$  into the blood. Neutrophil granulocytes, mobile phagocytic cells that defend against microbial  
20 invaders, undergo large depolarization subsequent to agonist stimulation, and IK channels have been implicated in repolarizing the stimulated granulocyte. The IK channel is sensitive to a number of peptide toxins as well as organic blockers. Charybdotoxin is the most potent peptide blocker ( $\text{IC}_{50} = 5\text{-}10\text{ nM}$  on cloned hKCa4 (SEQ ID NO:3, *infra*)) and clotrimazole is the most potent of the antimycotics ( $\text{IC}_{50} = 100\text{-}300\text{ nM}$ ) (Figs 5-6). Clotrimazole is the  
25 compound which is being evaluated for the treatment of sickle cell anemia.

### **Brain-derived calcium activated potassium channels**

Köhler, M. *et al.*, describe the molecular structures of members of a previously unidentified family of calcium activated potassium channels with six transmembrane domains cloned from  
30 rat and human brain. Three full-length coding sequences were isolated and translated into their respective amino acid residue sequences, one from human, hSK1 (SEQ ID NO:4) (561

residues), and two from rat brain, rSK2 (SEQ ID NO:5) (580 residues), and rSK3 (SEQ ID NO:6) (553 residues). Expression of the respective mRNAs in *Xenopus* oocytes resulted in calcium-activated, voltage-independent potassium channels. Science, 273:1709 (1996).

Hydrophobicity analysis predicts six transmembrane segments of each with the NH<sub>2</sub>- and  
5 COOH-termini residing inside the cell. The sequences are highly conserved across their transmembrane cores (80 to 90% homology) but diverge in sequence and length within their NH<sub>2</sub>- and COOH-terminal domains. The cloned channels show no significant amino acid homology to other cloned potassium channel subunits except for a 12-residue stretch within the putative pore domain. Hydrophobicity plots predicted that these subunits contain six  
10 transmembrane domains, a topology shared with members of the voltage-gated class of K<sup>+</sup> channels. Jan, L.Y., et al., Nature, 345:672 (1990). Although related in topology to voltage-dependent K<sup>+</sup> channels, including a P region and S4 segment, the Köhler, M. *et al.* clones reside on a distinct evolutionary branch within the K<sup>+</sup> channel superfamily. Science, 273:1709 (1996).

15

**Novel human lymph node derived calcium activated potassium channel (IK) polypeptide**

Potassium channels play a critical role in modulating calcium signaling of lymphocytes.

Verheugen, J. A., *et al.*, Cell Calcium, 17: 287 (1995). Human T lymphocytes express at least two types of potassium channels: those that open in response to changes in membrane  
20 potential (Kv channels) and those that are activated following elevations of intracellular calcium levels (KCa channels). Lewis, R. S., *et al.*, Annu. Rev. Immunol., 13:623 (1995). The predominant Kv channel in human T cells is encoded by Kv1.3, a *Shaker*-related voltage-gated potassium channel gene. Kv1.3 has been extensively characterized at the molecular and physiological level and plays a vital role in controlling T cell proliferation, mainly by  
25 maintaining the resting membrane potential of the cell. Chandy, K. G., *et al.*, Sem. Neurosci., 5:125 (1993). Upon T cell activation, there is at most a 2-fold increase in Kv1.3 currents. The predominant KCa channel in human T lymphocytes is of intermediate conductance, is voltage-insensitive and is potently blocked by the scorpion venom peptide, charybdotoxin (CTX; 4). Unlike Kv1.3, KCa channel currents are up-regulated dramatically (10-25 fold)  
30 following mitogenic or antigenic stimulation, and are thought to play a significant role in post-activation and secondary immune phenomena. Grissmer, S., *et al.*, 102:601 (1993);



Verheugen, J. A., *et al.*, Cell Calcium. 21:1 (1997). KCa channels with very similar biophysical and pharmacological properties to the T lymphocyte channel have also been identified in red blood cells, e.g., Gardos' channel, macrophages, neutrophils, B lymphocytes as well as in other non-immune peripheral tissues. Grygorczyk, R., *et al.*, Biophys. J., 45:693  
5 (1984); Gallin, E. K., Am. J. Physiol., (Cell Physiol. 26) 257:C77-C-85 (1989); Karl-Heinz, *et al.*, J. Clin. Invest., 85:491 (1990); Mahaut-Smith, M. P., *et al.*, J. Physiol., 415:69 (1989). However, the molecular identity of this channel type has hitherto been unknown. The cloning and characterization of an intermediate-conductance, CTX-sensitive KCa channel, called hKCa4 (SEQ ID NO:2; SEQ ID NO:3) from a human lymph node cDNA library is described  
10 herein. Convergent molecular, biophysical, and pharmacological evidence that the novel sequences presented herein (e.g., SEQ ID NO:2, SEQ ID NO:3, SEQ ID NO:1, and SEQ ID NO:7) correspond to the the predominant KCa channel in human T cells is presented.

Example full-length cDNAs from human lymph node which comprises the coding region for  
15 the novel calcium activated potassium channel polypeptide described herein is set forth in FIG.1 (SEQ ID NO:1) and FIG.12 (SEQ ID NO:7). The native structural coding region for the subject polypeptide is set forth in FIG.2 (SEQ ID NO:2). The amino acid residue sequence of the native calcium activated potassium channel polypeptide described herein is set forth in FIG.3 (SEQ ID NO:3). The translated protein SEQ ID NO:3 is comprised of 427  
20 amino acid residues FIG.3 and appears to have 6 transmembrane regions, i.e. residues 25-42, 64-79, 105-120, 150-174, 205-223, and 265-285, S1 through S6, respectively, revealed by hydropathy analysis FIG.8. Moreover the SEQ ID NO:3 polypeptide displays a pore region of the calcium activated potassium channel, i.e., residues 245-260 which contain the consensus signature sequence of potassium channels. Heigenbothams, L., *et al.*, Biophys. J., 66(4):1061  
25 (1994).

The amino acid sequence of the novel calcium activated potassium channel (SEQ ID NO:3) is about 41% homologous to the calcium activated potassium channels (hSK1, rSK2 and rSK3) cloned from brain. Kohler, M., *et al.*, Science, 273:1709 (1996). FIG.4 shows SEQ ID NO:4  
30 which is the amino acid residue sequence of the HSK1 brain derived small-conductance calcium activated potassium channel polypeptide. FIG.5 shows SEQ ID NO:5 which is the

amino acid residue sequence of the RSK2 brain derived small-conductance calcium activated potassium channel polypeptide. FIG.6 shows SEQ ID NO:6 which is the amino acid residue sequence of the RSK3 brain derived small-conductance calcium activated potassium channel polypeptide. FIG. 7 shows a comparison between the amino acid residue sequence of the  
5 novel human lymph node derived calcium activated potassium channel polypeptide described herein (SEQ ID NO:3) (designated HKCA4), and the amino acid residue sequences of the brain derived small-conductance calcium activated potassium channel polypeptides HSK1 (SEQ ID NO:4), RSK2 (SEQ ID NO:5), and RSK3 (SEQ ID NO:6). Conserved amino acid residues are boxed. Dashes represent gaps introduced to optimize the alignment. Sequences  
10 shown in this figure were produced using the multisequence alignment program of DNASTAR software (DNASTAR Inc, Madison WI).

The pore region, between residues 245-260 of SEQ ID NO:3 show significant amino acid differences from its homologs, HSK1 (SEQ ID NO:4), RSK2 (SEQ ID NO:5) or RSK3 (SEQ  
15 ID NO:6), see FIG.7, predicting different pharmacology and biophysics from the brain derived calcium activated potassium channels. Kohler, M., *et al.*, Science, 273:1709 (1996). A splice variant of the subject gene has also been identified as discussed *infra* at the region corresponding to the S2-S3 regions of SEQ ID NO:3.

20 Similar to the previous description by Grissmer, S., *et al.*, wherein expression of the charybdotoxin (CTX) -sensitive calcium activated potassium channel increases dramatically after T-cell activation. J. Gen. Physiol, 102:601 (1993). The novel human lymph node derived calcium activated potassium channel polypeptide described herein (SEQ ID NO:3) is shown have high level expression, demonstrated via Northern analysis FIG.9, in activated T-  
25 cells relative to resting T-cells. FIG.9 shows a strong mRNA signal of the the novel human lymph node derived calcium activated potassium channel described herein upon T-cell activation in Northern blot analysis.

The observation of dramatic up-regulation of the novel calcium activated potassium channel  
30 upon T-cell activation provides evidence that this channel plays a crucial role in lymphocyte signalling events. In addition, high expression of the calcium activated potassium channel in

prostate, colon and placenta suggests an important role in the physiology of these tissues as well. Transcripts of the novel calcium activated potassium channel described herein are not detected in the brain. The expression of the subject channel in cells of the myeloid, lymphoid and erythroid lineages indicates that they could play a major role in modulating immune responses as well as in red blood cell and platelet physiology. Accordingly, due to the evidence of very strong mRNA signal in activated T-cells and the fact that SEQ ID NO:1 (SEQ ID NO:7) was isolated as described *infra* from a lymph node library which is enriched for activated T-cells the polypeptide of the current invention appears to be an ideal pharmacological target for modulating the activity of human leukocytes and other cells of hematopoietic origin and T-cells in particular.

It is therefore probable that compounds (e.g., small-molecules, peptides, analogs, mimetics) that modulate the lymph node derived calcium activated potassium channel described herein could be used in the treatment of a variety of disease conditions including acute and chronic inflammation, asthma, allergies, graft rejection, proliferative disorders, anemias, neurodegenerative diseases with immunological components, as well as autoimmune diseases including rheumatoid arthritis, type-1 diabetes mellitus, multiple sclerosis, myasthenia gravis, systemic lupus erythematosus, Sjogren's syndrome, mixed connective tissue disease, and experimental allergic encephalomyelitis (EAE).

20

A gene sequence described herein (e.g., SEQ ID NO:7) which encodes an intermediate-conductance, calcium-activated potassium channel, has been isolated from a human lymph node library. The translated protein is 427 amino-acids long, has 6 transmembrane segments, S1-S6, and a pore-motif between S5 and S6. The novel protein (SEQ ID NO:3) shares 40-42% similarity at the amino acid level with three small-conductance calcium-activated potassium channels cloned from brain. SEQ ID NO:7 maps to human chromosome 19q13.1-13.2. Northern blot analysis of primary human T lymphocytes reveals a 2.2 Kb transcript that is highly upregulated in activated compared to resting cells, concomitant with an increase in KCa currents. Transcripts of the novel gene presented herein are detected by PCR and by Northern blots in placenta, prostate, thymus, spleen, colon and many cell lines of hematopoietic origin. Patch-clamp recordings of SEQ ID NO:2 -transfected HEK 293 cells reveal a large

voltage-independent, inwardly rectifying potassium current with a single channel conductance of  $33 \pm 2$  pS in symmetrical potassium solutions. The channel is activated by intracellular calcium ( $K_d = 270 \pm 8$  nM) with a highly cooperative interaction of  $\sim 3$  calcium ions per channel. hKCa4 currents are blocked by externally applied tetraethylammonium ( $K_d = 30 \pm 6$  mM),  
5 charybdotoxin ( $K_d = 10 \pm 1$  nM) and clotrimazole ( $K_d = 380 \pm 34$  nM), but are resistant to apamin, iberiotoxin, kaliotoxin and scyllatoxin ( $K_d > 1$   $\mu$ M). These properties of the cloned channel are very similar to those reported for the KCa channel in activated human T lymphocytes (Grissmer *et al*, J. Gen. Physiol. 102: 601, 1993), sequences described herein (e.g., SEQ ID NO:7; SEQ ID NO:1; SEQ ID NO:2; and SEQ ID NO:3) are example embodiments which  
10 encode this native channel type.

A 2.2 Kb cDNA clone described herein (e.g., SEQ ID NO: 7) has been identified from a human lymph node lambda library using a probe derived from an EST sequence that was identified from a database search for novel potassium channels. The clone has 400 bp of 5' UTR, 1.3 Kb of coding region and 540 bp of 3' UTR (fig. 17A). The entire transcript was  
15 mapped of the predominant 2.2 Kb band seen in Northern blots from various tissues (see fig 18). Four independent lambda clones were identified which started with the same 5' UTR sequence ( $\pm 15$  bp) and the 3' UTR ended with a poly-adenylation signal followed by a poly-A tail. We also detected an in-frame stop codon upstream of the initiator ATG in all the 4  
20 clones, ruling out the existence of alternate upstream initiator ATG's. The novel gene (e.g., SEQ ID NO:2) maps to chromosome 19q13.1-13.2. The translated protein comprises 427 amino acids. Hydropathy plots reveal a short intracellular N-terminus, followed by 6 transmembrane segments (S1-S6) and a long intracellular C-terminus (Fig. 17 B). The loop between S5 and S6 contains the highly conserved GYG sequence characteristic of all cloned  
25 potassium channel pores (10). The channel protein (SEQ ID NO:3) has one consensus N-linked glycosylation site between S5 and the pore and several sites for serine-threonine phosphorylation. There are no consensus E-F hand motifs in SEQ ID NO:3. A comparison of the amino acid sequence of the hKCa4 with other representative members of the  $K^+$  channel super-family reveals that it is most similar to the small-conductance KCa channels (hSK1,  
30 rSK2, rSK3; SEQ ID NO: 4, SEQ ID NO: 5 and SEQ ID NO: 6) recently cloned from brain,

although it shares only 40-42% amino acid identity to them; warranting placement within a distinct sub-family (fig.15 C).

The novel gene described herein (e.g., SEQ ID NO:2) maps to human chromosome 19q13.1-

- 5 13.2. This region of the chromosome 19 has several disease susceptibility loci, including B-cell lymphoma and Diamond Black Fan Syndrome. Abnormalities of the gene described herein may be associated with any disorders which map to this general area including disorders that are currently being mapped as part of the human genome project.
- 10 Northern blot analysis reveals a strong signal at ~2.2 Kb in activated human T lymphocytes, placenta, prostate, and colon (fig. 16A and B). Moderate signal was observed in spleen, thymus and peripheral blood leukocytes (fig. 16 A); there was no detectable signal in brain. Minor transcripts of larger sizes were also seen in some tissues (fig. 16 A).
- 15 The novel gene described herein (SEQ ID NO:2) when co-transfected with a reporter gene for GFP into HEK 293 cells, produced a calcium-dependent  $K^+$  current with strong inward rectification in symmetrical  $K^+$  solutions upon perfusion of 1  $\mu$ M ionomycin (Fig. 13 A). Current induced by ionomycin was completely blocked by 100 nM CTX, but unaffected (i.e., <10% change in amplitude) by 1  $\mu$ M apamin ( $n = 2$ ). Internal dialysis with a solution
- 20 containing 1  $\mu$ M  $[Ca^{2+}]_{free}$  also activated a stable large current ( $15.4 \pm 2.4$  nS slope conductance between -100 and -30 mV;  $n = 11$ ) with similar characteristics. Pharmacological evaluation (Fig. 13 B) revealed that this current was potently inhibited by CTX ( $K_d = 10 \pm 1$  nM) and also blocked by TEA ( $K_d = 30 \pm 7$  mM) and Clotrimazole ( $K_d = 387 \pm 34$  nM), but insensitive to apamin, iberiotoxin, and kaliotoxin, and scyllatoxin ( $K_d > 1$   $\mu$ M,  $n=3$  for each
- 25 experiment). Control current during ramps in symmetrical 160 mM  $K^+$  converged with current elicited during CTX perfusion close to 0 mV (e.g., Fig. 13 A), supporting  $K^+$  selectivity of the channel.  $K^+$  selectivity was evaluated further by examining reversal potentials of control vs. CTX currents over a range of different  $K^+_{(out)}$  concentrations (Fig. 13 C). Reversal potential shifted -57 mV per e-fold increase in  $K^+_{(out)}$  in close agreement with
- 30 the predicted Nernstian value for a  $K^+$ -selective channel. Untransfected HEK 293 cells, or cells transfected with GFP alone, showed small currents in response to voltage ramps during

dialysis with  $1 \mu\text{M Ca}^{2+}_{\text{(free)}}$  (conductance  $< 0.1 \text{ nS}$ ;  $n = 8$ ), and these currents were not measurably affected by  $100 \text{ nM CTX}$ .

Cell-attached recordings revealed single channel openings during perfusion with Ringer  
5 solution containing  $1 \mu\text{M}$  ionomycin. Channels showed a unitary conductance (measured during voltage ramps between  $-120$  and  $-30 \text{ mV}$  of  $33 \pm 2 \text{ pS}$  ( $n = 3$ ; Fig. 14A) with pipettes containing  $160 \text{ mM K}^+$ , and  $9 \pm 1 \text{ pS}$  with pipettes containing  $4.5 \text{ mM K}^+$ . Pipette potentials between  $-100$  and  $+20 \text{ mV}$  had no apparent effect on the probability of channel openings ( $P_o$ ) during ramps (e.g., Fig. 16A) or steps. The inside-out configuration in symmetrical  $\text{K}^+$  also  
10 revealed single channels of similar conductance with gating that clearly depended on the  $[\text{Ca}^{2+}]_{\text{free}}$  at the cytoplasmic face of the patch (Fig. 14 B). The  $P_o$  with different  $[\text{Ca}^{2+}]_{\text{free}}$  varied considerably between cells, but never exceeded  $\sim 0.5$ . Fitting the open probability vs.  $[\text{Ca}^{2+}]_{\text{free}}$  curve revealed an activation  $K_d = 270 \pm 8 \text{ nM Ca}^{2+}$  with a Hill-coefficient of  $2.7 \pm 0.2$ , indicative of a highly cooperative interaction between calcium ions and the channel (Fig.  
15 14 C).

The evidence demonstrated herein that the novel gene (e.g., SEQ ID NO:2) encodes the calcium-activated potassium channel in T lymphocytes is summarized as follows: (a) the clone was identified from a lymph node library which is enriched for activated T cells (b)  
20 Northern blot analysis of resting and activated T cells reveals  $\sim 10$ -fold up-regulation of transcript levels, corresponding to an increase in current levels (FIG.16 B) (c) the biophysical and pharmacological properties of the currents are very similar to those reported for T cells (Grissmer, S., et al., 102:601 (1993); Verheugen, J. A. H., *et al.*, Cell Calcium 21:1 (1997)), including inward rectification in high K solution, block by CTX and TEA, single-channel  
25 conductance, and  $K_d$  for calcium-dependence (FIG.14 D). Furthermore, SEQ ID NO:2 quite possibly may encode the KCa channels found in other cells of hematopoietic origin, including erythrocytes, monocytes, B cells, neutrophils and eosinophils. The reported electrophysiological properties of the channel from these cells is compatible with the data presented herein, for example SEQ ID NO:2 heterologously expressed in HEK 293 cells.  
30 Moreover, a 1Kb band was amplified using SEQ ID NO:7-specific primers, for example, and confirm its sequence from many cell-lines of hematopoietic origin (e.g.: U937, HL-60, fetal

liver cells, Jurkats). Subtle changes in reported properties is attributable to differences in post-translational modifications and/or association with accessory subunits in different cell-types.

- 5 The significant up-regulation of the novel transcript levels during T cell activation highlights the importance of these channels in early-activation and post-activation immune responses. Engagement of the T cell receptor by mitogen or antigen evokes an increase in intracellular calcium concentrations, leading to membrane hyperpolarizations (rather than depolarization) caused by opening of KCa channels. This event in turn maximizes the electrical driving force  
10 for calcium influx through calcium-release activated channels, facilitating sustained calcium oscillations required for cell proliferation, cytokine production and the expression of immune function. Blockade of both KCa and Kv channels have been reported to cause profound inhibition of T cell proliferation (Rader, R. K., *et al.*, J. Immunol. 156:1425 (1996)) while blockade of KCa channels is sufficient to prevent secondary immune responses (Verheugen, J.  
15 A. H., *et al.*, Cell Calcium, 21:1 (1997)).

**The calcium activated potassium channel subunit (SEQ ID NO:3) (hKCA4) transcript is present in microglial cells**

- Microglia are resident macrophages in the brain that play vital roles in neuroinflammatory  
20 responses and in neurodegenerative diseases (Wood, P.L., *Role of CNS macrophages in neurodegeneration in Neuroinflammation: mechanism and management*, edited by PL Wood, Human Press Inc., Totowa, NJ). In view of SEQ ID NO:2 expressed in several cells of hematopoietic origin, including monocytes and macrophages, we asked whether it could be present in microglia. A predicted 800 bp band was obtained from microglial cDNA, and  
25 sequence was confirmed, indicating the presence of the transcript in brain microglial cells.  
*See EXAMPLE XIX.*

**Calmodulin is an interaction partner for the novel calcium activated potassium channel subunit (hkca4) and is possibly the calcium sensor**

- 30 The small- and intermediate-conductance calcium-activated potassium channels KCa2, KCa3, KCa4 belong to a distinct sub-family of K<sup>+</sup> channels that play crucial roles in

hyperpolarization of excitable and non-excitable cells. Despite their exquisite calcium sensitivity, the protein coding region of these channels do not contain any consensus calcium-bowl or EF-hand motifs. To investigate whether an accessory protein is involved in calcium sensing, we used the 146-a.a. C-tail of hKCa4 as a bait to screen an activated T-cell library in  
5 a yeast-two-hybrid system (discussed *infra*). Several of the positive clones encoded for calmodulin. See, Vogel, H., The Merck Frosst Award Lecture 1994, *Calmodulin: a versatile calcium mediator protein.*, Biochem Cell Biol., 72(9-10):357 (1994).

GST-fusion pull-down experiments have confirmed the interaction and deletion analysis have  
10 identified the first 97 amino acids in the C-tail of hKCa4 as being required for calmodulin binding (FIG.21). This region is highly conserved among all SKCa channels and appears to contain hydrophobic and positively charged amino acids representative of calmodulin binding domains. The interaction of calmodulin with the C-tail does not appear to require calcium and suggests that calmodulin could be pre-bound to the channel on the membrane. The SEQ ID  
15 NO:3 calcium activated potassium channel (hKCa4) currents in HEK 293 cells are inhibited (50-60% at 10 uM) by the calmodulin antagonist W7. Taken together, our data point to the fact that calmodulin quite possibly may serve as the calcium sensor for this class of channels.

### Embodiments

- 20 The present invention relates to nucleic acid (e.g., SEQ ID NO:2) and amino acid sequences (SEQ ID NO:3) of a novel human lymph node derived calcium activated potassium channel and variations thereof and to the use of these sequences to identify compounds that modulate the activity of human leukocytes, as described *infra*.
- 25 The invention further relates to the use of the novel calcium activated potassium channel in expression systems as assays for agonists or antagonists of calcium-activated potassium channels. The invention also relates to the diagnosis, study, prevention, and treatment of disease related to dysfunctional leukocytes.
- 30 Polynucleotide sequences which encode the human leukocyte calcium activated potassium channel SEQ ID NO:3 or a functionally equivalent derivative thereof may be used in



accordance with the present invention which comprise deletions, insertions and/or substitutions of the SEQ ID NO:2 nucleic acid sequence. Biologically active and/or pharmacologically active variants of the calcium activated potassium channel subunit of the present invention may also be comprised of deletions, insertions or substitutions of SEQ ID  
5 NO:3 amino acid residues. A purified polynucleotide comprising a nucleic acid sequence encoding the polypeptide having the sequence substantially as depicted in SEQ ID NO:3 or a biologically active fragment thereof is a particularly preferred embodiment of the present invention.

10 Amino acid substitutions of SEQ ID NO:3 may be made, for instance, on the basis of similarity in polarity, charge, solubility, hydrophobicity, hydrophilicity, and/or the amphipathic nature of the residues as long as the biological activity of the calcium activated potassium channel is retained. For example, negatively charged amino acids include aspartic acid and glutamic acid; positively charged amino acids include lysine and arginine; and amino  
15 acids with uncharged polar head groups having similar hydrophilicity values include leucine, isoleucine, valine; glycine, alanine; asparagine, glutamine; serine, threonine phenylalanine, and tyrosine.

Nucleic acid sequences which encode the amino acid sequence of the calcium activated  
20 potassium channel described herein are of an exponential sum due to the potential substitution of degenerate codons (different codons which encode the same amino acid). The oligonucleotide sequence selected for heterologous expression is therefore preferably tailored to meet the most common characteristic tRNA codon recognition of the particular host expression system used as well known by those skilled in the art.

25

Suitable conservative substitutions of amino acids are known to those of skill in this art and may be made without altering the biological activity of the resulting polypeptide, regardless of the chosen method of synthesis. The phrase "conservative substitution" includes the use of a chemically derivatized residue in place of a non-derivatized residue provided that such  
30 polypeptide displays the desired binding activity. D-isomers may also be substituted for the

naturally occurring amino acids. Substitutions are preferably, although not exclusively, made in accordance with those set forth in TABLE 1 as follows:

5

TABLE 1

Original residue	Example conservative substitution embodiments
Ala (A)	Gly; Ser; Val; Leu; Ile; Pro
Arg (R)	Lys; His; Gln; Asn
Asn (N)	Gln; His; Lys; Arg
Asp (D)	Glu
Cys (C)	Ser
Gln (Q)	Asn
Glu (E)	Asp
Gly (G)	Ala; Pro
His (H)	Asn; Gln; Arg; Lys
Ile (I)	Leu; Val; Met; Ala; Phe
Leu (L)	Ile; Val; Met; Ala; Phe
Lys (K)	Arg; Gln; His; Asn
Met (M)	Leu; Tyr; Ile; Phe
Phe (F)	Met; Leu; Tyr; Val; Ile; Ala
Pro (P)	Ala; Gly
Ser (S)	Thr
Thr (T)	Ser
Trp (W)	Tyr; Phe
Tyr (Y)	Trp; Phe; Thr; Ser
Val (V)	Ile; Leu; Met; Phe; Ala

The nucleotide sequences of the present invention may also be engineered in order to alter a  
 10 coding sequence for a variety of reasons, including but not limited to, alterations which

modify the cloning, processing and/or expression of the gene product. For example, mutations may be introduced using techniques which are well known in the art, eg, site-directed mutagenesis to insert new restriction sites, to alter glycosylation patterns, to change codon preference, etc.

5

Dominant negative mutant embodiments of the channel subunit are contemplated herein (*see*; FIG.15 for site change embodiments (sites for glycosylation are indicated by a triangle, serine-threonine phosphorylation sites are circled) (e.g., change to D or R), *any* changes in the pore region are preferred embodiments as dominant negative mutants) as pharmacologically

10 valuable agents in order to attenuate the corresponding channel biological activity *in vivo*.

Generally mutated versions, or site directed mutagenesis species embodiments of the pore region of SEQ ID NO:2, for instance, can be used to produce non-functional calcium activated potassium channels to be subsequently over-expressed in heterologous cells, preferably ones  
15 of hematopoietic origin, to knock-out endogenous channels by dominant-negative suppression of the channel. Ribera, A.B., J. Neurosci., 16:1123, (1996); Babila, T., *et al.*, Neuron, 12:615 (1994).

Included within the scope of the present invention are alleles of the calcium activated  
20 potassium channel gene of the present invention. As used herein, an allele or allelic sequence is an alternative form of the subunit herein. Alleles result from nucleic acid mutations and mRNA splice-variants or longer transcripts which comprise SEQ ID NO:2 which produce polypeptides whose structure or function may or may not be altered.

25 Particularly preferred embodiments of the present invention are host cells transformed with a purified polynucleotide comprising a nucleic acid sequence encoding the polypeptide having the sequence substantially as depicted in SEQ ID NO:3 or a biologically active fragment thereof. Cells of this type or preparations made from them may be used to screen for pharmacologically active modulators of the novel calcium activated potassium channel  
30 activity using methods which are well known in the art. U.S. Patent No. 5,397,702, *Assay For and Treatment of Autoimmune Diseases*, issued March 14, 1995; U.S. Patent No. 5,637,470,

*Screening array using cells expressing recombinant alpha and beta subunits of the mammalian large-conductance (maxi-K) potassium channel*, issued June 10, 1997; U.S. Patent No. 5,607,843, *Nucleic Acid Sequence Encoding Apamin Binding Protein*, issued March 4, 1997; and U.S. Patent No. 5,602,169, *3-substituted oxindole derivatives as potassium channel modulators*, issued Feb. 11, 1997.

Such modulators are potentially useful in treating disease which is manifested by dysfunctional leukocytes including but not limited to acute and chronic inflammation, asthma, allergies, graft rejection, proliferative disorders, anemias, neurodegenerative diseases with immunological components, as well as autoimmune diseases including rheumatoid arthritis, type-1 diabetes mellitus, multiple sclerosis, myasthenia gravis, systemic lupus erythematosus, Sjogren's syndrome, mixed connective tissue disease, and experimental allergic encephalomyelitis (EAE).

As used herein, a "functional derivative" of the calcium activated potassium channel disclosed herein is a compound that possesses a biological activity (either functional or structural) that is substantially similar to SEQ ID NO:3. The term "functional derivatives" is intended to include the "fragments," "variants," "degenerate variants," "analogs" and "homologues", and to "chemical derivatives". The term "variant" is meant to refer to a molecule substantially similar in structure and function to either the entire calcium activated potassium channel molecule or to a fragment thereof. A molecule is "substantially similar" to the calcium activated potassium channel polypeptide if both molecules have substantially similar structures or if both molecules possess similar biological activity. The term "analog" refers to a molecule substantially similar in function to either the entire calcium activated potassium channel polypeptide, or to a fragment thereof.

The cloned calcium activated potassium channel DNA obtained through the methods described herein may be recombinantly expressed by molecular cloning into an expression vector containing a suitable promoter and other appropriate transcription regulatory elements, and transferred into prokaryotic or eukaryotic host cells to produce the calcium activated potassium channel polypeptide. Techniques for such manipulations are fully described in

Sambrook, J., et al., Molecular Cloning Second Edition, 1990, Cold Spring Harbor Press and are well known in the art.

Expression vectors are defined herein as DNA sequences that are required for the transcription  
5 of cloned copies of genes and the translation of their mRNAs in an appropriate host cell.  
Such vectors can be used to express nucleic acid sequences in a variety of hosts such as  
bacteria, bluegreen algae, plant cells, insect cells, fungal cells, human, and animal cells.  
Specifically designed vectors allow the shuttling of DNA between hosts such as bacteria-  
yeast, or bacteria-animal cells, or bacteria-fungal cells, or bacteria-invertebrate cells.

10

#### ***In vitro* synthesis of capped mRNA**

The full-length cDNA (SEQ ID NO:1), may be used in standard procedures to synthesize  
biologically active mRNA for functional expression in heterologous cells, in *Xenopus* oocytes,  
for instance, or in various types of mammalian cells. See, e.g., Goldin A, Methods Enzymol.  
15 207:279 (1992).

The coding region for the novel calcium activated potassium channel described herein (e.g., a  
region which comprises SEQ ID NO:2) which encodes the polypeptide having a bio-active  
sequence substantially as depicted in SEQ ID NO:3 or a biologically active fragment thereof.  
20 may be cloned 3', for example, to a bacteriophage promoter, e.g., an SP6, T7 or T3 promoter.  
The PGEM vectors from Promega, Madison, WI, are examples of preferred vectors which  
may be used with the present invention. Standard vectors known in the art, such as pSP64T or  
pBSTA, which enhance stability of the message and increase specific expression in *Xenopus*  
oocytes, are especially preferred. Kreig and Melton, Nucleic Acids Res., 12:7057 (1984).  
25 These particular preferred vectors contain the *Xenopus* beta-globin 5' and 3' untranslated  
mRNA regions flanking the 5'end and a poly-A tail on the 3' end of the gene.

A plasmid vector DNA construct which contains an insert which encodes the novel calcium  
activated potassium channel described herein or a biologically active fragment thereof (e.g., a  
30 SEQ ID NO:1 or SEQ ID NO:7 region which comprises SEQ ID NO:2, or a truncated version  
thereof) may be then cut with a restriction enzyme to linearize the construct 3' to the structural

coding region. The linearized vector should be extracted with phenol-chloroform-isoamyl alcohol, precipitated with ethanol and re-suspended in RNase-free water for use as a transcription template. The transcription reaction, for example, may be carried out as described *infra* to synthesize biologically active mRNA for subsequent *in vivo* or *in vitro* translation by methods which are well-known in the art. The mMessage mMachine™ kit, AMBION, Austin, TX, is especially preferred for synthesizing biologically active mRNA. Alternately, the mRNA transcripts may be used as probes for analysis of tissue distribution by Northern analyses and or RNase protection assays.

#### 10 **mRNA SYNTHESIS**

PGEM, Promega, Madison, WI

	10x SP6/T7 buffer	5 ul
	10X ATP, CTP, UTP (5 mM each)	5 ul
	10x GTP (5 mM)	1 ul
15	10x GpppG (Cap; 5 mM)	5 ul
	DTT (1M)	0.5 ul
	Rnase inhibitor (40U/ul)	1.25 ul
	Water	27 ul
	Linearized DNA (1 ug/ul)	1-5 ul
20	SP6 or T7 RNA polymerase (20U/ul)	1-3 ul

TOTAL REACTION VOLUME 50 ul

---

Incubate at 37 degrees for 1-2 hours. Add 1 ul of RNase-free Dnase to degrade template DNA. Digest for 10 minutes. Add 75 ul of water. Extract with phenol/CHCl<sub>3</sub>. Ethanol precipitate RNA. Store in aliquots at -70 degrees.

#### **Functional expression**

Biologically active mRNA can be introduced into heterologous cells for functional expression and analyses by methods well-known in the art. Synthetic mRNA from example constructs described *supra*, for example, may be injected into *Xenopus* oocytes for functional expression

and analyses. Goldin, A., *Methods Enzymol.*, 207:266, (1992). Heterologous calcium-activated potassium channels may be examined using standard two-electrode voltage clamp techniques. *See, e.g.*, Stuhmer, W., *Methods in Enzymol.*, 207:319, (1992); Kohler, *et al.*, *Science*, 273:1709 (1996). Calcium concentrations inside the cell may be altered, for example, by adding Ionomycin, a calcium-ionophore, co-expression with a G-protein coupled receptor that causes a rise in intracellular calcium. Calcium concentrations may alternately be altered by pulling inside-out patches and changing calcium concentrations in the bath medium. Grissmer, S., *et al.*, *Calcium-activated Potassium Channels in Resting and Activated Human T Lymphocytes*, *J. Gen. Physiol.*, 102:601 (1993). Standard biophysical parameters, such as activation, calcium dependence, single-channel conductance, inactivation, tail currents, potassium selectivity, and thorough pharmacology of various K channel blockers including TEA, Charybdotoxin, Apamin, and others may also be tested. Grissmer, S., *et al.*, *Calcium-activated Potassium Channels in Resting and Activated Human T Lymphocytes*, *J. Gen. Physiol.*, 102:601 (1993).

Alternatively, cRNA (synthetic mRNA from a cDNA construct) can be introduced into heterologous mammalian cells, for example, RBL cells (ATCC # CRL 1378), and 293 cells (ATCC # CRL 1573), may be transformed using standard art-methods. For example, the Eppendorf microinjection system may be used (Micromanipulator 5171 and Transjector 5242). Transformed cells may be analysed for calcium-activated K currents about 4 hours later using patch-clamp techniques which are well-documented. *E.g.*, Ikeda, *et al.*, *Pflug. Arch. Eur. J. Physiol.*, 422:201 (1992); Grissmer, *et al.*, *J. Gen. Physiol.*, 102:601 (1993).

#### **Over-expression of the novel channel in cell-lines**

Transient and/or stable eucaryotic transfectant cells comprised of the coding region(s) described herein are contemplated for high-level expression of the novel calcium activated potassium channel.

Eucaryotic transfectants are preferred embodiments of the present invention for employment in pharmacological target binding studies for the identification molecules which block the novel channel described herein *in vivo*. HEK cells are particularly preferred.

Transient expression of coding regions for the novel channel can be achieved by straight transfection into mammalian cells, by standard techniques. Omari, K. *et al.*, J. Physiol., 499:369, (1997); Panyi, G. *et al.*, J. Gen. Physiol., 107(3):409 (1996). High level transient  
5 expression may be achieved using standard viral systems, e.g., Baculovirus, Adenovirus, or Vaccinia virus. Channel numbers resulting from these systems are typically 5-500K per cell. Kamb, A., Methods Enzymol. 207:423 (1992); Sun, T. *et al.*, Biochemistry, 33(33):9992 (1994); Spencer, R.H., *et al.*, J. Biol. Chem., 272:2389 (1997). See, EXAMPLE XV, Baculovirus expression system for hKCa4 (SEQ ID NO:2).

10

Stable transfection of heterologous cells using sequences which encode the novel calcium activated potassium channel described herein (SEQ ID NO:3) or biologically active variations or fragments thereof can be generated using, for example, NIH-3t3, L929, COS, HEK, or CHO cells. See, e.g., EMBO, 11(6):2033 (1992); Grissmer, *et al.*, Mol. Pharm., 45:1227  
15 (1994).

A preferred vector for use with the present invention is pcDNA/Neo, which is commercially available from INVITROGEN, Carlsbad, CA.

20 Cells, NIH-3t3, for example, are grown to 50% confluency in 60mm plates (media, and conditions are according to requirements of the particular cell line) and transfected with 5 ug of pure DNA comprising a coding region for the novel calcium activated potassium channel, e.g. SEQ ID NO:2, in pcDNA/Neo using the Lipofection reagent, as described by the supplier (LIFE TECHNOLOGIES Gibco BRL, Bethesda, MD). After transfection, the cells  
25 are incubated at 37°C, conditions for 3 days in medium with 10% FCS. Cells are trypsinized seeded onto 100mm dishes, and then selected with 300ug/ml of G418 (Neomycin). Only cells that have stable integration of the heterologous coding region will grow in the presence of G418, which is conferred by the Neomycin-resistance gene in the plasmid. Isolated clones are processed for 2-3 rounds of purification and subjected to patch-clamp analysis for Kca  
30 currents.



- A variety of mammalian expression vectors may be used to express the recombinant calcium activated potassium channel polypeptide disclosed herein in mammalian cells. Commercially available mammalian expression vectors which are suitable for recombinant expression, include but are not limited to, pcDNA3 (Invitrogen), pMC1neo (Stratagene), pXT1
- 5 (Stratagene), pSG5 (Stratagene), EBO-pSV2-neo (ATCC 37593) pBPV-1(8-2) (ATCC 37110), pdBPV-MMTneo(342-12) (ATCC 37224), pRSVgpt (ATCC 37199), pRSVneo (ATCC 37198), pSV2-dhfr (ATCC 37146), pUCTag (ATCC 37460), and IZD35 (ATCC 37565).
- 10 A variety of insect cell expression vectors may be used to express recombinant calcium activated potassium channel polypeptide in insect cells along with the alpha subunit. Commercially available insect cell expression vectors which may be suitable for recombinant expression of beta subunit include but are not limited to pBlue Bac III (INVITROGEN).
- 15 Eukaryotic recombinant host cells are especially preferred. Examples include but are not limited to yeast, mammalian cells including but not limited to cell lines of human, bovine, porcine, monkey and rodent origin, and insect cells including but not limited to Drosophila and silkworm derived cell lines. Cell lines derived from mammalian species which may be suitable and which are commercially available, include but are not limited to, L cells L-
- 20 M(TK-) (ATCC CCL 1.3), L cells L-M (ATCC CCL 1.2), 293 (ATCC CRL 1573), Raji (ATCC CCL 86), CV-1 (ATCC CCL 70), COS-1 (ATCC CRL 1650), COS-7 (ATCC CRL 1651), CHO-K1 (ATCC CCL 61), 3T3 (ATCC CCL 92), NIH/3T3 (ATCC CRL 1658), HeLa (ATCC CCL 2), C1271 (ATCC CRL 1616), BS-C-1 (ATCC CCL 26) and MRC-5 (ATCC CCL 171).
- 25
- The calcium activated potassium channel may also be expressed as a recombinant protein with one or more additional polypeptide domains added to facilitate protein purification. Such purification facilitating domains include, but are not limited to, metal chelating peptides such as histidine-tryptophan modules that allow purification on immobilized metals (Porath, J.,
- 30 Protein Expr Purif., 3:263 (1992)), protein A domains that allow purification on immobilized immunoglobulin, and the domain utilized in the FLAGS extension/affinity purification system

(Immunex Corp, Seattle WA). The inclusion of a cleavable linker sequences such as Factor XA or enterokinase (Invitrogen, San Diego CA) between the purification domain and TMP is useful to facilitate purification.

**5 Glutathione S-Transferase (GST) Fusion Proteins, *in vitro* Translation of Calmodulin, and Binding of Calmodulin to GST Fusion Proteins**

GST fusion plasmids were based on PGEX-6P-1 (Pharmacia, Piscataway, NJ) and were constructed by conventional means using standard PCR with primer to the C-terminal of KCa4 sequence with the appropriate restriction site, EcoRI site and XhoI site in this case.

- 10 Positions 1252-1272 of SEQ ID NO:1 is the sense oligonucleotide and SEQ ID NO:1 Positions 1660-1680 is the antisense oligonucleotide for the full length hKCa4 C-terminal GST construct. After transfection into *Escherichia coli* BL21, synthesis of the fusion proteins was induced with 0.2 mM isopropyl b-D-thiogalactoside (IPTG) in a liquid culture grown to OD of 1.0. After 2-3 hr at 37°C, the cells were collected by centrifugation, resuspended in
- 15 NETN lysis buffer (0.5% Nonidet P-40/1 mM EDTA/20 mM Tris-HCl, pH 8.0/100 mM NaCl; 1.0 ml per 20 ml of culture) and lysed by sonication. The lysate was cleared by centrifugation at 10,000 x g for 10 min at 4°C. The GST fusion proteins in the supernatant were adsorbed for 30 min at room temperature to glutathione-agarose beads in NETN [1 volume of lysate:1 volume of 50% (vol/vol) slurry of agarose-GSH beads (Pharmacia) in
- 20 NETN]. The last wash was with binding buffer [1% (vol/vol) Polyoxyethylene 9 Lauryl Ether/2 mM EDTA/100 mM NaCl/20 mM Tris-HCl, pH 8.01 instead of NETN]. Slurries [50% (vol/vol)] slurries of washed agarose-GSH beads with GST fusion proteins adsorbed to them (agarose-GSH::GST fusions) were then incubated for 30 min at room temperature in an equal volume of binding buffer with <sup>35</sup>S-labeled calmodulin, synthesized by coupled
- 25 transcription-translation (Cat#L4610, Promega, Madison, WI), as described (Pragnell, *et al.*, Nature, 368:67 (1994). After a 50-fold dilution in binding buffer, the beads were centrifuged and washed three times with binding buffer and resuspended in an equal volume of 4x Laemmli's sample buffer. Proteins released from the beads were analyzed by 12% SDS/PAGE followed by autoradiography to detect retention of calmodulin by the fragments
- 30 fused to GST. *In vitro*-translated luciferase was used to test for specificity in the interaction of calmodulin with the KCa4 C-terminal fragments.

**Stable mammalian cell expression of hkca4 for assay development**

Heterologous expression of cloned channels is a pre-requisite for assay development, particularly preferred binding assays which require high density of pure channels. The first  
5 step towards assay development for the novel channel subunit (hKCa4) (SEQ ID NO:3) was therefore the generation of a high-level expression system. A stable transfectant of the novel gene described herein was developed in HEK 293 cells.

SEQ ID NO:2 was subcloned into mammalian cell expression vectors to generate both stably  
10 integrated and episomally replicating lines. For stable lines, the (pGEN IRES-neo) construct and the pcDNA3 construct is used (described in JBC, 272 (52): 32723). The pGEN IRES-neo vector was made by cleaving a 1.3 kb SmaI/ScaI fragment containing the coding region (SEQ ID NO:2) of the calcium activated potassium channel. This fragment was subcloned into the vector at the SmaI site. This cloning strategy results in a construct that uses the authentic  
15 initiator methionine.

For the episomal lines, SEQ ID NO:2 was first subcloned into pFastBac1 from the pcDNA3 construct using a strategy that eliminates the three upstream fusion amino acids. A 1101 bp  
20 NcoI fragment was cleaved from the pcDNA3 construct and cloned into the NcoI site of pBlueBac III (Invitrogen, Carlsbad, CA). A colony was identified whose orientation contained the BamHI site from the pBlueBac III multiple cloning region on the 5' end of the cDNA (hKCa4) fragment. The 5' end of the cDNA was then cleaved from pBlueBac III using BamHI (5') and EcoRV (internal to hKCa4). The 3' end of the gene was provided by cleaving the pcDNA3 construct with EcoRV and XhoI (3'). The reconstructed coding region was  
25 cloned into the BamHI (5') and XhoI (3') sites of the pFastBac1 vector (Life Technologies, Gaithersburg, MD). This cloning strategy eliminated the N-terminal MGA created by cloning into pcDNA3. The coding region of the calcium activated potassium channel subunit was then cleaved from this pFastBac1 construct using KpnI (5') and NotI (3'), and cloned into the KpnI and NotI sites of pCEP4 (Invitrogen, Carlsbad, CA).

HEK293 cells were cultured in DMEM medium containing 10% fetal bovine serum, 100 units per mL of penicillin and 100 µg per mL of streptomycin. CHO-K1 cells were cultured in Ham's F12 medium containing 10% fetal bovine serum. Transfections into HEK293 cells were achieved using the LipoTAXI mammalian transfection kit (cat # 204110) from

5 Stratagene (La Jolla, CA), and, separately, a modified calcium phosphate transfection protocol (see EXAMPLE XIII). LipoTAXI transfections were according to the protocol, with the following notes: 1 mL of serum-free DMEM was mixed with 70 µL transfection reagent, and 15 µg of either (pGEN IRES-neo) or pcDNA3 - hKCa4 construct was added. The DNA/lipid complexes were allowed to form for 20 minutes at room temperature. The cells were

10 incubated with the DNA/lipid complexes for four hours at 37°C, 5% CO<sub>2</sub>. Ten mL of DMEM containing 20% fetal bovine serum (FBS) were added and the cells were grown overnight. Calcium phosphate transfections followed the attached protocol (EXAMPLE XIII). The episomal lines in HEK293 cells were generated using the pCEP4/SEQ ID NO:2 construct and the modified calcium phosphate method according to the protocol (EXAMPLE XIII).

15 Transfections into CHO-K1 cells were achieved using Cellfectin reagent from Life Technologies, Gaithersburg, MD, (cat # 10362-010), according to EXAMPLE XIV.

For all lines, media was removed and replaced with fresh the morning after the transfection.

20 Forty eight hours later (i.e., three days after the transfection), cells were split as follows: HEK293 stable lines, 1:10; HEK293 episomal lines, 1:4; CHO-K1 stable lines, 1:4. Twenty four hours after the split (i.e., four days after the transfection), stable lines were subjected to medium containing 1 mg/mL G418 sulfate (MediaTech, Inc., Herndon, VA, cat # 30-234-CR; or Life Technologies, Gaithersburg, MD, cat # 10131-035) and the HEK293 episomal lines

25 were subjected to medium containing 250µg/mL hygromycin B (Calbiochem, San Diego, CA, cat # 400051). The HEK293 episomal lines were not cloned, but rather were maintained as four separate populations. For all stable lines, colonies were picked after 10 days in selective medium using sterile cloning discs (Research Products International, Mount Prospect, IL, cat. # 248710, according to their protocol (see attached protocol #3) Clones were expanded as per

30 protocol and characterized by Northern blot analysis and <sup>125</sup>I-CTX binding.

**RNA dot-blot assay for screening positive clones**

Initial characterizations were done by seeding HEK293 stable or episomal lines in two duplicate 96 well plates and growing to confluency. One plate was processed for crude total RNA by rinsing the cells with 250  $\mu$ L of phosphate buffered saline and lysing the cells

5 directly in the well with 100  $\mu$ L per well of RLT buffer plus  $\beta$ -mercaptoethanol from Qiagen, Inc., Chatsworth, CA (cat # 79216). The lysates were vacuumed onto Hybond N+ positively charged nylon membrane (Amersham Life Science, Inc., Arlington Heights, IL, cat # RPN203B) that had been soaked in 2 X concentrated SSC solution (20 X SSC is 0.3 M  $\text{Na}_3$  citrate and 3 M NaCl at pH 7.0). The wells were rinsed 2 to 3 times in 2 X SSC and the blot

10 was crosslinked on a Spectrolinker XL1500 (Spectronics Corp., Westbury, NY) under the optimal setting. The vacuum apparatus used was the Bio-dot 96 well apparatus, Bio-Rad Laboratories, Hercules, CA, cat. # 170-6545. The probe used to hybridize the blot was a PCR product corresponding to nucleotide # 658-1098 of hKCa4. Fifty ng of fragment was labeled using  $\alpha^{32}\text{P}$ -dCTP (NEN Life Science Products, Inc., Boston, MA, cat. # NEG-013H) and

15 Ready-To-Go labeling beads (-dCTP) (Pharmacia Biotech, Piscataway, NJ, cat # 27-9240-01) according to the manufacturer's recommendations. Labeled product was purified from unincorporated radioactive nucleotide using Probe Quant G-50 Micro Columns (Pharmacia Biotech, Piscataway, NJ, cat # 27-5335-01) according to the manufacturer's recommendations. The blot was prehybridized in ExpressHyb solution (Clontech, Palo Alto,

20 CA, cat # 8015-1) for 1.5 hours at 68°C. Twenty five ng ( $28 \times 10^6$  cpm) of probe was denatured by boiling and added to 10 mL of fresh prewarmed ExpressHyb. This replaced the prehybridization solution for one hour at 68°C. The blot was washed twice for twenty minutes at room temperature in 2 X SSC and 0.1% SDS, followed by two fifteen-minute washes at 50°C in 0.1 X SSC and 0.1% SDS. The blot was then exposed to film overnight

25 with two intensifying screens.

**CTX-Binding Assay For Screening Positive Clones**

30 The second plate was processed for  $^{125}\text{I}$ -CTX binding by rinsing each well in 250  $\mu$ L of phosphate buffered saline and detaching the cells using 50  $\mu$ L per well of versene (Life

Technologies, Gaithersburg, MD, cat # 15040-066). The cells were pipetted to break up clumps and to ensure they were detached, and the cells were then centrifuged five minutes at 1000 rpm to bring the cells to the bottom of the wells. Supernatant was carefully removed and the cell pellets were stored at -20°C until the binding assay was performed. To begin the  
5 experiment, plates were allowed to thaw to room temperature and 200 ul of assay buffer (5mM NaCl, 5mM KCl, 10mM HEPES, 6mM glucose, pH 8.4) containing 0.02% bovine serum albumin was added to each well. 50ul of <sup>125</sup>I-charybdotoxin (Dupont New England Nuclear, Boston, MA, Cat # NEX-276) was added to each well for a final concentration of 50 pM. The plates were incubated for 1 hour at room temperature with gentle agitation on a  
10 Titertek plate mixer. After one hour, bound ligand was separated from free by filtration onto GF/C Unifilters (Packard, Meriden, CT) that had been pre-soaked in 0.6% polyethylenimine. The samples were washed twice rapidly, using ice cold wash buffer (200mM NaCl, 20mM HEPES, pH 8.0). The Unifilter plates were air-dried overnight. When dry, 25 ul of Microscint-20 was added to the Unifilters, and the plates were counted in a Top Count  
15 scintillation counter (Packard, Meriden, CT).

Cell lines that appeared positive for both SEQ ID NO:2 RNA transcript and <sup>125</sup>I-CTX binding were characterized further by whole cell patch clamp analysis.

#### 20 Patch clamp recording of stable transfectants:

Currents were recorded with an Axopatch 200A amplifier (Axon Instr., Foster City, CA) using the whole cell and inside-out configurations (Hammill, *et al.*, Pflugers Arch., 391:85 (1981)). Thin-wall borosilicate glass pipettes were fabricated, sylgarded, and fire-polished to a DC resistance of 2 - 8 MΩ. The resistance of patch seals was >10 GΩ. Liquid junction potentials  
25 were corrected for in all experiments, and series resistance compensation of >70% was used where maximal current was >0.5 nA. Voltage clamp protocols were implemented and data acquisition performed with pClamp 6.0 software (Axon Instr., CA). Currents were low-pass filtered at (-3 db at 1 kHz) and then digitized at 3-8 kHz as computer files with a TL-1 interface (Scientific Solutions, Solon, OH). Currents were measured with p-Clamp software,  
30 and iterative curve fittings were performed with either p-Clamp or Origin software (Version 3.73; Microcal Inc., Northampton, MA).

For *whole cell recording*, the pipette solution contained (in mM) 160 K aspartate, 2 MgCl<sub>2</sub>, 5 HEPES, and 1.6 EGTA with either 0.8 CaCl<sub>2</sub> (calculated  $[Ca^{2+}]_{free} = 100 \text{ nM}$ ; Eqcal software, Biosoft Corp, Cambridge, UK) or 1.6 CaCl<sub>2</sub> (calculated  $[Ca^{2+}]_{free} = 1 \text{ }\mu\text{M}$ ) at pH 7.2 (by KOH) and osmolality ~315 mOsm (by sucrose). Cells were perfused locally with a solution containing 160 KCl, 2 CaCl<sub>2</sub>, 1 MgCl<sub>2</sub>, 5 HEPES, and 5 glucose at pH 7.4 (by KOH) and osmolality ~325 mOsm. In K<sup>+</sup> selectivity experiments, equimolar Na<sup>+</sup> was substituted for K<sup>+</sup>. Whole-cell voltage clamp protocol evoked current by 273 ms voltage ramp from -100 to +40 mV. Slope conductance was measured during voltage ramp between -99 to -50.1 mV. Slope conductance of Charybdotoxin block at 100 or 300 nM was subtracted from control current as leak.

Channel numbers were calculated using the formula :  
 (Whole-cell conductance / Unitary conductance) / (open probability of 0.5). Cell-attached voltage clamp protocol evoked single-channel current by 500 ms voltage ramp from -120 to +180 mV. Unitary conductance was measured during voltage ramp between -99 to +20 mV. Several stable transfectant clones expressed huge KCa currents. A current trace from an example clone (HC23) is shown in FIG.17. This current was blocked by CTX (<sup>125</sup>I charybdotoxin binding assay for KCa4 (SEQ ID NO:3) (*See examples I, II and III*)). An example of a dose-response curve for CTX-binding on KCa4 (SEQ ID NO:3) channels is shown in FIG.18.

#### **Human calcium activated potassium channel subunit (hKCa4) expression in anergic T-cells**

In the event a T-lymphocyte interacts with a specific antigen and results in signaling through the T-cell antigen receptor, yet does not lead to T cell proliferative response, the cell enters a state of nonresponsiveness, or *anergy*. Johnson, J.G., *et al.*, Life Sciences, 55: 1767 (1994). For normal T cell activation, IL-2 synthesis and proliferation, two signals are required: one through the T cell receptor (TCR) complex and another through a co-stimulatory molecule on the T cell known as CD28. Anergy induction can result from a number of different situations, including antigen presentation in the absence of costimulation, pharmacological blockade of

T-cell proliferation or chronic stimulation of the TCR by antigen. Anergy is a long-lived state characterized by a profound inability of the T-cell to produce IL-2. An understanding of the mechanism of anergy induction and maintenance will most certainly lead to beneficial clinical applications, including improving graft acceptance and avoidance of such deleterious immune responses such as autoimmunity and allergy. See EXAMPLE XVII.

We have shown dramatic upregulation of the novel calcium activated potassium channel subunit RNA and currents when T-cells are activated. There appears to be a correlation between induction of hKCa4 (SEQ ID NO:3) currents and the anergic state.

10

Resting (no stimulus), activated (anti-CD3+CD28) and T-cells induced to anergy (anti-CD3 alone) were patch clamped, to examine relative KCa4 channel activities. Dramatic downregulation of KCa4 current is herein demonstrated in T-cells receiving TCR occupancy alone compared to T cells receiving TCR occupancy plus costimulation. ( $p < 0.001$ ). Results are summarized in FIG.20. Our data point to the importance of the costimulatory signal in upregulation of this channel during T cell activation.

### Antibodies

Three peptides were designed to SEQ ID NO:3 (hKCa4) for antibody generation. These peptides were chosen based on specificity for KCa4 over KCa1-3 as well as high antigenicity index.

A: Positions 415-427 + C of SEQ ID NO:3 (COOH-terminus with an added cysteine)

B: Positions 415-427 of SEQ ID NO:3 (bare COOH-terminus)

25 C: Positions 135-148 of SEQ ID NO:3 (S3-S4 loop)

The peptides were ordered from Research Genetics, Huntsville, AL. The peptides were conjugated to KLH and to BSA using the Inject Immunogen EDC Conjugation kit (Cat#77101, Pierce Inc., IL). The KLH-conjugated peptides were sent to VETERINARY MEDICINE for antibody generation in rabbits. 300 ug of the KLH-peptide was used per injection. The immunization protocol was as follows:



	Aug 19	Prebleed
	Aug 20	1 <sup>st</sup> Injection(Complete adj.)
	Aug 27	2 <sup>nd</sup> Injection(Incomplete adj.)
5	Sept. 10	3 <sup>rd</sup> Injection(Incomplete adj.)
	Sept. 17	1 <sup>st</sup> Bleed
	Sept. 24	4 <sup>th</sup> Injection(Incomplete adj.)
	Oct. 1	2 <sup>nd</sup> Bleed
	Oct. 8	5 <sup>th</sup> Injection(Incomplete adj.)
10	Oct 15	3 <sup>rd</sup> Bleed
	Oct. 22	6 <sup>th</sup> Injection(Incomplete adj.)
	Oct. 29	Terminal Sample

An ELISA (ELISAmate kit, Kirkegaard and Perry Laboratories Inc. MD) using the BSA-  
 15 peptide as the antigen, and HRP-conjugated goat anti-rabbit antibody as the secondary  
 reagent, gave very good titers in immune sera for all three peptide antibodies tested. Western  
 blot analysis (using the ECL detection kit, Amersham, England) of a C-terminus tail GST  
 fusion protein (made in bacteria, see elsewhere in this document under *two yeast hybrid*)  
 picked up a strong band at the predicted size when peptide B antiserum was used but not  
 20 with peptide C antiserum (which is directed to S3-S4 loop). See, FIG.19.

In addition, recombinant the calcium activated potassium channel polypeptide can be  
 separated from other cellular proteins by use of an immunoaffinity column made with  
 monoclonal or polyclonal antibodies specific for full length nascent calcium activated  
 25 potassium channel polypeptide, or polypeptide fragments of the calcium activated potassium  
 channel.

Calcium activated potassium channel polypeptides described herein may be used to affinity  
 purify biological effectors from native biological materials, e.g., disease tissue. Affinity  
 30 chromatography techniques are well known to those skilled in the art. A calcium activated  
 potassium channel peptide described herein or an effective fragment thereof, is fixed to a solid

matrix, e.g. CNBr activated Sepharose according to the protocol of the supplier (Pharmacia, Piscataway, NJ), and a homogenized/buffered cellular solution containing a potential molecule of interest is passed through the column. After washing, the column retains only the biological effector which is subsequently eluted, e.g., using 0.5M acetic acid or a NaCl  
5 gradient.

Monospecific antibodies to the calcium activated potassium channel polypeptide are purified from mammalian antisera containing antibodies reactive against the polypeptide or are prepared as monoclonal antibodies reactive with calcium activated potassium channel  
10 polypeptide using the technique of Kohler and Milstein. Nature 256, 495-497 (1975). Monospecific antibody as used herein is defined as a single antibody species or multiple antibody species with homogenous binding characteristics for the novel calcium activated potassium channel. Homogenous binding as used herein refers to the ability of the antibody species to bind to a specific antigen or epitope, such as those associated with the novel calcium activated  
15 potassium channel, as described. Novel calcium activated potassium channel specific antibodies are raised by immunizing animals such as mice, rats, guinea pigs, rabbits, goats, horses and the like, with rabbits being preferred, with an appropriate concentration of calcium activated potassium channel polypeptide either with or without an immune adjuvant.

20 Preimmune serum is collected prior to the first immunization. Each animal receives between about 0.1 mg and about 1000 mg of calcium activated potassium channel polypeptide associated with an acceptable immune adjuvant. Such acceptable adjuvants include, but are not limited to, Freund's complete, Freund's incomplete, alum-precipitate, water in oil emulsion containing Corynebacterium parvum and tRNA. The initial immunization consists of calcium  
25 activated potassium channel polypeptide in, preferably, Freund's complete adjuvant at multiple sites either subcutaneously (SC), intraperitoneally (IP) or both. Each animal is bled at regular intervals, preferably weekly, to determine antibody titer. The animals may or may not receive booster injections following the initial immunization. Those animals receiving booster injections are generally given an equal amount of the antigen in Freund's incomplete  
30 adjuvant by the same route. Booster injections are given at about three week intervals until maximal titers are obtained. At about 7 days after each booster immunization or about weekly

after a single immunization, the animals are bled, the serum collected, and aliquots are stored at about - 20o C.

Monoclonal antibodies (mAb) reactive with the calcium activated potassium channel  
5 polypeptide are prepared by immunizing inbred mice, preferably Balb/c, with the calcium activated potassium channel polypeptide. The mice are immunized by the IP or SC route with about 0.1 mg to about 10 mg, preferably about 1 mg, of calcium activated potassium channel polypeptide in about 0.5 ml buffer or saline incorporated in an equal volume of an acceptable adjuvant, as discussed above. Freund's complete adjuvant is preferred. The mice receive an  
10 initial immunization on day 0 and are rested for about 3 to about 30 weeks. Immunized mice are given one or more booster immunizations of about 0.1 to about 10 mg of calcium activated potassium channel polypeptide in a buffer solution such as phosphate buffered saline by the intravenous (IV) route. Lymphocytes, from antibody positive mice, preferably splenic lymphocytes, are obtained by removing spleens from immunized mice by standard procedures  
15 known in the art. Hybridoma cells are produced by mixing the splenic lymphocytes with an appropriate fusion partner, preferably myeloma cells, under conditions which will allow the formation of stable hybridomas. Fusion partners may include, but are not limited to: mouse myelomas P3/NS1/Ag 4-1; MPC-11; S-194 and Sp 2/0, with Sp 2/0 being preferred. The antibody producing cells and myeloma cells are fused in polyethylene glycol, about 1000  
20 molecular weight, at concentrations from about 30% to about 50%. Fused hybridoma cells are selected by growth in hypoxanthine, thymidine and aminopterin supplemented Dulbecco's Modified Eagles Medium (DMEM) by procedures known in the art. Supernatant fluids are collected from growth positive wells on about days 14, 18, and 21 and are screened for antibody production by an immunoassay such as solid phase immunoradioassay (SPIRA)  
25 using calcium activated potassium channel polypeptide as the antigen. The culture fluids are also tested in the Ouchterlony precipitation assay to determine the isotype of the mAb. Hybridoma cells from antibody positive wells are cloned by a technique such as the soft agar technique of MacPherson, Soft Agar Techniques, in Tissue Culture Methods and Applications, Kruse and Paterson, Eds., Academic Press, 1973.

Monoclonal antibodies are produced in vivo by injection of pristane primed Balb/c mice, approximately 0.5 ml per mouse, with about  $2 \times 10^6$  to about  $6 \times 10^6$  hybridoma cells about 4 days after priming. Ascites fluid is collected at approximately 8-12 days after cell transfer and the monoclonal antibodies are purified by techniques known in the art.

5

In vitro production of anti- calcium activated potassium channel polypeptide mAb is carried out by growing the hybridoma in DMEM containing about 2% fetal calf serum to obtain sufficient quantities of the specific mAb. The mAb are purified by techniques known in the art.

10

Antibody titers of ascites or hybridoma culture fluids are determined by various serological or immunological assays which include, but are not limited to, precipitation, passive agglutination, enzyme-linked immunosorbent antibody (ELISA) technique and radioimmunoassay (RIA) techniques. Similar diagnostic assays are used to detect the presence  
15 of the novel calcium activated potassium channel polypeptide in body fluids or tissue and cell extracts.

Diagnostic assays using calcium activated potassium channel polypeptide specific antibodies are useful for the diagnosis of conditions, disorders or diseases characterized by abnormal  
20 expression of the calcium activated potassium channel polypeptide or expression of genes associated with abnormal cell growth. Diagnostic assays for calcium activated potassium channel polypeptide include methods utilizing the antibody and a label to detect calcium activated potassium channel polypeptide in human body fluids, cells, tissues or sections or extracts of such tissues. The polypeptides and antibodies of the present invention may be used  
25 with or without modification. Frequently, the polypeptides and antibodies will be labeled by joining them, either covalently or noncovalently, with a reporter molecule, a myriad of which are well-known to those skilled in the art.

A variety of protocols for measuring calcium activated potassium channel polypeptide, using  
30 either polyclonal or monoclonal antibodies specific for the respective protein are known in the art. Examples include enzyme-linked immunosorbent assay (ELISA), radioimmunoassay

(RIA) and fluorescent activated cell sorting (FACS). A two-site, monoclonal-based immunoassay utilizing monoclonal antibodies reactive to two non-interfering epitopes on calcium activated potassium channel polypeptide is preferred, but a competitive binding assay may be employed. These assays are described, among other places, in Maddox, DE et al  
5 (1983, J Exp Med 158:1211).

In order to provide a basis for the diagnosis of disease, normal or standard values for calcium activated potassium channel polypeptide expression must be established. This is accomplished by combining body fluids or cell extracts taken from normal subjects, either  
10 animal or human, with antibody to calcium activated potassium channel polypeptide under conditions suitable for complex formation which are well known in the art. The amount of standard complex formation may be quantified by comparing it with a dilution series of positive controls where a known amount of antibody is combined with known concentrations of purified calcium activated potassium channel polypeptide. Then, standard values obtained  
15 from normal samples may be compared with values obtained from samples from subjects potentially affected by a disorder or disease related to calcium activated potassium channel polypeptide expression. Deviation between standard and subject values establishes the presence of the disease state.

20 It is readily apparent to those skilled in the art that the above described methods for producing monospecific antibodies may be utilized to produce antibodies specific for calcium activated potassium channel polypeptide fragments, or the full-length nascent calcium activated potassium channel polypeptide. Specifically, it is readily apparent to those skilled in the art that monospecific antibodies may be generated which are specific for the fully functional  
25 receptor or fragments thereof.

Calcium activated potassium channel polypeptide antibody affinity columns are made by adding the antibodies to Affigel-10 (Biorad), a gel support which is activated with N hydroxysuccinimide esters such that the antibodies form covalent linkages with the agarose  
30 gel bead support. The antibodies are then coupled to the gel via amide bonds with the spacer arm. The remaining activated esters are then quenched with 1M ethanolamine HCl (pH 8).

The column is washed with water followed by 0.23M glycine HCl (pH 2.6) to remove any non-conjugated antibody or extraneous protein. The column is then equilibrated in phosphate buffered saline (pH 7.3) with appropriate detergent and the cell culture supernatants or cell extracts containing calcium activated potassium channel polypeptide made using appropriate  
5 membrane solubilizing detergents are slowly passed through the column. The column is then washed with phosphate buffered saline/detergent until the optical density falls to background, then the protein is eluted with 0.23M glycine-HCl (pH 2.6)/detergent. The purified calcium activated potassium channel polypeptide is then dialyzed against phosphate buffered saline/detergent.

10

### Screening methods

The present invention is also directed to methods for screening for compounds which modulate the expression of DNA or RNA encoding calcium activated potassium channel polypeptide, as well as the function of calcium activated potassium channel polypeptide *in*  
15 *vivo*. Compounds which modulate these activities may be DNA, RNA, peptides, proteins, or non-proteinaceous organic molecules. Compounds may modulate by increasing or attenuating the expression of DNA or RNA encoding calcium activated potassium channel polypeptide, or the function of calcium activated potassium channel polypeptide. Compounds that modulate the expression of DNA or RNA encoding calcium activated potassium channel polypeptide or  
20 the function of calcium activated potassium channel polypeptide may be detected by a variety of assays. The assay may be a simple "yes/no" assay to determine whether there is a change in expression or function. The assay may be made quantitative by comparing the expression or function of a test sample with the levels of expression or function in a standard sample.

25 The calcium activated potassium channel polypeptide described herein, its immunogenic fragments or oligopeptides can be used for screening therapeutic compounds in any of a variety of drug screening techniques. The fragment employed in such a test may be free in solution, affixed to a solid support, borne on a cell surface, or located intracellularly. The abolition of activity or the formation of binding complexes, between calcium activated  
30 potassium channel polypeptide and the agent being tested, may be measured. Accordingly, the present invention provides a method for screening a plurality of compounds for specific

binding affinity with calcium activated potassium channel polypeptide or a fragment thereof, comprising providing a plurality of compounds; combining the calcium activated potassium channel polypeptide of the present invention or a fragment thereof with each of a plurality of compounds for a time sufficient to allow binding under suitable conditions; and detecting  
5 binding of the calcium activated potassium channel polypeptide, or fragment thereof, to each of the plurality of compounds, thereby identifying the compounds which specifically bind the calcium activated potassium channel polypeptide.

Since the novel gene, e.g., SEQ ID NO:1 is highly expressed in activated T-cells, and since  
10 T-cell calcium activated potassium channels are potently blocked by the scorpion toxin Charybdotoxin (CTX), various cell-lines that heterologously over-expresses the novel channel structural coding regions described herein can be used in radio-labeled binding assays to screen for pharmacologically active molecules which block the novel channel using  $^{125}\text{I}$ -CTX as the ligand.

15

**CTX-binding assay for high-throughput screening for modulators of the novel channel**

Particularly preferred embodiments of the present invention are cell-lines that heterologously over-expresses the calcium activated potassium channel described herein (e.g., SEQ ID NO:3 or a biologically active truncated version thereof or a chimeric fusion) and their use in binding  
20 assays to identify pharmacologically active molecules which block or otherwise modulate biological activity of the novel channel.

For example, a radio-labeled binding assay using  $^{125}\text{I}$ -CTX as the ligand may be used as previously described by Hill, R.J., Mol. Pharm., 48:98 (1995), and Deutsch, C., *et al.*, J. Biol.  
25 Chem., 266:3668 (1991). Membrane preparations of cell-lines which over-express the novel calcium activated potassium channel described herein are made by homogenizing the cells using a Polytron for 25 seconds at 13,000 RPM and spun at low speed (100 g) for 2 minutes. The supernatant is spun at high speed (50,000 g) for 10 minutes. The pellet is suspended in 1 ml of assay buffer (5 mM NaCl, 5 mM KCL, 10 mM HEPES, 6 mM glucose, pH 8.4) and  
30 diluted to 50 ug/ml.

To each well of a 96-well microtiter plate, 130  $\mu$ l of assay buffer is added, as well as 20  $\mu$ l of test molecule compound (test drug may be, for instance, a small molecule, peptide, analog or mimetic compound) /control assay buffer /non-specific (10 nM cold CTX) (not-labelled) 50  $\mu$ l of membranes from cells which over-express the novel calcium activated potassium channel  
5 at 50  $\mu$ g/ml and 50  $\mu$ l of  $^{125}$ I-CTX (25 pM; NEN, 2200 Ci/mmol) are incubated for 20 minutes at 21°C with mixing. Bound radiolabeled CTX is separated from free radiolabeled CTX in solution by filtering over pre-soaked GF/C Unifilters (Packard Instruments) and washing rapidly in ice-cold wash buffer. Upon drying, the filter plates are scintillation counted. Data from saturation experiments are subject to Scatchard analysis and linear  
10 regression. Deutsch, *et al.*, J. Biol. Chem., 266:3668 (1991). Compounds that compete with the radio-labeled CTX for binding the novel calcium activated potassium channel are identified which produce a reduction in specific counts.

#### Scintillation proximity assay

15 In another embodiment, a scintillation proximity assay (SPA) that eliminates the need for filters, can be easily adapted for high throughput screening (HTS) assays. Hoffman, R., *et al.*, Anal. Biochem., 20370 (1992); Kienhuis, C.B.M., *et al.*, J. Recept. Res., 12:389 (1992). See, EXAMPLE XVI.

#### 20 Fluorometric imaging plate reader (FLIPR)

The FLIPR system is most useful as a high-throughput primary screen for ion channels, which influence the cell membrane potential (e.g. KCa4 (SEQ ID NO:3) described herein). FLIPR is a commercial device (Molecular Devices, Sunnyvale, CA) which is designed to measure fluorescent signals simultaneously from 96 wells in a microtiter plate using an argon ion laser  
25 and a semi-confocal optical system. It is equipped with a 96-tip pipettor allowing simultaneous drug administration to all wells. The time resolution is down to 1 second. The fluorescent dyes which are most useful are DiBAC4(3) for membrane potential (Epps, D.E., *et al.*, Chem Phys Lipids., 69(2):137 (1994)) and Fluo-3 for calcium (Rijkers, G.T., *et al.*, Cytometry, 11(8):923 (1990). Cells stably expressing SEQ ID NO:2 are plated onto 96-well  
30 plates at a density of  $0.5 \times 10^4$  cells/well and allowed to form a uniform monolayer. Next day, cell media is removed and incubated with EBSS buffer containing 20 mM HEPES buffer and



5uM DiBAC (cat# B-438, Molecular Probes) with or without test-compounds and incubated at 37 degrees for 1.5 hours. Cells are then washed in a non-fluorescent salt-buffer such as EBSS. Fluorescence is measured using FLIPR laser optics. (Molecular Devices, Sunnyvale, CA)

5

Methods of identifying compounds that modulate the activity of a human leukocyte calcium-activated potassium channel polypeptide are generally preferred, which comprise

- (a) combining a candidate compound modulator of a human leukocyte calcium-activated potassium channel activity with a polypeptide of a human calcium-activated potassium
- 10 channel having the sequence substantially as depicted in SEQ ID NO:3, and
- (b) measuring an effect of the candidate compound modulator on the channel.

Methods of identifying compounds that modulate the activity of a human leukocyte calcium-activated potassium channel, are especially preferred which comprise

- 15 (a) combining a candidate compound modulator of a human leukocyte calcium-activated potassium channel activity with a host-cell expressing the polypeptide of a human calcium-activated potassium channel having the sequence substantially as depicted in SEQ ID NO:3, and
- (b) measuring an effect of the candidate compound modulator on the channel.

20

Compounds which are identified generally according to methods described and referenced herein that modulate the activity of a calcium-activated potassium channel comprised of the sequence substantially as depicted in SEQ ID NO:3 are especially preferred embodiments of the present invention.

25

### **Two hybrid assays**

- In another embodiment of the invention, a nucleic acid sequence which encodes a calcium activated potassium channel substantially as depicted in SEQ ID NO:3 or a biologically active fragment thereof may be ligated to a heterologous sequence to encode a fusion protein. For
- 30 example, for screening compounds for modulation of biological activity, it may be useful to encode a chimeric calcium activated potassium channel as described herein for expression in

hererologous host cells. Chimeric constructs may also be used to express a 'bait', according to methods well known using a yeast two-hybrid system, to identify accessory native peptides that may be associated with the calcium activated potassium channel polypeptide described herein. Fields, S., *et al.*, Trends Genet., 10:286 (1994); Allen, J.B., *et al.*, TIBS, 20:511  
5 (1995). A yeast two-hybrid system has been described wherein protein:protein interactions can be detected using a yeast-based genetic assay via reconstitution of transcriptional activators. Fields, S., Song, O., Nature 340:245 (1989). The two-hybrid system used the ability of a pair of interacting proteins to bring a transcription activation domain into close proximity with a DNA-binding site that regulates the expression of an adjacent reporter gene.  
10 See also, Mendelsohn, A.R., Brent, R., Curr. Op. Biotech., 5:482 (1994); Phizicky, E.M., Fields, S., Microbiological Rev., 59(1):94 (1995); Yang, M., *et al.*, Nucleic Acids Res., 23(7):1152 (1995); Fields, S., Sternglanz, R., TIG, 10(8):286 (1994); and US Patents 5,283,173, *System to Detect Protein-Protein Interactions*, and 5,468,614, which are incorporated herein by reference.

15

#### **Yeast 2-Hybridization experiments (implication of calmodulin)**

The MATCHMAKER two-hybrid system from Clontech, Palo Alto CA, (Cat # K1604-1) from is a complete GAL4-based system that provides transcriptional assay for detecting protein-protein interactions in yeast. The two-hybrid system utilizes a powerful growth  
20 selection-the conditional expression of a nutritional reporter gene-to screen large numbers of yeast transformants with a specially constructed fusion library for interacting proteins. hKCa4 (SEQ ID NO:2) C-terminal was subcloned into GAL4 DNA-BD vector (pAS2-1, Clontech, Palo Alto, CA) by conventional means using standard PCR (Clontech, PT3061-1). The restriction sites used in the subcloning are EcoRI and XhoI. This construct was then  
25 confirmed by sequencing and used as the bait in the library screening. The library used is human leukocyte MATCHMAKER cDNA (HL4021AB, Clontech). All the screening procedures were performed based on manufacturer's recommendations (Clontech manual #PT3061-1). Ten thousand putative positive clones were identified after first round screening. Tow thousand of these were then subjected to colony-lift lacZ assay, which  
30 identified seven positives. Positive colonies are then sequenced for further analysis, such as comparison with the sequence database. Six of the seven clones coded for calmodulin.

**Antisense**

The cDNA sequences SEQ ID NO:1 and SEQ ID NO:7 provided herein, may be used in another embodiment of the invention to study the physiological relevance of the novel calcium  
5 activated potassium channel in cells, especially cells of hematopoietic origin, by knocking out the endogenous gene by use of anti-sense constructs that target, for example, the region of the initiator methionine or the pore-region described *supra*. See, e.g., Shirihai, O. *et al.*, *Pfugers Arch.* 431:632 (1996); Chung, S., *et al.*, *PNAS*, 92:5955 (1995).

10 U.S. Patent No. 5,639,595, *Identification of Novel Drugs and Reagents*, issued Jun. 17, 1997, wherein methods of identifying oligonucleotide sequences that display *in vivo* activity are thoroughly described, is herein incorporated by reference. Expression vectors containing random  
oligonucleotide sequences derived from previously known polynucleotides are transformed  
15 into cells.

The cells are then assayed for a phenotype resulting from the desired activity of the oligonucleotide.

Once cells with the desired phenotype have been identified, the sequence of the oligonucleotide having the desired activity can be identified. Identification may be  
20 accomplished by recovering the vector or by polymerase chain reaction (PCR) amplification and sequencing the region containing the  
inserted nucleic acid material.

Nucleotide sequences that are complementary to the calcium activated potassium channel  
25 polypeptide encoding polynucleotide sequence can be synthesized for antisense therapy. These antisense molecules may be DNA, stable derivatives of DNA such as phosphorothioates or methylphosphonates, RNA, stable derivatives of RNA such as 2'-O-alkylRNA, or other oligonucleotide mimetics. Calcium activated potassium channel polypeptide antisense molecules may be introduced into cells by microinjection, liposome encapsulation or by  
30 expression from vectors harboring the antisense sequence. Calcium activated potassium

channel polypeptide antisense therapy may be particularly useful for the treatment of diseases where it is beneficial to reduce calcium activated potassium channel polypeptide activity.

Calcium activated potassium channel polypeptide gene therapy may be used to introduce the polypeptide into the cells of target organs. Conversely, calcium activated potassium channel polypeptide *antisense* gene therapy may be used to reduce the expression of the polypeptide in the cells of target organs. The calcium activated potassium channel polypeptide coding region can be ligated into viral vectors which mediate transfer of the calcium activated potassium channel polypeptide DNA by infection of recipient host cells. Suitable viral vectors include retrovirus, adenovirus, adeno-associated virus, herpes virus, vaccinia virus, polio virus and the like. Alternatively, calcium activated potassium channel polypeptide DNA can be transferred into cells for gene therapy by non-viral techniques including receptor-mediated targeted DNA transfer using ligand-DNA conjugates or adenovirus-ligand-DNA conjugates, lipofection membrane fusion or direct microinjection. These procedures and variations thereof are suitable for *ex vivo*, as well as *in vivo* calcium activated potassium channel polypeptide gene therapy.

### Compositions

Pharmaceutically useful compositions comprising calcium activated potassium channel polypeptide DNA, calcium activated potassium channel polypeptide RNA, or calcium activated potassium channel polypeptide, or modulators of calcium activated potassium channel polypeptide activity, may be formulated according to known methods such as by the admixture of a pharmaceutically acceptable carrier. Examples of such carriers and methods of formulation may be found in *Remington's Pharmaceutical Sciences*. To form a pharmaceutically acceptable composition suitable for effective administration, such compositions will contain an effective amount of the protein, DNA, RNA, or modulator.

Therapeutic or diagnostic compositions of the invention are administered to an individual in amounts sufficient to treat or diagnose calcium activated potassium channel polypeptide related disorders. The effective amount may vary according to a variety of factors such as the individual's condition, weight, sex and age. Other factors include the mode of administration.

The pharmaceutical compositions may be provided to the individual by a variety of routes such as subcutaneous, topical, oral and intramuscular.

- 5 The term "chemical derivative" describes a molecule that contains additional chemical moieties which are not normally a part of the base molecule. Such moieties may improve the solubility, half-life, absorption, etc. of the base molecule. Alternatively the moieties may attenuate undesirable side effects of the base molecule or decrease the toxicity of the base molecule. Examples of such moieties are described in a variety of texts, such as *Remington's*  
10 *Pharmaceutical Sciences*.

Compounds identified according to the methods disclosed herein may be used alone at appropriate dosages defined by routine testing in order to obtain optimal modulation of a calcium activated potassium channel, or its activity while minimizing any potential toxicity. In  
15 addition, co-administration or sequential administration of other agents may be desirable.

The present invention also has the objective of providing suitable topical, oral, systemic and parenteral pharmaceutical formulations for use in the novel methods of treatment of the present invention. The compositions containing compounds identified according to this  
20 invention as the active ingredient for use in the modulation of calcium activated potassium channels can be administered in a wide variety of therapeutic dosage forms in conventional vehicles for administration. For example, the compounds can be administered in such oral dosage forms as tablets, capsules (each including timed release and sustained release formulations), pills, powders, granules, elixirs, tinctures, solutions, suspensions, syrups and  
25 emulsions, or by injection. Likewise, they may also be administered in intravenous (both bolus and infusion), intraperitoneal, subcutaneous, topical with or without occlusion, or intramuscular form, all using forms well known to those of ordinary skill in the pharmaceutical arts. An effective but non-toxic amount of the compound desired can be employed as a calcium activated potassium channel modulating agent.

The daily dosage of the products may be varied over a wide range from 0.01 to 1,000 mg per adult human/per day. For oral administration, the compositions are preferably provided in the form of scored or unscored tablets containing 0.01, 0.05, 0.1, 0.5, 1.0, 2.5, 5.0, 10.0, 15.0, 25.0, and 50.0 milligrams of the active ingredient for the symptomatic adjustment of the dosage to the patient to be treated. An effective amount of the drug is ordinarily supplied at a dosage level of from about 0.0001 mg/kg to about 100 mg/kg of body weight per day. The range is more particularly from about 0.001 mg/kg to 10 mg/kg of body weight per day. Even more particularly, the range varies from about 0.05 to about 1 mg/kg. Of course the dosage level will vary depending upon the potency of the particular compound. Certain compounds will be more potent than others. In addition, the dosage level will vary depending upon the bioavailability of the compound. The more bioavailable and potent the compound, the less compound will need to be administered through any delivery route, including but not limited to oral delivery. The dosages of the calcium activated potassium channel modulators are adjusted when combined to achieve desired effects. On the other hand, dosages of these various agents may be independently optimized and combined to achieve a synergistic result wherein the pathology is reduced more than it would be if either agent were used alone.

## EXAMPLES

### 20 EXAMPLE I

#### Membrane Preparation of Transfected HEK Cells for Screening Compounds

Membranes are prepared when the HEK cells in each T-150 (or T-175 or T-225) flask are  
25 confluent.

The preparation of membranes from each flask of cells is as follows:

1. Pour off the cell culture media and save (since it always contains some cells that have just been dislodged from handling the flask).

- 57 -

2. Rinse the confluent layer of cells with about 8 mls of phosphate buffered saline (PBS), 37°C. This removes residual fetal calf serum which can inhibit the trypsin used in step 3. Save the PBS from the rinse, too, to capture any other cells that are dislodged.
- 5 3. Add about 5 ml of trypsin-EDTA, 37°C. Slowly rock the flask back and forth to let the trypsin wash over the whole layer of cells. In about 30 sec, some of the cells will come loose; tap the flask firmly on the lab bench to dislodge the rest of the cells. Quickly add fresh medium to the flask to stop the trypsin reaction. Add about 10ml medium for every 5ml of trypsin. Pipet the trypsin/cells/medium into a centrifuge tube. Rinse the flask with about 5 ml  
10 more medium and combine with other tubes.
4. Wash the cells by spinning at 250xg, 8 min. Resuspend the individual cell pellets in a small volume of assay buffer, and combine them into one tube.
- 15 5. Homogenize the cells using a Polytron, 25 sec., 13000 RPM
6. LOW SPEED SPIN: 800 RPM (~100xg) 2 min.; save the supernatant; discard the pellet
- 20 7. HIGH SPEED SPIN: transfer the sup into high speed tubes, and spin 20,000 RPM (50,000xg), 10 min. Discard the supernatant; add 1 ml of assay buffer on top of pellet
8. Freeze on dry ice; cap the tubes and store at -80 C.
- 25 Dilute the membranes to 50ug/ml, and use 2.5ug in each well (50ul of 50ug/ml)

## EXAMPLE II

- 30 <sup>125</sup>-I Charybdotoxin binding assay for (SEQ ID NO:3) HKca4

## (A) FILTRATION METHOD

**Materials:**Assay Buffer:

5mM NaCl

5mM KCl

10mM HEPES

6mM glucose

230mM sucrose

bring to pH 8.4 with TRIS base

add BSA to a [final] of 0.02%

Wash Buffer:

200mM NaCl

20mM HEPES

pH 8.0 with Tris base

SEQ ID NO:2 stably transfected whole MEL cells 500,000 cell per well - or - SEQ ID  
 5 NO:2 stably transfected HEK cells

misc. 96 well microtiter plates, Falcon, round bottom, non tissue culture treated  
 Phosphate buffered saline (PBS) no Ca<sup>++</sup> or Mg<sup>++</sup>  
 96 well Unifilters, GF/C, (Packard, Meriden, CT).  
 10 Microscint-20 scintillation fluor, Packard  
 Top Count scintillation counter, Packard  
 Polyethylenimine, 0.6% for soaking Unifilters  
 125I-Charybdotoxin, (NEN, Boston, MA, NEX-276) 2200Ci/ mmol  
 nonspecific binding defined with 100nM cold ChTX (Bachem, King of  
 15 Prussia, PA).

**Assay:** 130ul assay buffer  
 20ul test compound or buffer or cold CTX (nonspecific)  
 20 50ul cells at 1E<sup>6</sup>/ml (500,000 per well)  
50ul 125I-ChTX [final] 50 pM



250ul total volume

**Assay Method:**

Place 130ul assay buffer per well in 96well plate. Add 20ul of test compound, or control  
5 solution (more buffer or test compound diluent), or nonspecific (defined with 100nM cold  
CTX). Wash the cells twice in PBS, resuspend in assay buffer, and add 50ul of cells to the  
assay plate. Add 50ul  $^{125}\text{I}$ -ChTX, [final] 50 pM, and incubate for 45 min. 21° C, mixing on  
the plate mixer. Separate bound from free  $^{125}\text{I}$ -ChTX by filtering over the presoaked GF/C  
Unifilters; wash twice rapidly, using ice cold wash buffer. Let Unifilter plates dry overnight.  
10 When dry, add 25ul Microscint-20 to Unifilters, and count in Top Count scintillation counter.

An example of a dose-response curve for CTX-binding/novel calcium activated potassium  
channel subunit (hKCa4) is shown in FIG.18.

15

**EXAMPLE III**

**$^{125}\text{I}$ -ChTX Binding Assays**

Lymphocytes are incubated in 12 x 75-mm polystyrene tubes with  $^{125}\text{I}$ -ChTX ( $1-2 \times 10^3$   
20 Ci/mmol). Unless otherwise noted, cells are suspended in isotonic sucrose medium (Medium  
I) containing 10 mM NaHepes, 5mM KCl, 5 mM NaCl, and 6 mM glucose, pH 8.4, and  
incubated with  $^{125}\text{I}$ -ChTX for 1 h at room temperature on a rotary shaker. Nonspecific  
binding is determined in the presence of 10 nM native ChTX. Stock cell suspensions are  
diluted to give a final cell concentration of  $2.5 \times 10^5$  cells/ml in a total volume of 400  $\mu\text{l}$ . At  
25 the end of the incubation period, samples are diluted with 4 ml of ice-cold Quench solution,  
which contain 200 mM NaCl, 20 mM Hepes (free acid), titrated to pH 8.0 with Trisbase.  
Quenched samples are filtered through GF/C glass microfiber filters, that had been presoaked  
in 0.6% polyethylenimine, and washed twice with ice-cold Quench solution. Triplicate  
samples are run for each experimental point. Standard deviation of the mean is typically less  
30 than 5%. Different cell preparations may produce somewhat different ratios of

nonspecific/total  $^{125}\text{I}$ -ChTX binding. Stock solutions of ChTX are prepared in 100 mM NaCl, 20 mM Tris-HCl, pH 7.4, 0.1% bovine serum albumin.

*Analysis of Data*—Data from saturation experiments are subjected to Scatchard analysis, and linear regression is performed to yield to equilibrium dissociation constant ( $K_d$ ) and maximum receptor concentration ( $B_{\max}$ ). Correlation coefficients for these determinations are typically greater than 0.95. Data from competition experiments are analyzed by the standard method of Cheng and Prusoff to determine  $K_i$  values. The dissociation rate constant ( $k_{-1}$ ) is determined directly from a first order plot of ligand dissociation *versus* time. The rate of ligand association ( $k_1$ ) is determined from the equation  $k_1 = k_{\text{obs}}([\text{LR}]_e/([\text{L}][\text{LR}]_{\max}))$ , where  $[\text{L}]$  is the concentration of ligand,  $[\text{LR}]_e$  is the concentration of the complex at equilibrium,  $[\text{LR}]_{\max}$  is the maximum number of receptors present, and  $k_{\text{obs}}$  is the slope of the pseudo-first order plot  $\ln([\text{LR}]_e/([\text{L}]_e - [\text{LR}]_t))$  *versus* time. Association and dissociation rate of  $^{125}\text{I}$ -ChTX are also determined by measuring the kinetics of radiolabeled toxin binding at different ligand concentrations, determining  $k_{\text{obs}}$  at each concentration of toxin from semilogarithmic representations of these data, and determining  $k_1$  and  $k_{-1}$  from the slope and y intercept, respectively, of the plot of  $k_{\text{obs}}$  *versus* ChTX concentration.

$^{125}\text{I}$ -ChTX Binding Assays may be used with the present invention substantially as described by Deutsch, C., *et al.*, J. of Biological Chemistry, 266: No. 6, 3668-3674 (1991).

25

#### EXAMPLE IV

##### Database search

30 A pore motif of voltage-gated potassium channels, SEQ ID NO:10: PASFWWATITMTTVGYGDIYP, was derived from hKv2.1, a Shab-related K channel

(Chandy and Gutman, Handbook of Receptors and Channels, *Ligand and Voltage-gated Ion Channels*, Edited by R. Alan North, CRC Press, Inc., Chapter 1 (1995)). SEQ ID NO:10 was used to query via *tblastn* search of an un-annotated proprietary database. *Basic Local Alignment Search Tool*, Altschul, S.F., *et al.*, J. Mol. Biol. 215:403 (1990). SEQ ID NO:8 (FIG.8) was identified in this manner. Seq ID NO: 8 was then used as a query sequence in a *blastn* search against a proprietary database to identify SEQ ID NO:9 (FIG.9) through overlapping clones.

Both SEQ ID NO:8 and SEQ ID NO:9 clones were isolated from a male and female donor cDNA library of adherent mononuclear cells. Library was constructed using 2 micrograms of polyA RNA isolated from adherent mononuclear cells, which came from a pool of male and female donors. The cells were cultured for 24 hours following Ficoll Hypaque centrifugation. CDNA synthesis was initiated using a NotI-oligo(dT) primer. Double-stranded cDNA was blunted, ligated to EcoRI adaptors, digested with NotI, size-selected, and cloned into the NotI and EcoRI sites of the pSPORT vector. The SEQ ID NO:8 and SEQ ID NO:9 clones were sequenced (ABI PRISM™ Dye Terminator Cycle sequencing on ABI PRISM™377 automated sequencer). SEQ ID NO:8 is 384 bp shorter than SEQ ID NO:9 and is identical in all but the 5'- end, wherein they differ by fifty five (55) nucleotides. This evidence suggests the existence of an alternate splice variant of the native gene. SEQ ID NO:9 was used to identify the full-length clone SEQ ID NO:1 described herein, by screening a human lymph-node cDNA library, as described below.

25

## EXAMPLE V

### cDNA library screening

30 A human lymph-node lambda gt10 cDNA library, oligo-dT plus random primed, from mRNA obtained from 6 donor male/female Caucasians, aged 26-29 whole lymph-nodes was used.

Clontech Laboratories, Palo Alto, CA (Cat # HL5000a; Lot # 46005). Standard methodology was followed as described in the Clontech Lambda Library protocol handbook (PT1010-1).

### Primary screen

- 5 About 30,000 plaque-forming units (pfu's) of lambda phage were plated onto 150 mm plates containing LB agar + 10 mM MgSO<sub>4</sub> + 200 ul of an overnight culture of the bacterial host strain E.coli C600 Hfi. Twenty such plates were made, a total of 600,000 plaques for primary screening.
- 10 Plates were inverted and incubated at 37°C overnight. The plaque diameters had just begun to make contact with one another. Plates were chilled for 1 hour.

- A 150 mm nylon membrane filter (Amersham) was marked in pencil and placed onto a LB agarose plate containing the plaques. The filter was further marked in 3 asymmetric
- 15 peripheral locations by puncturing the filter into the agar with a 16 gauge needle. The bottom of the plate corresponding to the hole was marked with a marker. After 1 minute, the filter was carefully peeled off and floated on top of DNA-denaturing solution (1.5 M NaCl, 0.5 M NaOH), plaque-side up, for 30 sec. The filter was removed and immersed in Neutralizing solution (1.5 M NaCl, 0.5 M Tris-HCl, pH 8) for 5 minutes, rinsed in 2x SSC (0.3 M NaCl,
- 20 0.03M Na-citrate) and placed on a Whatman 3M paper to dry. The DNA was fixed onto the filter by UV-crosslinking for 45 seconds. This procedure was repeated for all 20 plates. Second, duplicate filters were placed onto the original plates and similarly processed.

25

### Hybridization

The 40 filters prepared were pre-hybridized in 6X SSPE, 5X Denhart's and 0.25% SDS and 100 ug/ml of salmon-sperm DNA, for 4 hour at 42°C in a water bath shaker.

- 30 *Probe for library screening*

Primers were designed based on SEQ ID NO:9 to generate a probe via PCR. The primer sequences used were: 5' (forward) corresponding to SEQ ID NO:1 nucleotides 658-676, and 5' (reverse) corresponding to SEQ ID NO:1 nucleotides 1643-1661. The target sequence was amplified from 1 ng of SEQ ID NO:9 plasmid (in pSport1 vector, Life Technologies, Gaithersburg, MD) in a 50 µL reaction using Advantage™ KlenTaq polymerase (Clontech # 8417-1, lot 7020348) according to the manufacturer's recommendations. Cycling parameters were employed according to the Clontech Marathon-Ready kit conditions (touchdown PCR with three minute extensions). The target nucleic acid product was purified using Qiagen's Qiaquick PCR Purification kit, cat # 28104) according to the manufacturer's protocol and subsequently used to screen the human lymph node cDNA library. The probe was radiolabeled using <sup>32</sup>P-dCTP (Pharmacia Inc., Piscataway, NJ, cat# 27-9240 B).

The radiolabelled probe was added to the pre-hybridization mix containing the filters and allowed to hybridize at 42°C overnight. The filters were washed twice in buffer 2 (1X SSC, 0.5% SDS) at 65°C for 1 hour, placed on a Whatman filter paper and exposed to X-ray film overnight. Thirty eight positive plaques were identified. The plaques were picked and eluted in 1ml of dilution buffer using a cut P1000 pipet tip.

#### Secondary and tertiary screens

Of the 38 positive clones, 12 were [randomly] chosen for secondary screening. The clones were plated at dilutions of 1: 10,000, grown overnight, lifted on nylon filters, denatured, neutralized and probed with the 1 Kb SEQ ID NO:9 probe using a non-radioactive ECL method (Amersham, Inc., Cat# RPN 3001). Ten positive clones were picked (one from each plate) and processed for tertiary screening as described *supra*. At least 80% to 100% of the tertiary plaques were positive. One positive plaque was picked from each of the 10 plates.

#### Amplifying Lambda clone inserts

The 10 lambda inserts were amplified from 1 µl of phage supernatant in a 50 µl reaction using Advantage™ KlenTaq polymerase (Clontech # 8417-1, lot 7020348) according to the manufacturer's recommendations, with the addition of 5% final concentration DMSO. The primers were vector primers corresponding to λgt10 nucleotides 201-231 (forward) and 267-

298 (reverse), which surround the EcoRI cloning site (Appendix C, Clontech Lambda Library Protocol Handbook, manual # PT1010-1). Cycling parameters were according to Clontech Marathon-Ready cDNA User Manual (PT1156-1) p. 19, program 1 (briefly, touchdown PCR with three minute extensions). Product was purified using Qiagen's Qiaquick PCR  
5 Purification kit (cat # 28104) according to the manufacturer's protocol. Four PCR products were sequenced directly in the presence of DMSO (ABI PRISM™ Dye Terminator Cycle sequencing on ABI PRISM™377 automated sequencer). One clone was identified which contained the full-length cDNA: SEQ ID NO:1.

10

#### EXAMPLE VI

##### Northern blots

In order to analyze tissue-distribution and size of the transcripts corresponding to SEQ ID  
15 NO:1, Northern blot analyses were performed using two separate probes. Pre-made mRNA blots were used from Clontech which contained approximately 2ug poly A+ RNA per lane.

Probe 1 was generated by digesting SEQ ID NO:8 (in pSport1, Life Technologies) with PstI and ScaI restriction enzymes according to manufacturer's recommendations (Promega), and  
20 gel-isolating the 444 bp fragment using standard molecular biology techniques. This fragment corresponds to 3' UTR, SEQ ID NO:1 nucleotides 1716-2163, an area of exact match between SEQ ID NO:8 and SEQ ID NO:9. The restriction fragment was purified from agarose using Qiagen's Qiaquick Gel Extraction kit according to the manufacturer's recommendations and quantitated by 260 nm OD measurement. Twenty ng of fragment was labeled using  $\alpha^{32}\text{P}$ -  
25 dCTP (Dupont NEN) and Pharmacia's Ready-To-Go™ labeling beads (-dCTP) (cat # 27-9240-01) according to the manufacturer's recommendations. Labeled product was purified from unincorporated radioactive nucleotide using Pharmacia's ProbeQuant™ G-50 Micro Columns (cat # 27-5335-01) according to the manufacturer's recommendations.

30 Twenty ng ( $17 \times 10^6$  cpm) of labeled restriction fragment was used to probe Clontech Multiple Tissue Northern (catalogue # 7760-1 and 7759-1) and a prepared blot (described

*infra*) according to Clontech's recommendations (User Manual #PT1200-1). The blots were prehybridized for 45 minutes at 68°C in Clontech's ExpressHyb™ solution (cat # 8015-1). Denatured probe was added to 10 ml fresh pre-warmed ExpressHyb™. This replaced the prehybridization solution for 1.5 hours at 68°C. The blots were washed three times, ten minutes each, in 2X SSC + 0.1% SDS at room temperature, followed by two washes, twenty minutes each, in 0.1X SSC + 0.1% SDS at 50°C. Blots were exposed to film overnight using two intensifying screens.

Probe 2 was generated by PCR using the primers: 5' (forward) corresponding to SEQ ID NO:1 nucleotides 658-676, and 5' (reverse) corresponding to SEQ ID NO:1 nucleotides 1079-1098. The 440 bp product is unique to the SEQ ID NO:9 clone except for the reverse primer sequence, which is also found in the SEQ ID NO:8 clone. PCR reagents, template, and cycling conditions were the same as those used for generating cDNA library probe except that extension times were 1.5 minutes. Product was purified in the same manner as the library probe. Forty ng of Northern probe 2 was labeled exactly as Northern probe 1 and used to screen the Clontech Multiple Tissue Northern as described *supra* for Northern probe 1.

20

## 25 EXAMPLE VII

**Northern blot II**

Since calcium activated potassium channels are known to be dramatically up-regulated upon T-cell activation (e.g., Grissmer, *et al.*, J. Gen. Physiol. 102:601 (1993)), the expression pattern of SEQ ID NO:1 (SEQ ID NO:7) in resting versus activated human peripheral blood T-lymphocytes was determined.

30

Mononuclear cells were isolated from whole blood (obtained from healthy donors) on Ficoll-Hypaque density gradients as in Current Protocols in Immunology, vol.1, ed. by John E. Coligan, *et al.*, John Wiley and Sons, (1996). Red blood cells were removed by hypotonic lysis for 45 seconds in 0.2% saline followed by an equal volume of 1.6% saline to bring the saline back to physiological concentration. After an additional wash in Hank's Balanced Salt Solution (Life Technologies cat # 14175-095), CD14+ (monocyte and macrophage) and CD19+ (B cell) populations were removed using Dynal's Dynabeads® M450 Pan-B (CD19, cat # 111.03) and M450 CD14 (cat # 111.11) according to the standard Dynal protocol.

10 Remaining cells were T cells 80-90% as assessed by FACS analysis.

The T cells were cultured overnight at 37°C, 5% CO<sub>2</sub>, in complete RPMI 1640 medium (RPMI 1640, Life Technologies, cat # 11875-093, containing: 10% fetal bovine serum, Life Technologies, cat # 26140-087, heat inactivated 30 minutes at 56°C; 1X penicillin/streptomycin, Life Technologies, cat # 15140-122). The cells were then activated with 10 µg/mL final concentration phytohemagglutinin (PHA-P, Sigma cat # L9132) in complete RPMI for 48-72 hours at 37°C, 5% CO<sub>2</sub>. Messenger RNA was isolated using Invitrogen's FastTrack® 2.0 kit (cat # K1593-02). 2.8 µg of each mRNA was run on a 1% agarose/formaldehyde/formamide gel according to Sambrook *et al.*, Molecular Cloning Lab Manual, Second Edition, Cold Spring Harbor Press (1989), section 7:43, with 3µg of RNA ladder as size markers (Life Technologies, cat # 15620-016). RNA was transferred overnight to positively charged nylon membrane (Hybond™-N+, Amersham # RPN203B), according to the standard Amersham protocol, and then UV crosslinked to the membrane using Stratagene's Stratalinker on the automatic setting.

20

25

#### EXAMPLE VIII

##### PCR of SEQ ID NO:1 from various human cDNAs of hematopoietic origin

SEQ ID NO:1 (SEQ ID NO:7) message is determined to be present in other cells of hematopoietic origin by PCR with cDNA's isolated from different cell-lines of hematopoietic origin.

30



Target product was amplified from 2-2.5 ng of reverse transcribed mRNAs in a 20  $\mu$ L reaction using Advantage™ KlenTaq polymerase (Clontech # 8417-1, lot 7020348) according to the manufacturer's recommendations. Primer set 1 sequences are: 5' (forward) corresponding to  
5 SEQ ID NO:1 nucleotides 1462-1479, and 5' (reverse), corresponding to SEQ ID NO:1 nucleotides 1904-1922. Primer set 2 was the same set used to generate probe for the cDNA library screening. Primer set 2 reactions included DMSO at a final concentration of 5%. Cycling parameters were according to Clontech Marathon-Ready cDNA User Manual (PT1156-1) p. 19, program 1 (briefly, touchdown PCR with 1.5-2 minute extensions). 10  $\mu$ L  
10 of each reaction was analyzed on a 1% agarose gel containing 0.5  $\mu$ g/mL final concentration ethidium bromide in TAE buffer as in Sambrook *et al.*, Molecular Cloning Lab Manual, Second Edition, Cold Spring Harbor Press (1989), using 600 ng of 1kb DNA ladder as markers (Life Technologies, cat # 15615-016).

15

#### EXAMPLE IX

##### SEQ ID NO:7

*Library screening:* A  $^{32}$ P-dCTP labeled fragment from a clone (corresponding to SEQ ID  
20 NO:7 positions 658-1661) was used as a probe to screen ~ 600,000 recombinant plaques from a human lymph-node lambda gt10 cDNA library (Clontech Laboratories, Inc., Palo Alto, CA). Hybridizations were at 42°C overnight. Filters were washed twice in 1X SSC, 0.5% SDS at 65°C for 1 hour, and exposed to X-ray film overnight. Of 38 doubly-positive clones, 10 were subjected to two rounds of plaque purification and rescreening. Inserts were amplified using  
25 lambda-specific primers and amplicons were sequenced directly by automated sequencing as above. Six clones had SEQ ID NO:7 sequence information, and three were full-length. One full-length clone was subcloned and sequenced entirely in both directions, and used for subsequent expression constructs.

30 Computer analysis of the structural sequence was done using Lasergene software (DNASTar, Inc, Madison, WI). Alignments with other potassium channel sequences were performed

using the CLUSTAL algorithm (Megalign program, Lasergene) and these were utilized to create a dendrogram. The gap penalty and the gap length penalty were 10 each. Hydropathy plots were according to Kyte-Doolittle criteria, averaging over a 9-residue window (Protean program, Lasergene). Post-translational modification sites were identified using pattern  
5 searches within the Protean program (Lasergene). Patterns were derived from the Prosite database, and the threshold for matching was 100 percent.

#### EXAMPLE X

10

##### **Transient transfections**

A ~1.3 Kb SmaI/ScaI fragment containing the coding region (Nucleotides 390-1718 of **SEQ ID NO:7** and containing entire sequence for **SEQ ID NO:2**) was cloned into pcDNA3 vector (Invitrogen) at the EcoRV site. This cloning strategy introduced an additional methionine and  
15 two amino acids (G, A) upstream and in-frame with the authentic initiator methionine.

Approximately  $3 \times 10^5$  HEK 293 cells (ATCC, Rockville, MD) were transfected with 5  $\mu$ g of the novel gene in pcDNA3 vector along with 1  $\mu$ g of green fluorescent protein (GFP) in pEGFP-C1 vector (Clontech) using the LipoTAXI™ transfection kit (Stratagene, La Jolla, CA) as per manufacturer's instructions. Currents were recorded 24-72 hours later. GFP was  
20 used to track fluorescent cells for patch-clamping.

#### EXAMPLE XI

##### **Chromosomal Localization of the Novel Gene**

25 Two (2) ng of human genomic DNA (Sigma #D-4642 lot 41019) was used to amplify a 2 KB genomic DNA PCR fragment in a 20  $\mu$ L reaction using 1X Clontech's Advantage KlenTaq polymerase in 1X buffer (Clontech # 8417-1) and 200  $\mu$ M dNTPs (final concentration) plus 0.2  $\mu$ M each primer (final concentration). Forward primer corresponds to nucleotides 1462-1479 of **SEQ ID NO: 7**; Reverse primer corresponds to nucleotides 1904-1922 of **SEQ ID**  
30 **NO: 7**. The cycling parameters were as follows: one cycle of 94°C, 1 min. 30 sec.; five

cycles of 94°C, 30 sec., 72°C, 3 min.; five cycles of 94°C, 30 sec., 70°C, 3 min.; twenty five cycles of 94°C, 30 sec., 68°C, 3 min.

The 2 kb product obtained from the above PCR reaction was subcloned into pT7Blue  
5 (Novagen # 69820-1). Two clones were sequenced with vector-specific primers. PCR  
primers (Forward primer: 5' GTGGACATCTCCAAGGTACTAGGA 3' (SEQ ID NO:11) &  
Reverse primer: 5' TTTACTGACAGGCTGTGTGTGCCA 3' (SEQ ID NO:12)) were  
designed from the 2 Kb genomic sequence described *supra*, to reproducibly amplify a single  
139bp product from 25ng of Genomic DNA via the polymerase chain reaction (PCR) under  
10 the following conditions; 2mM MgCl<sub>2</sub>, 100µM dNTPs, & 1µM of each primer, for 35 cycles  
of 94°C for 30s, 65°C for 30s and 72°C for 30s. These primers and conditions were  
subsequently used to screen the "Research Genetics" BAC library in a 3D fashion by PCR  
(Bentley, *et al.*, 1992), using 10ng of pool DNA in a 50µl reaction volume. BAC 107 D 14  
was identified and streaked to give single colonies after overnight growth @ 37°C on LB-agar  
15 plates +chloramphenicol (12.5µg/ml). Colonies were screened by PCR (Huxley, *et al.*, 1990)  
to confirm the correct BAC had been isolated. One of the colonies was used as an inoculum  
for 500mls of LB+ chloramphenicol which was grown @37°C overnight with shaking @  
200rpm. This culture was used to prepare DNA by the alkaline-lysis method (Birnboim &  
Doly, 1979) and purified by banding on caesium chloride. The DNA was used as a probe for  
20 chromosomal-localization by fluroscent in-situ hybridization (Trask B, 1995), described  
below.

Chromosome spreads were prepared from PHA stimulated human lymphocytes which had  
been cultured in BrdU (200mg/ml) supplemented medium for 16 hours and in BrdU-free  
25 medium for a further 5 hours prior to colcemid treatment. Slides were stained for 30 minutes  
at room temperature in Hoechst 33258 (5mg/ml in 2xSSC), then placed face down in 2xSSC  
on a uv transilluminator for 30 minutes and dehydrated through an ethanol gradient. The  
chromosomes were denatured by incubation in 70% formamide, 2xSSC at 75°C for 5 minutes.

30 DNA from BAC 107d14 was biotinylated by nick translation using Bionick Kit according to  
the manufacturers instructions (BRL). 200ng labelled BAC and 10mg Cot-1 DNA (BRL) were

dissolved in 20ml hybridization buffer (50% formamide, 2xSSC, 10% dextran sulphate), denatured, and incubated at 37°C for 15 minutes before being hybridised to the denatured chromosomes at 37°C for 72 hours. Slides were washed for 3x5 minutes in 50% formamide/2xSSC, 2xSSC at 42°C and 0.1 xSSC at 60°C. Slides were blocked for one hour  
5 by incubation at room temperature in 4xSSC/ 5% dried milk/0.05% Triton X-100. Probe detection was with FITC-avidin (5mg/ml) with one round of signal amplification using biotinylated anti-avidin (5mg/ml) (Vector Laboratories) in 4xSSC/ 5% dried milk/ 0.05% Triton X-100. Washes between treatments were 3x 3minutes in 4xSSC/0.05% Triton X-100. Slides were counterstained with 0.1mg/ml DAPI and analysed on an Axioskop microscope  
10 (Zeiss) fitted with a Nu200 CCD camera (Photometrics) and SmartCapture software (Digital Scientific).

For the colocalization of BAC 107d14 and YAC 798E5, DNA (total yeast plus YAC) was labelled with digoxigenin dUTP by nick translation (Boehringer Mannheim kit). 300mg  
15 labelled DNA and 3mg Cot-1 DNA were added to the above mix and processed as above. Detection of the BAC was as described above but with rhodamine-avidin. YAC detection used 1mg/ml anti digoxigenin monoclonal antibody (Boehringer Mannheim), followed by 1/100 anti mouse IgG FITC conjugate, followed by 1/500 anti rabbit IgG FITC conjugate (Sigma). Signal was visible by eye. Paired signals were obtained on Chromosome 19 with  
20 BAC 107d14. Co-localization of BAC 107d14 and YAC 798E5 gave signals superimposed on each other indicating a localization to 19q13.1-13.2.

#### References for chromosomal localization

- Birnboim, H.C. & Doly, J. (1979). A rapid alkaline extraction procedure for screening  
25 recombinant plasmid DNA. *Nucleic Acids Research* 7 1513.
- Huxley, C., Green, E.D., & Dunham, I. (1990). Rapid assessment of *S.cerevisiae* mating type by PCR. *Trends Genet.* 6 236.
- Bentley, D.R., Todd, C., Collins, J., Holland, J., Dunham, I., Hassock, S., Bankier, A. & Gianelli, F. (1992). The development and application of automated gridding for efficient  
30 screening of yeast and bacterial ordered libraries. *Genomics* 12, 534 - 541.

Trask, B (1995) In genome analysis: A laboratory manual. Eds. B. Birren, E. Green and R. Myers, CSHL Press Inc.

## 5 EXAMPLE XII

### Patch Clamp Recording

Currents were recorded with an Axopatch 200A amplifier (Axon Instr., Foster City, CA) using the conventional whole cell, cell-attached, and inside-out configurations. Christian, E.P., *et al.*, J. Membr. Biol., 150:63 (1996). Transfected HEK cells were selected for recording by the presence of GFP labeling during epifluorescent illumination with a Nikon mercury arc light source and FITC filter block (excitation: 485/22 nm, emission: 505 nm). Thin-wall borosilicate glass pipettes were fabricated, sylgarded and fire-polished to a DC resistance of 2 - 8 MΩ. The resistance of patch seals was >10 GΩ. Liquid junction potentials (14) were corrected for in all experiments examining reversal potentials. Series resistance compensation of >70% was used in all experiments where maximal current was >0.5 nA. Voltage clamp protocols were implemented and data acquisition performed with pClamp 6.0 software (Axon Instr.). Currents were low-pass filtered at (-3 db at 1 or 2 kHz), and then digitized at 3-8 kHz as computer files with a TL-1 interface (Scientific Solutions, Solon, OH).

For seal formation, cells were bathed in normal Ringer solution containing (in mM) 160 NaCl, 4.5 KCl, 1 MgCl<sub>2</sub>, 5 HEPES and 2 CaCl<sub>2</sub> (except for excised inside-out patch formation, in which the bathing solution contained 0 CaCl<sub>2</sub>) at pH 7.4 (by NaOH) and osmolality ~325 mOsm (by sucrose). For *whole cell recording*, the pipette solution contained (in mM): 160 K aspartate, 2 MgCl<sub>2</sub>, 5 HEPES and 1.6 EGTA with either 0.8 CaCl<sub>2</sub> (calculated [Ca<sup>2+</sup>]<sub>free</sub> = 100 nM; Eqcal software, Biosoft Corp, Cambridge, UK), or 1.6 CaCl<sub>2</sub> (calculated [Ca<sup>2+</sup>]<sub>free</sub> = 1 μM) at pH 7.2 (by KOH) and osmolality ~315 mOsm (by sucrose). Cells were perfused locally with a solution containing 160 KCl, 2 CaCl<sub>2</sub>, 1 MgCl<sub>2</sub>, 5 HEPES, and 5 glucose at pH 7.4 (by KOH) and osmolality ~325 mOsm. In K<sup>+</sup> selectivity experiments equimolar Na<sup>+</sup> was substituted for K<sup>+</sup>. For *cell-attached patches*, the pipette solution contained 160 KCl, 2 CaCl<sub>2</sub>, 1 MgCl<sub>2</sub>, and 10 HEPES at pH 7.4 (by NaOH) and osmolality ~325 mOsm (for

determination of unitary conductance at physiological  $[K^+]_o$  the pipette contained Ringer).

Cells were perfused with Ringer solution +/- 1  $\mu$ M ionomycin. For *excised inside-out patches*, the pipette solution contained 160 KCl, 2  $CaCl_2$ , 1  $MgCl_2$  and 10 HEPES at pH 7.4 (by KOH) and osmolality ~325 mOsm. The local perfusion solution (cytoplasmic face)

- 5 contained 160 K aspartate, 2  $MgCl_2$ , 5 HEPES and 10 EGTA at pH 7.2 (by KOH) and osmolality ~329 mOsm with  $CaCl_2$  added at concentrations that yielded the calculated  $[Ca^{2+}]_{free}$  specified. Perfusion solutions were delivered locally using a solenoid valve system (13). Toxins were purchased from Sigma Chemicals (St. Louis, MO) and were diluted in solution containing 0.01% bovine serum albumin.

## EXAMPLE XIII

5

**CaPO<sub>4</sub> transfection protocol for HEK293 cells****REAGENTS**

10

Prepare 1) 2 x HBS

Reagent	Final	grams, dry
---------	-------	------------

HEPES pH7.05	50mM	5.0
--------------	------	-----

KCl	10mM	0.37
-----	------	------

15 dextrose	12mM	1.0
-------------	------	-----

NaCl	280mM	8.0
------	-------	-----

Na <sub>2</sub> HPO <sub>4</sub> (MW 141.96)	1.5mM	0.1065g
--	-------	---------

500ml total, pH to 7.05 (**NOTE pH is important**).

20

2) CaCl<sub>2</sub> 2M 29.4g100ml total

Stock solutions may be frozen at -20°C. Once opened stocks may be stored at 4°C for 6

25 weeks.

30

**METHOD** (for each T75 plate):

1. Plate cells 2-3 days prior to transfection to ensure ~75% confluence immediately prior to transfection. This is really important, it is vital that there are still some gaps on the plate before you start.
2. On day of transfection, prepare DNA in water to a final vol. of 876 µl in 15ml (polystyrene) Falcon tube. 20 mg of hKCa4 plasmid was used (pcDNA3 or pGEN IRES-neo) and 10 µg of pEGFP-C1.
3. Add 124 µl of 2M CaCl<sub>2</sub> to DNA and mix.
4. Remove media from cells very carefully.
5. Add 8 ml of fresh media very slowly down the side of the plate to ensure cells *are* not dislodged.
6. Add chloroquine to final concentration of 25 mM (0.1 ml of 100x stock (2.5mM)) and swirl gently.
7. Add 1 ml of 2X HBS onto DNA dropwise while bubbling (or mix the solution) for 30 secs.
8. Immediately, add mixture onto cells dropwise covering entire plate, and leave at 37°C for 6-10 hours (usually 6)
9. **NOTE: Leave chloroquine on cells for a MAX of 10 hours.**
10. At the end of the day, replace with fresh media and again the following morning. The cells should be ready to check for GFP at that point, and up to 2 days after the transfection.

For stables: On day 4, (48 hours after the last media change), split cells into 150 mm dishes (or T150s for episomal, if you're not going to clone individual colonies) at ~1:10 in complete medium.

On day 5, (24 hours after splitting), put cells into selective media; HEK293 take 1 mg/mL G418 and 250 µg/mL hygromycin B.

30

EXAMPLE XIV



**Transfection of cho-k1 cells with cellfectin**

1. Split cells 2-3 days ahead of time so that they are 50-70% confluent in a T75 flask ( e.g., grow CHO-K1 cells in Ham's F12 plus 10% FBS, no antibiotics).
- 5 2. Place 30 µg hKCa4 DNA (*either* pGEN IRES-neo *or* pcDNA3) plus 15 µg pEGFP-C1 in a 15 mL Falcon tube (polystyrene).
3. Add 2 mL of media without serum or antibiotics (in this case Ham's F12 (called SFM (for short))).
4. In a separate polystyrene tube, mix 2 mL of SFM with Cellfectin reagent (50 µL per T75  
10 (LTI# 10362-010)).\*
5. Combine the contents of the two tubes, mix, and let sit at room temperature 15 minutes.
6. Wash the cells in SFM, then add 5 mL SFM per T75.
7. Add 11 mL SFM to the DNA/Cellfectin mix (bringing the volume to 15 mL).
8. Immediately overlay the DNA/Cellfectin mix onto the cells.
- 15 9. Incubate at 37°C in 5% CO<sub>2</sub> for 5 hours.
10. Remove DNA/Cellfectin mix and replace with Ham's F12 plus 10% FBS.
11. Next day change media and check for EGFP.

For stables: On day 4, (48 hours after the last media change), split cells into 150 mm dishes  
20 (or T150s for episomal, if clone individual colonies are not required) at ~1:10 in complete medium.

On day 5, (24 hours after splitting), put cells into selective media; CHO-K1 take 1 mg/mL G418 and 300 µg/mL hygromycin B.

25

**EXAMPLE XV****Baculovirus expression system for hKCa4 (SEQ ID NO:2)**

A recombinant virus bearing the coding region of hKCa4 was constructed using the Bac-To-  
30 Bac Baculovirus Expression System (Life Technologies, Gaithersburg, MD, cat. # 10359-016). The coding region of hKCa4 was subcloned into the donor vector pFastBac1 by

- 76 -

cleaving the pcDNA3 construct (described in JBC 272 (52): 32723-32726) with EcoRI (5') and XhoI (3') sites, releasing the 1356 bp insert, and cloned into these sites within pFastBac1. This strategy retains the additional methionine and two amino acids (G, A) upstream and in-frame with the authentic initiator methionine. Transposition was as detailed in the protocol  
5 from Life Technologies, Gaithersburg, MD, with the following notes: transformation plates were incubated 48 hours at 37°C for transposition to occur. Restreaked colony plates were also incubated 48 hours at 37°C to ensure white phenotype. For transfecting Sf9 cells, 8 µL of bacmid miniprep DNA and 5 µL Cellfectin was used. The DNA plus Cellfectin mixes were allowed to incubate 20 minutes before applying to cells.

10

Passage two stocks of the virus were obtained by infecting  $3 \times 10^6$  Sf9 cells in a T25 flask with 100 µL of the passage one virus (the transfection supernatant produced above) in a total volume of 5 mL. The infection was allowed to proceed for three days at 27°C in a humidified chamber, and the cell culture supernatants were collected. This passage two virus was  
15 quantitated using the BacPAK Baculovirus Rapid Titer Kit (Clontech, Palo Alto, CA, cat # K1599-1) or by conventional plaque assay (Life Technologies, Gaithersburg, MD, Baculovirus Expression Training course manual, 3.54-3.58).

High-titer passage three viral stocks were generated from passage two virus by infecting Sf9 cells in Sf900II medium in suspension at a MOI of 0.1 and a cell density of  $4.0 \times 10^6$ /mL.  
20 After allowing virus to adhere to cells for one hour at room temperature with constant agitation, cells were harvested and resuspended in fresh medium to a cell density of  $1.0 \times 10^6$ /mL. The infection was allowed to continue for two days in suspension at 27°C and 110 rpm. The supernatant (passage three virus) was harvested and quantitated as for the passage 2 virus. Infections for protein production were performed at various MOIs and times of harvest  
25 as per the Bac-To-Bac manual. Cells were examined by patch clamp analysis and/or  $^{125}$ I-charybdotoxin (CTX) binding.

30

## EXAMPLE XVI

**Scintillation proximity assay (SPA)**

This method has been successfully applied to receptor binding assays by immobilizing receptors directly on to SPA beads via a number of coupling methods (Bosworth N and Towers P, Nature, 342: 167-168, 1989; Undenfriend S et. al., Anal. Biochem. 161: 494-500, 1987). The most commonly used beads are those coupled to Wheat Germ Agglutinin (WGA), which bind to N-acetyl b-D glucosamine residues in membrane glycolipids and glycoproteins. (Nagata Y and Burger MM, JBC, 249: 3116-3122, 1974). Once immobilized, the receptor is close enough to the bead so that, should a suitably radiolabelled ligand bind to the receptor, it will be in close enough proximity to stimulate the bead to emit light. Any unbound radioligand is too distant from the bead to transfer energy and goes undetected. The bead therefore only detects the population of ligand molecules which are receptor-bound. The discrimination of binding by proximity means that no separation of bound and free ligand is required, as in traditional methods. SPA is an extremely powerful technique when large numbers of samples are required to be assayed simultaneously, as during a high-throughput assay for compounds (News, 27: November 1996, Amersham Life Sciences).

All reagents for the SPA assay (except buffer solutions), were ordered from Amersham Life Sciences, Arlington Heights, IL.

Assay Buffer:

5mM NaCl  
5mM KCl  
10mM HEPES  
6mM glucose  
bring to pH 8.4 with TRIS base

Wash Buffer:

200mM NaCl  
20mM HEPES  
pH 8.0 with Tris base

20

assay run in 96-well SPA microtiter plates

0.5-3 mg per well of wheat germ agglutinin-coated SPA beads (Amersham, Cat# RPNQ0001).

5-30 ug of KCa4-HEK cell membranes per well

- 78 -

50 pM of <sup>125</sup>I charybdotoxin (CTX, NEN, Boston, MA)

non specific defined with 50 nM cold charybdotoxin (Bachem, King of Prussia, PA).

**ASSAY:**

5 50 ul beads + membranes

50 ul assay buffer

50 ul experimental compound or nonspecific or buffer (defines total binding)

50 ul ligand

10 Membranes are allowed to pre-incubate with SPA beads for 1-2 hours; unbound membranes are washed away from beads with centrifugation.

Assay buffer is added to each well, followed by compounds, bead-membranes and ligand

15 The assay incubates with shaking for 1 hour at room temperature

The beads are allowed to settle overnight, then the plate is counted in a Top Count scintillation counter (Packard, Meriden, CT).

20

**EXAMPLE XVII**

**Anergic T-cell line**

3926 alloT is one of a number of HLA-DR4Dw4-specific alloreactive T cell lines which were

25 produced according to standard methodology. Twenty-five million peripheral blood mononuclear cells from a donor not carrying the HLA-DR4 allotype were mixed with the same number of irradiated blood mononuclear cells from a donor carrying the HLA-DR4Dw4 allotype in 5ml of RPMI 1640 cell growth medium supplemented with glutamine and 5% human serum. After 14 days culture, viable cells were isolated by centrifugation over a  
30 Lymphopaque density gradient, and 5 million cells added to a further twenty-five million irradiated blood mononuclear cells from the same HLA-DR4Dw4-positive donor.

Recombinant hIL-2 (2 Units/ml) was added on day 4, and every 3 days afterwards. The line was maintained by periodic restimulation and IL-2 addition as described. The HLA-DR4Dw4 allospecificity of the lines was confirmed by demonstrating that they were able to proliferate in response to HLA-DR4Dw4 homozygous EBV B cells (JAH cell line) and HLA-DR4Dw4-  
5 expressing L cell transfectants.

#### **maintenance**

A frozen vial of the 3926 T cell line from UK was thawed and maintained in RPMI medium  
10 with 10% FCS and glutamine. The cells were maintained on a "rest-stimulation" protocol by stimulating them with soluble anti-CD3 (2.5 mg/ml catalog # 30100D plus 2.5 mg/ml catalog # 30110D from PharMingen, San Diego, CA) plus mitomycin-c-treated untyped human PBMC to provide the costimulatory signal. The PBMC would also act to crosslink the soluble CD3 on the TCR. After two days, 5 Units/ml of human recombinant IL-2 (Boehringer  
15 Mannheim, Indianapolis, IN) was added to the cells. During this stimulation step, the T cells were seen to clump and proliferate. T cells used in experiments were rested at least one week after addition of IL-2 and were separated from debris and dead PBMC by Histopaque 1077 (Sigma, St. Louis, MO) density gradient centrifugation.

#### **20 pretreatment conditions for resting, activated and anergic cells**

Purified 3926 allo T cells ( $4 \times 10^5$  cells/well) were incubated in 24-well plates which were uncoated ("resting" cells) or coated with anti-CD3 plus anti-CD28 (2.5 mg/ml each of the two anti-CD3 antibodies plus 20 mg/ml catalog # 33740D- "two signal"); or in order to induce anergy, with 2.5 mg/ml each of the two human anti-CD3 antibodies in a 6 well plate ("one  
25 signal"). All antibodies were diluted in PBS. The cells were washed away from the antibodies after one day and rested in fresh wells for two more days (total of three days following exposure to the antibodies) at which time the T cells were patch-clamped for KCa channels as described elsewhere in the document.

30

#### **EXAMPLE XVIII**

### Cloning of murine KCa4

An effort to clone a rodent KCa4 gene was initiated because there are several rodent models of autoimmunity available, and it was important to verify that the rodent genes are essentially the same as the human gene, for these animal models to be used. An earlier publication had already shown that the KCa channel in rodent lymphocytes behaves very similar to its human counterpart by electrophysiology; it is low in resting cells and is dramatically upregulated upon lymphocyte activation (Mahaut-Smith and Mason J. Physiol. 1991).

- 10 The amino acid sequence of human KCa4 (SEQ ID NO:3) was used as a query sequence in a TBLASTN search against the embl est database. One mouse EST, accession number W30402, showed 79% identity to hKCa4 at the amino terminus (N-term.) (over an 89 amino acid region) and 63% identity to hKCa4 at the carboxy terminus (C-term.) (over only a 20 amino acid region). The conceptual translation of the EST resulted in a truncated protein compared to that expected for a 6 transmembrane ion channel subunit, wherein the portion of the molecule corresponding to the third transmembrane region through all but the last 20 amino acids are deleted. However, because of the high homology to the human sequence, primers were designed upstream of the putative start methionine and downstream of the putative stop codon to PCR "full-length" mouse KCa4. The forward primer sequence is
- 20 corresponds to W30402 nucleotides 20-42. The reverse primer sequence corresponds to W30402 nucleotides 413-436.

- cDNA was made from the mouse erythroleukemic cell line MEL-C88 (Deisseroth, A., *et al.*, Cell 15: 55 (1978)). Cells were cultured in  $\alpha$ MEM medium containing 10% fetal bovine serum, 100 units per mL of penicillin and 100  $\mu$ g per mL of streptomycin. Total RNA was isolated using the RNeasy Midi prep Kit from Qiagen, Inc., Chatsworth, CA, cat. # 75142 and treated with DNase I, Life Technologies, Gaithersburg, MD, cat. # 18068-015, according to manufacturer's instructions. cDNA synthesis was with SuperScript II reverse transcriptase (Life Technologies, Gaithersburg, MD, cat. # 18064-014) using the Marathon cDNA synthesis primer (Marathon cDNA Amplification kit, Clontech, Palo Alto, CA, cat. # K1802-1),
- 25
- 30

according to the SuperScript II protocol. First strand cDNA was digested with RNase H (Life Technologies, Gaithersburg, MD, cat. # YD1220) and used at ~1 ng per PCR.

mKCa4 products were amplified from the cDNA in 25  $\mu$ L reactions using Advantage™  
5 KlenTaq polymerase (Clontech, Palo Alto, CA # 8417-1) according to the manufacturer's recommendations, with the addition of 5% final concentration DMSO. Cycling parameters were as follows: one cycle at 94°C for 1.5 minutes; thirty cycles of: 94°C thirty seconds, 65°C forty five seconds, 68°C one minute. Ten  $\mu$ L of each reaction was analyzed on a 1% agarose gel containing 0.5  $\mu$ g/mL final concentration ethidium bromide in TAE buffer as in  
10 Sambrook et al, Molecular Cloning Lab Manual, Second Edition, Cold Spring Harbor Press 1989, using 600 ng of 1kb DNA ladder as marker (Life Technologies, Gaithersburg, MD, cat # 15615). The remaining 15  $\mu$ L of the reactions were purified using the Qiaquick PCR purification kit according to the manufacturer's recommendations (Qiagen, Inc., Chatsworth, CA, cat # 28104), quantitated by Optical Density measurements at 260 nm, and subjected to  
15 automated sequencing (ABI PRISM Dye Terminator Cycle sequencing on ABI PRISM 377 automated sequencer).

The full-length coding region nucleotide sequence, translated amino acid sequence and amino acid alignment of the mouse ortholog and human KCa4 sequences are shown in FIG.22 (SEQ  
20 ID NO:13), FIG.23 (SEQ ID NO:14) (murine homolog (ortholog of human Kca4, SEQ ID NO:3)), and FIG.24, respectively.

#### EXAMPLE XIX

25

##### **Isolation of microglia & transcripts**

Human fetal brain tissue (from 16- to 22-week-old fetuses) was obtained commercially from the Anatomic Gift Foundation (Woodbine, GA). Cerebral hemispheres were placed in chilled, sterile, calcium-free, magnesium-free Hanks Balanced Salt Solution (HBSS, # 14170-112,  
30 Gibco-BRL, Gaithersburg, MD). Meninges were removed with sterile forceps and tissue was dissociated initially by repeated trituration through sterile pipets. Tissue was incubated with

0.05% Trypsin/0.53 mM EDTA (# 25300-054, Gibco-BRL, Gaithersburg, MD) and 0.15 mg/ml DNase (# D-5025, Sigma Chemical Co., St. Louis, MO) at 37°C for 45 minutes with gentle shaking. 10% Fetal Bovine Serum (FBS, # 160001-036, Gibco-BRL, Gaithersburg, MD) was added to the suspension to stop trypsinization. The cells were then passed through  
5 210 µm and 149 µm polypropylene mesh (# 08-670-188 and # 08-670-189, Fisher Scientific, Pittsburgh, PA) sequentially. The filtrate was washed twice and resuspended in complete media [DMEM (with high glucose, L-glutamine and HEPES), Gibco # 12430-054; 100 U-µg/ml penicillin/streptomycin, Gibco # 15070-030; 10% FBS, Gibco # 160001-036]. Cells were plated at density of 80 million/75 cm<sup>2</sup> flask (# 12-565-52, Fisher Scientific, Pittsburgh,  
10 PA). Cultures were incubated at 37°C, 5% CO<sub>2</sub> for 2 weeks (with one media change after one week in culture to remove cellular debris). Mature cultures consisted of astrocytes, neurons, microglia and oligodendrocytes. The microglial population is enriched in the supernatant, which was harvested (astrocytes and oligodendrocytes stick to the plate). Cells were lysed in Buffer RLT (QIAGEN) with 0.01% b-ME (Sigma Cat.# M7154) and frozen at -80°C until  
15 use. Lysed cells were thawed on ice, and processed through a QIAshredder (QIAGEN Cat.# 79654 ) to homogenize.

Total RNA was then prepared with the RNeasy Mini Kit (QIAGEN Cat.#74103) according to manufacturer's directions. Concentration was determined by absorbance at 260 nm. CDNA  
20 was prepared from 3.2 ug of microglial RNA by standard reverse transcription using SuperScript II reverse transcriptase kit (Life Technologies, Gaithersburg, MD, cat. # 18064-014). PCR was performed at 65 degrees annealing. Human KCa4 specific primers used were: Forward primer: positions 862-883 of SEQ ID NO:1; and reverse primer corresponding to positions 1659-1680 of SEQ ID NO:1.

25

\* \* \*

All publications and patents mentioned in the above specification are herein incorporated by  
30 reference. Various modifications and variations of the described methods and system of the invention will be apparent to those skilled in the art without departing from the scope and



spirit of the invention. Although the invention has been described in connection with specific preferred embodiments, it should be understood that the invention as claimed should not be unduly limited to such specific embodiments. Indeed, various modifications of the described modes for carrying out the invention which are obvious to those skilled in molecular biology  
5 or related fields are intended to be within the scope of the following claims.

- 84 -

## SEQUENCE LISTING

## (I) GENERAL INFORMATION

## (i) APPLICANT:

(ii) TITLE OF THE INVENTION: HUMAN LYMPH NODE DERIVED  
CALCIUM ACTIVATED POTASSIUM CHANNEL

(iii) NUMBER OF SEQUENCES: 14

## (iv) CORRESPONDENCE ADDRESS:

(A) ADDRESSEE: ZENECA Pharmaceuticals, Inc.  
(B) STREET: 1800 Concord Pike  
(C) CITY: Wilmington  
(D) STATE: DE  
(E) COUNTRY: USA  
(F) ZIP: 19850-5437

## (v) COMPUTER READABLE FORM:

(A) MEDIUM TYPE: Diskette  
(B) COMPUTER: IBM Compatible  
(C) OPERATING SYSTEM: DOS  
(D) SOFTWARE: FastSEQ for Windows Version 2.0

## (vi) CURRENT APPLICATION DATA:

(A) APPLICATION NUMBER:  
(B) FILING DATE:  
(C) CLASSIFICATION:

## (vii) PRIOR APPLICATION DATA:

(A) APPLICATION NUMBER: (UK)9714760.7  
(B) FILING DATE: 15-JUL-1997

(A) APPLICATION NUMBER: (UK)9721366.4  
(B) FILING DATE: 09-OCT-1997

## (viii) ATTORNEY/AGENT INFORMATION:

(A) NAME: Higgins, Patrick H  
(B) REGISTRATION NUMBER: 39,709

- 85 -

(C) REFERENCE/DOCKET NUMBER: PHM.70235

(ix) TELECOMMUNICATION INFORMATION:

(A) TELEPHONE: 302.886.4889

(B) TELEFAX: 302.886.8221

(C) TELEX:

(2) INFORMATION FOR SEQ ID NO:1:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 2261 base pairs

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:1:

```
GTCTTCGGT GTCTGGGTGT GGTGAGTAGA GGTGTGTGTC ACAAAGTACA GACCATTGTG 60
TGTGACAAAG CCCATCGTGT GTCTGTGTGT GTCTTTATCC ACGTGGATGG ACGTCTCTTT 120
CTTGCTCTGC CCCAAGACAC ACCCTAGCCC CTCCTTATTC TCAAAAGGGG GAGCTGGGGA 180
GCCTCCCCCT ACCCTGGGGC CTCCCCTGCC CCTCCCCGCC CTGCCTGGCC GTCACCACTC 240
CCCAGAGGGC ACAGGGCTCT GCTGTGCCTC AGAGCAAAAG TCCCAGAGCC AGCAGAGCAG 300
GTGACGACCT GCAAGCCACA GTGGCTGCCC TGTGCGTGCT GCGAGGTGGG GGACCCTGGG 360
CAGGAAGCTG GCTGAGCCCC AAGACCCCGG GGGCCATGGG CGGGGATCTG GTGCTTGGCC 420
TGGGGGCCCT GAGACGCCGA AAGCGCTTGC TGGAGCAGGA GAAGTCTCTG GCCGGCTGGG 480
CACTGGTGCT GGCAGGAAct GGCATTGGAC TCATGGTGCT GCATGCAGAG ATGCTGTGGT 540
TCGGGGGGTG CTCGTGGGCG CTCTACCTGT TCCTGGTTAA ATGCACGATC AGCATTTCCT 600
CCTTCTTACT CCTCTGCCTC ATCGTGGCCT TTCATGCCAA AGAGGTCCAG CTGTTCATGA 660
CCGACAACGG GCTGCGGGAC TGGCGCGTGG CGCTGACCGG GCGGCAGGCG GCGCAGATCG 720
TGCTGGAGCT GGTGGTGTGT GGGCTGCACC CGGCGCCCGT GCGGGGCCCC CCGTGCGTGC 780
AGGATTTAGG GCGCGCGCTG ACCTCCCCGC AGCCCTGGCC GGGATTCTTG GGCCAAGGGG 840
AAGCGCTGCT GTCCCTGGCC ATGCTGCTGC GTCTCTACCT GGTGCCCCGC GCCGTGCTCC 900
TGCGCAGCGG CGTCCTGCTC AACGCTTCCT ACCGCAGCAT CGGCGCTCTC AATCAAGTCC 960
GCTTCCGCCA CTGGTTCGTG GCCAAGCTTT ACATGAACAC GCACCCTGGC CGCTGCTGC 1020
TCGGCCTCAC GCTTGGCCTC TGGCTGACCA CCGCCTGGGT GCTGTCCGTG GCCGAGAGGC 1080
AGGCTGTAA TGCCACTGGG CACCTTTCAG ACACACTTTG GCTGATCCCC ATCACATTCC 1140
TGACCATCGG CTATGGTGAC GTGGTGCCGG GCACCATGTG GGGCAAGATC GTCTGCCTGT 1200
GCACTGGAGT CATGGGTGTC TGCTGCACAG CCCTGCTGGT GGCCGTGGTG GCCCGGAAGC 1260
TGGAGTTTAA CAAGGCAGAG AAGCACGTGC ACAACTTCAT GATGGATATC CAGTATACCA 1320
AAGAGATGAA GGAGTCCGCT GCCCGAGTGC TACAAGAAGC CTGGATGTTC TACAAACATA 1380
```

- 86 -

CTCGCAGGAA GGAGTCTCAT GCTGCCCCGA GGCATCAGCG CAAGCTGCTG GCCGCCATCA 1440  
 ACGCGTTCGG CCAGGTGCGG CTGAAACACC GGAAGCTCCG GGAACAAGTG AACTCCATGG 1500  
 TGGACATCTC CAAGATGCAC ATGATCCTGT ATGACCTGCA GCAGAATCTG AGCAGCTCAC 1560  
 ACCGGGGCCCT GGAGAAACAG ATTGACACGC TGGCGGGGAA GCTGGATGCC CTGACTGAGC 1620  
 TGCTTAGCAC TGCCCTGGGG CCGAGGCAGC TTCCAGAACC CAGCCAGCAG TCCAAGTAGC 1680  
 TGGACCCACG AGGAGGAACC AGGCTACTTT CCCCAGTACT GAGGTGGTGG ACATCGTCTC 1740  
 TGCCACTCCT GACCCAGCCC TGAACAAAGC ACCTCAAGTG CAAGGACCAA AGGGGGCCCT 1800  
 GGCTTGGAGT GGGTTGGCTT GCTGATGGCT GCTGGAGGGG ACGCTGGCTA AAGTGGGTAG 1860  
 GCCTTGGCCC ACCTGAGGCC CCAGGTGGGA ACATGGTCAC CCCCCTCTG CATAACCTCA 1920  
 TCAAAAACAC TCTCACTATG CTGCTATGGA CGACCTCCAG CTCTCAGTTA CAAGTGCAGG 1980  
 CGACTGGAGG CAGGACTCCT GGGTCCCTGG GAAAGAGGGT AATAGGGGCC CGGATXCAGG 2040  
 ATTCTGGGAG GCTTCAGTTA CCGCTGGCCG AGCTGAAGAA CTGGGTATGA GGCTGGGGCG 2100  
 GGGCTGGAGG TGGCGCCCCC TGGTGGGACA ACAAAGAGGA CACCATTTT CCAGAGCTGC 2160  
 AGAGAGCAXC TGGTGGGGAG GAAGAAGTGT AACTCACCAG CCTTTGCTCT TATTTTGTGA 2220  
 ATAAATGTGA AAGCCAGAAA AAAAAAAAAA AAAAAAAAAA A 2261

## (2) INFORMATION FOR SEQ ID NO:2:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 1284 base pairs  
 (B) TYPE: nucleic acid  
 (C) STRANDEDNESS: single  
 (D) TOPOLOGY: linear

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO:2:

ATGGGCGGGG ATCTGGTGCT TGGCCTGGGG GCCTTGAGAC GCCGAAAGCG CTTGCTGGAG 60  
 CAGGAGAAGT CTCTGGCCGG CTGGGCACTG GTGCTGGCAG GAACTGGCAT TGGACTCATG 120  
 GTGCTGCATG CAGAGATGCT GTGGTTCGGG GGGTGCTCGT GGGCGCTCTA CCTGTTCCTG 180  
 GTTAAATGCA CGATCAGCAT TTCCACCTTC TTA CTCTCT GCCTCATCGT GGCCTTTCAT 240  
 GCCAAAGAGG TCCAGCTGTT CATGACCGAC AACGGGCTGC GGGACTGGCG CGTGGCGCTG 300  
 ACCGGGCGGC AGGCGGCGCA GATCGTGCTG GAGCTGGTGG TGTGTGGGCT GCACCCGGCG 360  
 CCCGTGCGGG GCCCGCCGTG CGTGCAGGAT TAGGGGCGC CGCTGACCTC CCCGCAGCCC 420  
 TGGCCGGGAT TCCTGGGCCA AGGGGAAGCG CTGCTGTCCC TGGCCATGCT GCTGCGTCTC 480  
 TACCTGGTGC CCCGCGCCGT GTCCTGCGC AGCGGCGTCC TGCTCAACGC TTCCTACCGC 540  
 AGCATCGGCG CTCTCAATCA AGTCCGCTTC CGCCACTGGT TCGTGGCCAA GCTTTACATG 600  
 AACACGCACC CTGGCCGCTT GCTGCTCGGC CTCACGCTTG GCCTCTGGCT GACCACCGCC 660  
 TGGGTGCTGT CCGTGGCCGA GAGGCAGGCT GTTAATGCCA CTGGGCACCT TTCAGACACA 720  
 CTTTGGCTGA TCCCCATCAC ATTCCTGACC ATCGGCTATG GTGACGTGGT GCCGGGCACC 780  
 ATGTGGGGCA AGATCGTCTG CCTGTGCACT GGAGTCATGG GTGTCTGCTG CACAGCCCTG 840

- 87 -

CTGGTGGCCG TGGTGGCCCG GAAGCTGGAG TTTAACAAGG CAGAGAAGCA CGTGCACAAC 900  
 TTCATGATGG ATATCCAGTA TACCAAAGAG ATGAAGGAGT CCGCTGCCCC AGTGCTACAA 960  
 GAAGCCTGGA TGTTCTACAA ACATACTCGC AGGAAGGAGT CTCATGCTGC CCGCAGGCAT 1020  
 CAGCGCAAGC TGCTGGCCGC CATCAACGCG TTCCGCCAGG TCGGGCTGAA ACACCGGAAG 1080  
 CTCCGGGAAC AAGTGAATC CATGGTGGAC ATCTCCAAGA TGCACATGAT CCTGTATGAC 1140  
 CTGCAGCAGA ATCTGAGCAG CTCACACCGG GCCCTGGAGA AACAGATTGA CACGCTGGCG 1200  
 GGAAGCTGG ATGCCCTGAC TGAGCTGCTT AGCACTGCCC TGGGGCCGAG GCAGCTTCCA 1260  
 GAACCCAGCC AGCAGTCCAA GTAG

1284

## (2) INFORMATION FOR SEQ ID NO:3:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 427 amino acids  
 (B) TYPE: amino acid  
 (C) STRANDEDNESS: single  
 (D) TOPOLOGY: linear

## (ii) MOLECULE TYPE: peptide

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO:3:

Met	Gly	Gly	Asp	Leu	Val	Leu	Gly	Leu	Gly	Ala	Leu	Arg	Arg	Arg	Lys
1				5					10					15	
Arg	Leu	Leu	Glu	Gln	Glu	Lys	Ser	Leu	Ala	Gly	Trp	Ala	Leu	Val	Leu
			20					25					30		
Ala	Gly	Thr	Gly	Ile	Gly	Leu	Met	Val	Leu	His	Ala	Glu	Met	Leu	Trp
			35				40					45			
Phe	Gly	Gly	Cys	Ser	Trp	Ala	Leu	Tyr	Leu	Phe	Leu	Val	Lys	Cys	Thr
			50				55				60				
Ile	Ser	Ile	Ser	Thr	Phe	Leu	Leu	Leu	Cys	Leu	Ile	Val	Ala	Phe	His
65				70					75					80	
Ala	Lys	Glu	Val	Gln	Leu	Phe	Met	Thr	Asp	Asn	Gly	Leu	Arg	Asp	Trp
				85					90					95	
Arg	Val	Ala	Leu	Thr	Gly	Arg	Gln	Ala	Ala	Gln	Ile	Val	Leu	Glu	Leu
			100					105					110		
Val	Val	Cys	Gly	Leu	His	Pro	Ala	Pro	Val	Arg	Gly	Pro	Pro	Cys	Val
			115					120					125		
Gln	Asp	Leu	Gly	Ala	Pro	Leu	Thr	Ser	Pro	Gln	Pro	Trp	Pro	Gly	Phe

- 88 -

130	135	140
Leu Gly Gln Gly Glu Ala Leu Leu Ser Leu Ala Met Leu Leu Arg Leu		
145	150	155
Tyr Leu Val Pro Arg Ala Val Leu Leu Arg Ser Gly Val Leu Leu Asn		160
	165	170
Ala Ser Tyr Arg Ser Ile Gly Ala Leu Asn Gln Val Arg Phe Arg His		175
	180	185
Trp Phe Val Ala Lys Leu Tyr Met Asn Thr His Pro Gly Arg Leu Leu		190
	195	200
Leu Gly Leu Thr Leu Gly Leu Trp Leu Thr Thr Ala Trp Val Leu Ser		205
	210	215
Val Ala Glu Arg Gln Ala Val Asn Ala Thr Gly His Leu Ser Asp Thr		220
225	230	235
Leu Trp Leu Ile Pro Ile Thr Phe Leu Thr Ile Gly Tyr Gly Asp Val		240
	245	250
Val Pro Gly Thr Met Trp Gly Lys Ile Val Cys Leu Cys Thr Gly Val		255
	260	265
Met Gly Val Cys Cys Thr Ala Leu Leu Val Ala Val Val Ala Arg Lys		270
	275	280
Leu Glu Phe Asn Lys Ala Glu Lys His Val His Asn Phe Met Met Asp		285
	290	295
Ile Gln Tyr Thr Lys Glu Met Lys Glu Ser Ala Ala Arg Val Leu Gln		300
305	310	315
Glu Ala Trp Met Phe Tyr Lys His Thr Arg Arg Lys Glu Ser His Ala		320
	325	330
Ala Arg Arg His Gln Arg Lys Leu Leu Ala Ala Ile Asn Ala Phe Arg		335
	340	345
Gln Val Arg Leu Lys His Arg Lys Leu Arg Glu Gln Val Asn Ser Met		350
	355	360
Val Asp Ile Ser Lys Met His Met Ile Leu Tyr Asp Leu Gln Gln Asn		365
	370	375
Leu Ser Ser Ser His Arg Ala Leu Glu Lys Gln Ile Asp Thr Leu Ala		380
385	390	395
Gly Lys Leu Asp Ala Leu Thr Glu Leu Leu Ser Thr Ala Leu Gly Pro		400
	405	410
Arg Gln Leu Pro Glu Pro Ser Gln Gln Ser Lys		415

- 89 -

420

425

## (2) INFORMATION FOR SEQ ID NO:4:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 561 amino acids

(B) TYPE: amino acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

## (ii) MOLECULE TYPE: peptide

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO:4:

```

Met  Pro  Gly  Pro  Arg  Ala  Ala  Cys  Ser  Glu  Pro  Asn  Pro  Cys  Thr  Gln
 1              5              10              15
Val  Val  Met  Asn  Ser  His  Ser  Tyr  Asn  Gly  Ser  Val  Gly  Arg  Pro  Leu
          20          25          30
Gly  Ser  Gly  Pro  Gly  Ala  Leu  Gly  Arg  Asp  Pro  Pro  Asp  Pro  Glu  Ala
          35          40          45
Gly  His  Pro  Pro  Gln  Pro  Pro  His  Ser  Pro  Gly  Leu  Gln  Val  Val  Val
          50          55          60
Ala  Lys  Ser  Glu  Pro  Ala  Arg  Pro  Ser  Pro  Gly  Ser  Pro  Arg  Gly  Gln
65          70          75          80
Pro  Gln  Asp  Gln  Asp  Asp  Asp  Glu  Asp  Asp  Glu  Glu  Asp  Glu  Ala  Gly
          85          90          95
Arg  Gln  Arg  Ala  Ser  Gly  Lys  Pro  Ser  Asn  Val  Gly  His  Arg  Leu  Gly
          100          105          110
His  Arg  Arg  Ala  Leu  Phe  Glu  Lys  Arg  Lys  Arg  Leu  Ser  Asp  Tyr  Ala
          115          120          125
Leu  Ile  Phe  Gly  Met  Phe  Gly  Ile  Val  Val  Met  Val  Thr  Glu  Thr  Glu
          130          135          140
Leu  Ser  Trp  Gly  Val  Tyr  Thr  Lys  Glu  Ser  Leu  Tyr  Ser  Phe  Ala  Leu

```

145					150				155				160			
Lys	Cys	Leu	Ile	Ser	Leu	Ser	Thr	Ala	Ile	Leu	Leu	Gly	Leu	Val	Val	
				165				170				175				
Leu	Tyr	His	Ala	Arg	Glu	Ile	Gln	Leu	Phe	Met	Val	Asp	Asn	Gly	Ala	
				180				185				190				
Asp	Asp	Trp	Arg	Ile	Ala	Met	Thr	Cys	Glu	Arg	Val	Phe	Leu	Ile	Ser	
				195				200				205				
Leu	Glu	Leu	Ala	Val	Cys	Ala	Ile	His	Pro	Val	Pro	Gly	His	Tyr	Arg	
				210				215				220				
Phe	Thr	Trp	Thr	Ala	Arg	Leu	Ala	Phe	Thr	Tyr	Ala	Pro	Ser	Val	Ala	
225					230				235				240			
Glu	Ala	Asp	Val	Asp	Val	Leu	Leu	Ser	Ile	Pro	Met	Phe	Leu	Arg	Leu	
				245				250				255				
Tyr	Leu	Leu	Gly	Arg	Val	Met	Leu	Leu	His	Ser	Lys	Ile	Phe	Thr	Asp	
				260				265				270				
Ala	Ser	Ser	Arg	Ser	Ile	Gly	Ala	Leu	Asn	Lys	Ile	Thr	Phe	Asn	Thr	
				275				280				285				
Arg	Phe	Val	Met	Lys	Thr	Leu	Met	Thr	Ile	Cys	Pro	Gly	Thr	Val	Leu	
				290				295				300				
Leu	Val	Phe	Ser	Ile	Ser	Ser	Trp	Ile	Ile	Ala	Ala	Trp	Thr	Val	Arg	
305					310				315				320			
Val	Cys	Glu	Arg	Tyr	His	Asp	Lys	Gln	Glu	Val	Thr	Ser	Asn	Phe	Leu	
				325				330				335				
Gly	Ala	Met	Trp	Leu	Ile	Ser	Ile	Thr	Phe	Leu	Ser	Ile	Gly	Tyr	Gly	
				340				345				350				
Asp	Met	Val	Pro	His	Thr	Tyr	Cys	Gly	Lys	Gly	Val	Cys	Leu	Leu	Thr	
				355				360				365				
Gly	Ile	Met	Gly	Ala	Gly	Cys	Thr	Ala	Leu	Val	Val	Ala	Val	Val	Ala	
				370				375				380				
Arg	Lys	Leu	Glu	Leu	Thr	Lys	Ala	Glu	Lys	His	Val	His	Asn	Phe	Met	
385					390				395				400			
Met	Asp	Thr	Gln	Leu	Thr	Lys	Arg	Val	Lys	Asn	Ala	Ala	Ala	Asn	Val	



- 91 -

```

          405              410              415
Leu Arg Glu Thr Trp Leu Ile Tyr Lys His Thr Arg Leu Val Lys Lys
          420              425              430
Pro Asp Gln Ala Arg Val Arg Lys His Gln Arg Lys Phe Leu Gln Ala
          435              440              445
Ile His Gln Ala Gln Lys Leu Arg Ser Val Lys Ile Glu Gln Gly Lys
          450              455              460
Leu Asn Asp Gln Ala Asn Thr Leu Thr Asp Leu Ala Lys Thr Gln Thr
          465              470              475              480
Val Met Tyr Asp Leu Val Ser Glu Leu His Ala Gln His Glu Glu Leu
          485              490              495
Glu Ala Arg Leu Ala Thr Leu Glu Ser Arg Leu Asp Ala Leu Gly Ala
          500              505              510
Ser Leu Gln Ala Leu Pro Gly Leu Ile Ala Gln Ala Ile Arg Pro Pro
          515              520              525
Pro Pro Pro Leu Pro Pro Arg Pro Gly Pro Gly Pro Gln Asp Gln Ala
          530              535              540
Ala Arg Ser Ser Pro Cys Arg Trp Thr Pro Val Ala Pro Ser Asp Cys
          545              550              555              560
Gly

```

## (2) INFORMATION FOR SEQ ID NO:5:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 579 amino acids

(B) TYPE: amino acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO:5:

```

Met Ser Ser Cys Arg Tyr Asn Gly Gly Val Met Arg Pro Leu Ser Asn
 1              5              10              15

```

- 92 -

Leu Ser Ser Ser Arg Arg Asn Leu His Glu Met Asp Ser Glu Ala Gln  
                   20                                  25                                  30  
 Pro Leu Gln Pro Pro Ala Ser Val Val Gly Gly Gly Gly Gly Ala Ser  
                   35                                  40                                  45  
 Ser Pro Ser Ala Ala Ala Ala Ala Ser Ser Ser Ala Pro Glu Ile Val  
                   50                                  55                                  60  
 Val Ser Lys Pro Glu His Asn Asn Ser Asn Asn Leu Ala Leu Tyr Gly  
 65                                  70                                  75                                  80  
 Thr Gly Gly Gly Gly Gly Ser Thr Gly Gly Gly Gly Gly Gly Gly Gly Gly  
                                   85                                  90                                  95  
 Gly Gly Gly Ser Gly His Gly Ser Ser Ser Gly Thr Lys Ser Ser Lys  
                   100                                  105                                  110  
 Lys Lys Asn Gln Asn Ile Gly Tyr Lys Leu Gly His Arg Arg Ala Leu  
                   115                                  120                                  125  
 Phe Lys Arg Lys Arg Leu Ser Asp Tyr Ala Leu Ile Phe Gly Met Phe  
                   130                                  135                                  140  
 Gly Ile Val Val Met Val Ile Glu Thr Glu Leu Ser Trp Gly Ala Tyr  
 145                                  150                                  155                                  160  
 Asp Lys Ala Ser Leu Tyr Ser Leu Ala Leu Lys Cys Leu Ile Ser Leu  
                                   165                                  170                                  175  
 Ser Thr Ile Ile Leu Leu Gly Leu Ile Ile Val Tyr His Ala Arg Glu  
                   180                                  185                                  190  
 Ile Gln Leu Phe Met Val Asp Asn Gly Ala Asp Asp Trp Arg Ile Ala  
                   195                                  200                                  205  
 Met Thr Tyr Glu Arg Ile Phe Phe Ile Cys Leu Glu Ile Leu Val Cys  
                   210                                  215                                  220  
 Ala Ile His Pro Ile Pro Gly Asn Tyr Thr Phe Thr Trp Thr Ala Arg  
 225                                  230                                  235                                  240  
 Leu Ala Phe Ser Tyr Ala Pro Ser Thr Thr Thr Ala Asp Val Asp Ile  
                                   245                                  250                                  255  
 Ile Leu Ser Ile Pro Met Phe Leu Arg Leu Tyr Leu Ile Ala Arg Val  
                   260                                  265                                  270

- 93 -

Met Leu Leu His Ser Lys Leu Phe Thr Asp Ala Ser Ser Arg Ser Ile  
 275 280 285  
 Gly Ala Leu Asn Lys Ile Asn Phe Asn Thr Arg Phe Val Met Lys Thr  
 290 295 300  
 Leu Met Thr Ile Cys Pro Gly Thr Val Leu Leu Val Phe Ser Ile Ser  
 305 310 315 320  
 Leu Trp Ile Ile Ala Ala Trp Thr Val Arg Ala Cys Glu Arg Tyr His  
 325 330 335  
 Asp Gln Gln Asp Val Thr Ser Asn Phe Leu Gly Ala Met Trp Leu Ile  
 340 345 350  
 Ser Ile Thr Phe Leu Ser Ile Gly Tyr Gly Asp Met Val Pro Asn Thr  
 355 360 365  
 Tyr Cys Gly Lys Gly Val Cys Leu Leu Thr Gly Ile Met Gly Ala Gly  
 370 375 380  
 Cys Thr Ala Leu Val Val Ala Val Val Ala Arg Lys Leu Glu Leu Thr  
 385 390 395 400  
 Lys Ala Glu Lys His Val His Asn Phe Met Met Asp Thr Gln Leu Thr  
 405 410 415  
 Lys Arg Val Lys Asn Ala Ala Ala Asn Val Leu Arg Glu Thr Trp Leu  
 420 425 430  
 Ile Tyr Lys Asn Thr Lys Leu Val Lys Lys Ile Asp His Ala Lys Val  
 435 440 445  
 Arg Lys His Gln Arg Lys Phe Leu Gln Ala Ile His Gln Leu Arg Ser  
 450 455 460  
 Val Lys Met Glu Gln Arg Lys Leu Asn Asp Gln Ala Asn Thr Leu Val  
 465 470 475 480  
 Asp Leu Ala Lys Thr Gln Asn Ile Met Tyr Asp Met Ile Ser Asp Leu  
 485 490 495  
 Asn Glu Arg Ser Glu Asp Phe Glu Lys Arg Ile Val Thr Leu Glu Thr  
 500 505 510  
 Lys Leu Glu Thr Leu Ile Gly Ser Ile His Ala Leu Pro Gly Leu Ile  
 515 520 525

- 94 -

Ser Gln Thr Ile Arg Gln Gln Gln Arg Asp Phe Ile Glu Thr Gln Met  
 530 535 540  
 Glu Asn Tyr Asp Lys His Val Thr Tyr Asn Ala Glu Arg Ser Arg Ser  
 545 550 555 560  
 Ser Ser Arg Arg Arg Arg Ser Ser Ser Thr Ala Pro Pro Thr Ser Ser  
 565 570 575  
 Glu Ser Ser

## (2) INFORMATION FOR SEQ ID NO:6:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 551 amino acids
- (B) TYPE: amino acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO:6:

Met Ser Ser Cys Lys Tyr Ser Gly Gly Val Met Lys Pro Leu Ser Arg  
 1 5 10 15  
 Leu Ser Ala Ser Arg Arg Asn Leu Ile Glu Ala Glu Pro Glu Gly Gln  
 20 25 30  
 Pro Leu Gln Leu Phe Ser Pro Ser Asn Pro Pro Glu Ile Ile Ile Ser  
 35 40 45  
 Ser Arg Glu Asp Asn His Ala His Gln Thr Leu Leu His His Pro Asn  
 50 55 60  
 Ala Thr His Asn His Gln His Ala Gly Thr Thr Ala Gly Ser Thr Thr  
 65 70 75 80  
 Phe Pro Lys Ala Asn Lys Arg Lys Asn Gln Asn Ile Gly Tyr Lys Leu  
 85 90 95  
 Gly His Arg Arg Ala Leu Phe Glu Lys Arg Lys Arg Leu Ser Asp Tyr  
 100 105 110

- 95 -

Ala Leu Ile Phe Gly Met Phe Gly Ile Val Val Met Val Ile Glu Thr  
 115 120 125  
 Glu Leu Ser Trp Gly Leu Tyr Ser Lys Asp Ser Met Phe Ser Leu Ala  
 130 135 140  
 Leu Lys Cys Leu Ile Ser Leu Ser Thr Ile Ile Leu Leu Gly Leu Ile  
 145 150 155 160  
 Ile Ala Tyr His Thr Arg Glu Val Gln Leu Phe Val Ile Asp Asn Gly  
 165 170 175  
 Ala Asp Asp Trp Arg Ile Ala Met Thr Tyr Glu Arg Ile Leu Tyr Ile  
 180 185 190  
 Ser Leu Glu Met Leu Val Cys Ala Ile His Pro Ile Pro Gly Glu Tyr  
 195 200 205  
 Lys Phe Phe Trp Thr Ala Arg Leu Ala Phe Ser Tyr Thr Pro Ser Arg  
 210 215 220  
 Ala Glu Ala Asp Val Asp Ile Ile Leu Ser Ile Pro Met Phe Leu Arg  
 225 230 235 240  
 Leu Tyr Leu Ile Ala Arg Val Met Leu Leu His Ser Lys Leu Phe Thr  
 245 250 255  
 Asp Ala Ser Ser Arg Ser Ile Gly Ala Leu Asn Lys Ile Asn Thr Arg  
 260 265 270  
 Phe Val Met Lys Thr Leu Met Thr Ile Cys Pro Gly Thr Val Leu Leu  
 275 280 285  
 Met Phe Ser Ile Ser Leu Trp Ile Ile Ala Ala Trp Thr Val Arg Val  
 290 295 300  
 Cys Glu Arg Tyr His Asp Gln Gln Asp Val Thr Ser Asn Phe Leu Gly  
 305 310 315 320  
 Ala Met Trp Leu Ile Ser Ile Thr Phe Leu Ser Ile Gly Tyr Gly Asp  
 325 330 335  
 Met Val Pro His Thr Tyr Cys Gly Lys Gly Val Cys Leu Leu Thr Gly  
 340 345 350  
 Ile Met Gly Ala Gly Cys Thr Ala Leu Val Val Ala Val Val Ala Arg  
 355 360 365

- 96 -

Lys Leu Glu Leu Thr Lys Ala Glu Lys His Val His Asn Phe Met Met  
 370 375 380  
 Asp Thr Gln Leu Thr Lys Arg Ile Lys Asn Ala Ala Ala Asn Val Leu  
 385 390 395 400  
 Arg Glu Thr Trp Leu Ile Tyr Lys His Thr Lys Leu Leu Lys Lys Ile  
 405 410 415  
 Asp His Ala Lys Val Arg Lys His Gln Arg Lys Phe Leu Gln Ala Ile  
 420 425 430  
 His Gln Leu Arg Gly Val Lys Met Glu Gln Arg Lys Leu Ser Asp Gln  
 435 440 445  
 Ala Asn Thr Leu Val Asp Leu Ser Lys Met Gln Asn Val Met Tyr Asp  
 450 455 460  
 Leu Ile Thr Glu Leu Asn Asp Arg Ser Glu Asp Leu Glu Lys Gln Ile  
 465 470 475 480  
 Gly Ser Leu Glu Ser Lys Leu Glu His Leu Thr Ala Ser Phe Asn Ser  
 485 490 495  
 Leu Pro Leu Leu Ile Ala Asp Thr Leu Arg Gln Gln Gln Gln Gln Leu  
 500 505 510  
 Leu Thr Ala Phe Val Glu Ala Arg Gly Ile Ser Val Ala Val Gly Thr  
 515 520 525  
 Ser His Ala Pro Pro Ser Asp Ser Pro Ile Gly Ile Ser Ser Thr Ser  
 530 535 540  
 Phe Pro Glu Phe Leu Ile Phe  
 545 550

(2) INFORMATION FOR SEQ ID NO:7:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 2238 base pairs

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:7:

GTCTTCGGT GTCTGGGTGT GGTGAGTAGA GGTGTGTGC ACAAAGTACA GACCATTGTG 60  
TGTGACAAAG CCCATCGTGT GTCTGTGTGT GTCTTTATCC ACGTGGATGG ACGTCTCTTT 120  
CTTGCTCTGC CCCAAGACAC ACCCTAGCCC CTCCTTATTC TCAAAAGGGG GAGCTGGGGA 180  
GCCTCCCCCT ACCCTGGGGC CTCCCTGCC CCTCCCCGCC CTGCCTGGCC GTCACCACTC 240  
CCCAGAGGGC ACAGGGCTCT GCTGTGCCTC AGAGCAAAAG TCCCAGAGCC AGCAGAGCAG 300  
GCTGACGACC TGCAAGCCAC AGTGGCTGCC CTGTGCGTGC TGCGAGGTGG GGGACCCTGG 360  
GCAGGAAGCT GGCTGAGCCC CAAGACCCCG GGGGCCATGG GCGGGGATCT GGTGCTTGGC 420  
CTGGGGGCCT TGAGACGCCG AAAGCGCTTG CTGGAGCAGG AGAAGTCTCT GGCCGGCTGG 480  
GCACTGGTGC TGGCAGGAAC TGGCATTGGA CTCATGGTGC TGCATGCAGA GATGCTGTGG 540  
TTCGGGGGT GCTCGTGGG GCTCTACCTG TTCTGGTTA AATGCACGAT CAGCATTTC 600  
ACCTTCTTAC TCCTCTGCCT CATCGTGGCC TTTCATGCCA AAGAGGTCCA GCTGTTTCATG 660  
ACCGACAACG GGCTGCGGGA CTGGCGCGTG GCGCTGACCG GCGGGCAGGC GCGCGAGATC 720  
GTGCTGGAGC TGGTGGTGTG TGGGCTGCAC CCGGCGCCCG TGCGGGGCCC GCCGTGCGTG 780  
CAGGATTTAG GGGCGCCGCT GACCTCCCCG CAGCCCTGGC CGGGATTCTT GGGCCAAGGG 840  
GAAGCGCTGC TGTCCTGGC CATGCTGCTG CGTCTCTACC TGGTGCCCCG CGCCGTGCTC 900  
CTGCGCAGCG GCGTCTGCT CAACGCTTCC TACCGCAGCA TCGGCGCTCT CAATCAAGTC 960  
CGCTTCCGCC ACTGGTTCGT GGCCAAGCTT TACATGAACA CGCACCTGG CCGCTGCTG 1020  
CTCGGCTCA CGCTTGGCT CTGGCTGACC ACCGCTGGG TGCTGTCCGT GGCCGAGAGG 1080  
CAGGCTGTTA ATGCCACTGG GCACCTTTCA GACACACTTT GGCTGATCCC CATCACATTC 1140  
CTGACCATCG GCTATGGTGA CGTGGTGCCG GGCACCATGT GGGGCAAGAT CGTCTGCCTG 1200  
TGCACTGGAG TCATGGGTGT CTGCTGCACA GCCCTGCTGG TGGCCGTGGT GGCCCGGAAG 1260  
CTGGAGTTA ACAAGGCAGA GAAGCACGTG CACAACCTCA TGATGGATAT CCAGTATACC 1320  
AAAGAGATGA AGGAGTCCCG TGCCCCAGTG CTACAAGAAG CCTGGATGTT CTACAAACAT 1380  
ACTCGCAGGA AGGAGTCTCA TGCTGCCCCG AGGCATCAGC GCAAGCTGCT GGCCGCCATC 1440  
AACGCGTCC GCCAGGTGCG GCTGAAACAC CGGAAGCTCC GGGAACAAGT GAACTCCATG 1500  
GTGGACATCT CCAAGATGCA CATGATCCTG TATGACCTGC AGCAGAATCT GAGCAGCTCA 1560  
CACCGGGCCC TGGAGAAACA GATTGACACG CTGGCGGGGA AGCTGGATGC CCTGACTGAG 1620  
CTGCTTAGCA CTGCCCTGGG GCCGAGGCAG CTTCCAGAAC CCAGCCAGCA GTCCAAGTAG 1680  
CTGGACCAC GAGGAGGAAC CAGGCTACTT TCCCAGTAC TGAGGTGGTG GACATCGTCT 1740  
CTGCCACTCC TGACCCAGCC CTGAACAAAG CACCTCAAGT GCAAGGACCA AAGGGGGCCC 1800  
TGGCTTGGAG TGGGTTGGCT TGCTGATGGC TGCTGGAGGG GACGCTGGCT AAAGTGGGTA 1860  
GGCCTTGGCC CACCTGAGGC CCCAGGTGGG AACATGGTCA CCCCCACTCT GCATACCCTC 1920  
ATCAAAAACA CTCTACTAT GCTGCTATGG ACGACCTCCA GCTCTCAGTT ACAAGTGCAG 1980  
GCGACTGGAG GCAGGACTCC TGGGTCCCTG GGAAAGAGGG TACTAGGGGC CCGGATCCAG 2040  
GATTCTGGGA GGCTTCAGTT ACCGCTGGCC GAGCTGAAGA ACTGGGTATG AGGCTGGGGC 2100  
GGGGCTGGAG GTGGCGCCCC CTGGTGGGAC AACAAAGAGG ACACCATTTT TCCAGAGCTG 2160  
CAGAGAGCAC CTGGTGGGGA GGAAGAAGTG TAACTACCA GCCTCTGCTC TTATCTTTGT 2220  
AATAAATGTT AAAGCCAG 2238

- 98 -

## (2) INFORMATION FOR SEQ ID NO:8:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 1250 base pairs  
 (B) TYPE: nucleic acid  
 (C) STRANDEDNESS: single  
 (D) TOPOLOGY: linear

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO:8:

```

CAGACCTACT GAATCCTGGA GGGTAGGGCC AGCAATCCAC TTAGCAAGC CTTCCAGGCA   60
GGCTGTTAAT GCCACTGGGC ACCTTTCAGA CACACTTTGG CTGATCCCCA TCACATTCCT  120
GACCATCGGC TATGGTGACG TGGTGCCGGG CACCATGTGG GGCAAGATCG TCTGCCTGTG  180
CACTGGAGTC ATGGGTGTCT GCTGCACAGC CTTGCTGGTG GCCGTGGTGG CCCGGAAGCT  240
GGAGTTTAAC AAGGCAGAGA AGCACGTGCA CAACTTCATG ATGGATATCC AGTATACCAA  300
AGAGATGAAG GAGTCCGCTG CCCGAGTGCT ACAAGAAGCC TGGATGTTCT ACAAACATAC  360
TCGCAGGAAG GAGTCTCATG CTGCCCCGAG GCATCAGCGC AAGCTGCTGG CCGCCATCAA  420
CGCGTTCCGC CAGGTGCGGC TGAAACACCG GAAGCTCCGG GAACAAGTGA ACTCCATGGT  480
GGACATCTCC AAGATGCACA TGATCCTGTA TGACCTGCAG CAGAATCTGA GCAGCTCACA  540
CCGGGGCCCTG GAGAAACAGA TTGACACGCT GCGGGGAAG CTGGATGCCC TGA CTGAGCT  600
GCTTAGCACT GCCCTGGGGC CGAGGCAGCT TCCAGAACCC AGCCAGCAGT CCAAGTAGCT  660
GGACCCACGA GGAGGAACCA GGCTACTTTC CCCAGTACTG AGGTGGTGGA CATCGTCTCT  720
GCCACTCCTG ACCCAGCCCT GAACAAAGCA CCTCAAGTGC AAGGACCAAA GGGGGCCCTG  780
GCTTGGAGTG GGTGCTGCTG CTGATGGCTG CTGGAGGGGA CGCTGGCTAA AGTGGGTAGG  840
CCTTGGCCCC CTTGAGGCCC CAGGTGGGAA CATGGTCACC CCCACTCTGC ATACCCTCAT  900
CAAAAACACT CTCACTATGC TGCTATGGAC GACCTCCAGC TCTCAGTTAC AAGTGCAGGC  960
GACTGGAGGC AGGACTCCTG GGTCCCTGGG AAAGAGGGTA CTAGGGGCCC GGATCCAGGA 1020
TTCTGGGAGG CTTCAATTAC CGCTGGCCGA GCTGAAGAAC TGGGTATGAG GCTGGGGCGG 1080
GGCTGGAGGT GGCGCCCCCT GGTGGGACAA CAAAGAGGAC ACCATTTTTC CAGAGCTGCA 1140
GAGAGCAACT GGTGGGGAGG AAGAAGTGTA AATCACCAGC CTTTGCTTTT ATTTTGTAA 1200
TAAATGTAA AGCCAGXAAA AAAAAAAAAA AAAAAAAAAA AAAAAAAAAA 1250

```

## (2) INFORMATION FOR SEQ ID NO:9:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 1624 base pairs  
 (B) TYPE: nucleic acid  
 (C) STRANDEDNESS: single  
 (D) TOPOLOGY: linear



- 99 -

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:9:

```

CAAAGAGGTC CAGCTGTTCA TGACCGACAA CGGGCTGCGG GACTGGCGCG TGGCGCTGAC   60
CGGGCGGCAG GCGGCGCAGA TCGTGCTGGA GCTGGTGGTG TGTGGGCTGC ACCCGGCGCC   120
CGTGCGGGGC CCGCCGTGCG TGCAGGATTT AGGGGCGCCG CTGACCTCCC CGCAGCCCTG   180
GCCGGGATTC CTGGGCCAAG GGGAAGCGCT GCTGTCCCTG GCCATGCTGC TGCCTCTCTA   240
CCTGGTGCCC CGCGCCGTGC TCCTGCGCAG CGGCGTCCTG CTCAACGCTT CCTACCGCAG   300
CATCGGCCTG CTCAATCAAG TCCGCTTCCG CCACTGGTTC GTGGCCAAGC TTTACATGAA   360
CACGCACCTT GGCCGCTGCG TGCTCGGCCT CACGCTTGGC CTCTGGCTGA CCACCGCCTG   420
GGTGCTGTCC GTGGCCGAGA GGCAGGCTGT TAATGCCACT GGGCACCTTT CAGACACACT   480
TTGGCTGATC CCCATCACAT TCCTGACCAT CGGCTATGGT GACGTGGTGC CGGGCACCAT   540
GTGGGGCAAG ATCGTCTGCC TGTGCACTGG AGTCATGGGT GTCTGCTGCA CAGCCCTGCT   600
GGTGGCCGTG GTGGCCCGGA AGCTGGAGTT TAACAAGGCA GAGAAGCACG TGCACAACCT   660
CATGATGGAT ATCCAGTATA CCAAAGAGAT GAAGGAGTCC GCTGCCCCGAG TGCTACAAGA   720
AGCCTGGATG TTCTACAAAC ATACTCGCAG GAAGGAGTCT CATGCTGCCC GCAGGCATCA   780
GCGCAAGCTG CTGGCCGCCA TCAACGCGTT CCGCCAGGTG CGGCTGAAAC ACCGGAAGCT   840
CCGGGAACAA GTGAACTCCA TGGTGGACAT CTCCAAGATG CACATGATCC TGTATGACCT   900
GCAGCAGAAT CTGAGCAGCT CACACCGGGC CCTGGAGAAA CAGATTGACA CGCTGGCGGG   960
GAAGCTGGAT GGCCTGACTG AGCTGCTTAG CACTGCCCTG GGGCCGAGGC AGCTTCCAGA   1020
ACCCAGCCAG CAGTCCAAGT AGCTGGACCC ACGAGGAGGA ACCAGGCTAC TTTCCCCAGT   1080
ACTGAGGTGG TGGACATCGT CTCTGCCACT CCTGACCCAG CCCTGAACAA AGCACCTCAA   1140
GTGCAAGGAC CAAAGGGGGC CCTGGCTTGG AGTGGGTTGG CTTGCTGATG GCTGCTGGAG   1200
GGGACGCTGG CTAAAGTGGG TAGGCCTTGG CCCACCTGAG GCCCCAGGTG GGAACATGGT   1260
CACCCCCACT CTGCATACCC TCATCAAAAA CACTCTCACT ATGCTGCTAT GGACGACCTC   1320
CAGCTCTCAG TTACAAGTGC AGGCGACTGG AGGCAGGACT CCTGGGTCCC TGGGAAAGAG   1380
GGTAATAGGG GCCCGGATXC AGGATTCTGG GAGGCTTCAG TTACCGCTGG CCGAGCTGAA   1440
GAACTGGGTA TGAGGCTGGG GCGGGGCTGG AGGTGGCGCC CCCTGGTGGG ACAACAAAGA   1500
GGACACCATT TTTCCAGAGC TGCAGAGAGC AXCTGGTGGG GAGGAAGAAG TGTAACCTAC   1560
CAGCCTTTGC TCTTATTTTT GTAATAAATG TAAAGCCAG AAAAAAAAAA AAAAAAAAAA   1620
AAAA

```

1624

(2) INFORMATION FOR SEQ ID NO:10:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 amino acids

(B) TYPE: amino acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

- 100 -

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:10:

Pro Ala Ser Phe Trp Trp Ala Thr Ile Thr Met Thr Thr Val Gly Tyr  
1                      5                      10                      15  
Gly Asp Ile Tyr Pro  
                    20

(2) INFORMATION FOR SEQ ID NO:11:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 24 base pairs
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:11:

GTGGACATCT CCAAGGTACT AGGA

24

(2) INFORMATION FOR SEQ ID NO:12:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 24 base pairs
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:12:

TTTACTGACA GGCTGTGTGT GCCA

24

(2) INFORMATION FOR SEQ ID NO:13:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 1381 base pairs
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

- 101 -

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:13:

CAAGCTGGCT GCCATTTGCG TACTGAGAGG TGGGGGACCT GGGCAGGAAG CTGGCTGAGC 60  
CCCAAGACCT CAGGGGCCAT GGGCGGGGAG CTGGTGA CTG GCCTGGGGGC CCTGAGACGG 120  
AGAAAGCGCC TGCTGGAGCA GGAGAAGAGG GTGGCCGGCT GGGCGTTGGT GCTGGCGGGA 180  
ACTGGCATCG GACTCATGGT TCTGCACGCT GAGATGTTGT GGTTCTGGG CTGCAAGTGG 240  
GTGCTGTACC TGCTCCTGGT TAAGTGT TTG ATCACCTGT CCACTGCCTT CCTCCTTTGT 300  
CTTATTGTGG TCTTCCATGC CAAGGAGGTC CAGCTGTTCA TGA CTGACAA CGGGCTCCGG 360  
GACTGGCGCG TGGCGCTGAC CCGCGGCAG GTGGCGCAGA TCCTGCTGGA GCTGTTGGTG 420  
TGCGGGGTGC ACCCGGTGCC CCTACGGAGC CCGCACTGCG CCCTGGCGGG GGAGGCCACC 480  
GACGCGCAGC CCTGGCCGGG TTCTCTGGGC GAAGGCGAGG CGTTGCTGTC CCTGGCCATG 540  
CTCCTGCGTC TCTACCTGGT GCGCCGCGCG GTGCTGCTGC GCAGCGGGGT CCTGCTCAAC 600  
GCGTCTTACC GCAGCATCGG GCGCTCAAC CAAGTCCGCT TCCGCCACTG GTTCGTGGCC 660  
AAGCTGTACA TGAACACGCA CCCGGGTCGC CTGCTGCTGG GCCTACGCT GGGTCTCTGG 720  
CTCACACAG CTTGGGTGCT GTCTGTGGCT GAGAGGCAGG CTGTCAATGC CACGGGGCAC 780  
CTCACAGACA CACTGTGGCT GATTCCGATC ACATTCCTGA CCATTGGCTA TGGGGACGTG 840  
GTACCTGGCA CCATGTGGGG CAAGATTGTC TGCCTGTGCA CCGGAGTCAT GGGGGTCTGC 900  
TGCACAGCTC TCCTGGTGGC TGTGGTGGCT CGGAAGCTGG AGTTCAACAA GGCGGAGAAA 960  
CACGTGCACA ACTTCATGAT GGACATCCAT TATGCCAAAG AGATGAAGGA GTCAGCGGCG 1020  
CGGCTGCTGC AGGAAGCCTG GATGTACTAC AAGCACACTC GAAGGAAGGA CTCCCGGGCT 1080  
GCCCCGAGAC ATCAGCGCAA GATGCTGGCC GCCATCCACA CGTTCCGCCA GGTACGGCTG 1140  
AAACACCGGA AGTCCGGGA ACAAGTGAAT TCCATGGTGG ACATCTCAA GATGCACATG 1200  
ATCCTGTGCG ACCTGCAGCT GGGTCTCAGC TCCTCGCACC GTGCCCTGGA GAAGAGAATC 1260  
GACGGTCTGG CAGGAAAGCT GGATGCCCTG ACAGAGCTGC TCGGCACTGC TCTGCAGCAA 1320  
CAGCAGCTAC CAGAACCCAG TCAGGAGGCC ACATAGCTCC ACATGAACTC ACAGAAGAAC 1380  
A 1381

(2) INFORMATION FOR SEQ ID NO:14:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 425 amino acids

(B) TYPE: amino acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:14:

- 102 -

Met Gly Gly Glu Leu Val Thr Gly Leu Gly Ala Leu Arg Arg Arg Lys  
 1 5 10 15  
 Arg Leu Leu Glu Gln Glu Lys Arg Val Ala Gly Trp Ala Leu Val Leu  
 20 25 30  
 Ala Gly Thr Gly Ile Gly Leu Met Val Leu His Ala Glu Met Leu Trp  
 35 40 45  
 Phe Leu Gly Cys Lys Trp Val Leu Tyr Leu Leu Leu Val Lys Cys Leu  
 50 55 60  
 Ile Thr Leu Ser Thr Ala Phe Leu Leu Cys Leu Ile Val Val Phe His  
 65 70 75 80  
 Ala Lys Glu Val Gln Leu Phe Met Thr Asp Asn Gly Leu Arg Asp Trp  
 85 90 95  
 Arg Val Ala Leu Thr Arg Arg Gln Val Ala Gln Ile Leu Leu Glu Leu  
 100 105 110  
 Leu Val Cys Gly Val His Pro Val Pro Leu Arg Ser Pro His Cys Ala  
 115 120 125  
 Leu Ala Gly Glu Ala Thr Asp Ala Gln Pro Trp Pro Gly Phe Leu Gly  
 130 135 140  
 Glu Gly Glu Ala Leu Leu Ser Leu Ala Met Leu Leu Arg Leu Tyr Leu  
 145 150 155 160  
 Val Pro Arg Ala Val Leu Leu Arg Ser Gly Val Leu Leu Asn Ala Ser  
 165 170 175  
 Tyr Arg Ser Ile Gly Ala Leu Asn Gln Val Arg Phe Arg His Trp Phe  
 180 185 190  
 Val Ala Lys Leu Tyr Met Asn Thr His Pro Gly Arg Leu Leu Leu Gly  
 195 200 205  
 Leu Thr Leu Gly Leu Trp Leu Thr Thr Ala Trp Val Leu Ser Val Ala  
 210 215 220  
 Glu Arg Gln Ala Val Asn Ala Thr Gly His Leu Thr Asp Thr Leu Trp  
 225 230 235 240  
 Leu Ile Pro Ile Thr Phe Leu Thr Ile Gly Tyr Gly Asp Val Val Pro  
 245 250 255

- 103 -

Gly Thr Met Trp Gly Lys Ile Val Cys Leu Cys Thr Gly Val Met Gly  
                   260                                  265                                  270  
 Val Cys Cys Thr Ala Leu Leu Val Ala Val Val Ala Arg Lys Leu Glu  
                   275                                  280                                  285  
 Phe Asn Lys Ala Glu Lys His Val His Asn Phe Met Met Asp Ile His  
                   290                                  295                                  300  
 Tyr Ala Lys Glu Met Lys Glu Ser Ala Ala Arg Leu Leu Gln Glu Ala  
 305                                  310                                  315                                  320  
 Trp Met Tyr Tyr Lys His Thr Arg Arg Lys Asp Ser Arg Ala Ala Arg  
                                   325                                  330                                  335  
 Arg His Gln Arg Lys Met Leu Ala Ala Ile His Thr Phe Arg Gln Val  
                   340                                  345                                  350  
 Arg Leu Lys His Arg Lys Leu Arg Glu Gln Val Asn Ser Met Val Asp  
                   355                                  360                                  365  
 Ile Ser Lys Met His Met Ile Leu Cys Asp Leu Gln Leu Gly Leu Ser  
                   370                                  375                                  380  
 Ser Ser His Arg Ala Leu Glu Lys Arg Ile Asp Gly Leu Ala Gly Lys  
 385                                  390                                  395                                  400  
 Leu Asp Ala Leu Thr Glu Leu Leu Gly Thr Ala Leu Gln Gln Gln Gln  
                                   405                                  410                                  415  
 Leu Pro Glu Pro Ser Gln Glu Ala Thr  
                   420                                  425

-104-

**WHAT IS CLAIMED IS**

1. A purified polynucleotide comprising a nucleic acid sequence encoding the polypeptide having the sequence substantially as depicted in SEQ ID NO:3 or a biologically active fragment thereof.
2. The polynucleotide of Claim 1 wherein the polynucleotide sequence substantially comprises SEQ ID NO:2.
3. An expression vector comprising the polynucleotide of Claim 1.
4. An antisense molecule comprising the complement of the polynucleotide of Claim 2 or a biologically-effective portion thereof.
5. A host cell transformed with the expression vector of Claim 3.
6. A purified polypeptide comprising the amino acid sequence substantially as depicted in SEQ ID NO:3.
7. An antibody specific for the polypeptide of Claim 6.
8. A method for producing cells which express a biologically active polypeptide substantially as depicted in SEQ ID NO:3, said method comprising culturing a host cell according to Claim 5 under conditions suitable for the expression of said polypeptide.
9. A method for producing a polypeptide having the amino acid sequence substantially as depicted in SEQ ID NO:3, said method comprising the steps of: culturing a host cell according to Claim 5 under conditions suitable for the expression of said polypeptide, and recovering said polypeptide from the host cell culture.

-105-

10. A method of identifying compounds that modulate the activity of a potassium channel, comprising:
  - (a) combining a candidate compound modulator of a potassium channel activity with a polypeptide of a potassium channel having the sequence substantially as depicted in SEQ ID NO:3, and
  - (b) measuring an effect of the candidate compound modulator on the channel.
11. A method of identifying compounds that modulate the activity of a potassium channel, comprising:
  - (a) combining a candidate compound modulator of a potassium channel activity with a host-cell expressing the polypeptide of a potassium channel having the sequence substantially as depicted in SEQ ID NO:3, and
  - (b) measuring an effect of the candidate compound modulator on the channel.
12. A compound that modulates the activity of a potassium channel identified by the method of Claim 11.
13. A pharmaceutical composition comprising a compound that modulates the activity of a potassium channel identified by the method of Claim 11.
14. A method of treatment of a patient in need of such treatment for a condition which is mediated by a potassium channel, comprising administration of a potassium channel modulating compound according to Claim 12.
15. A method for inhibiting the expression of a potassium channel in a cell comprising administering an effective amount of the antisense molecule of Claim 4 to said cell.

-106-

16. A diagnostic composition for the identification of a polypeptide sequence comprising the amino acid sequence substantially as depicted in SEQ ID NO:3 comprising the antibody of Claim 7.
17. A method of screening a plurality of compounds for binding affinity with the potassium channel subunit polypeptide of Claim 6, said method comprising the steps of:
- a) providing a plurality of compounds;
  - b) combining the compounds with the potassium channel subunit for a time sufficient for the compound to bind the subunit; and
  - c) detecting and recovering the compound which binds the potassium channel subunit.
18. A method of treatment of a patient in need of such treatment for a condition which is mediated by the biological activity of a potassium channel, comprising administration of a subunit substantially as depicted in SEQ ID NO:3 or a biologically active fragment thereof.
19. A method of treatment of a patient in need of such treatment for a condition which is mediated by the biological activity of a potassium channel, comprising administration of a nucleic acid selected from the group consisting essentially of (SEQ ID NO:1, SEQ ID NO:2, and SEQ ID NO:7) or a biologically active fragment thereof.
20. A method of treatment of a patient in need of such treatment for a condition which is mediated by the biological activity of a potassium channel, comprising administration of an antisense molecule according to Claim 4 or a biologically-effective fragment thereof.
21. A compound that modulates the biological activity of native hKCa4 (SEQ ID NO:3) identified by the method of Claim 17.



-107-

22. A method of treatment of a patient in need of such treatment for a condition which is mediated by the biological activity of a human potassium channel, comprising administration of a modulating compound according to Claim 21.

23. A purified polynucleotide comprising a nucleic acid sequence encoding the polypeptide having the sequence substantially as depicted in SEQ ID NO:14 or a biologically active fragment thereof.

24. A purified polypeptide comprising the amino acid sequence substantially as depicted in SEQ ID NO:14.

25. A method of identifying compounds that modulate the activity of a potassium channel, comprising:

(a) combining a candidate compound modulator of a potassium channel activity with a polypeptide of a potassium channel having the sequence substantially as depicted in SEQ ID NO:14, and

(b) measuring an effect of the candidate compound modulator on the channel.

26. A method of identifying compounds that modulate the activity of a potassium channel, comprising:

(a) combining a candidate compound modulator of a potassium channel activity with a host-cell expressing the polypeptide of a potassium channel having the sequence substantially as depicted in SEQ ID NO:14, and

(b) measuring an effect of the candidate compound modulator on the channel.

27. A compound that modulates the activity of a potassium channel identified by the method of Claim 25.

-108-

28. A compound that modulates the activity of a calcium-activated potassium channel identified by the method of Claim 26.

1/29

FIG.1

5  
5'  
GTCCTTCGGTGTCTGGGTGTGGTGAGTAGAGGTGTGTGTGCACAAAGTACAGACCATTTGTG  
TGTGACAAAGCCCATCGTGTGTCTGTGTGTGTCTTTATCCACGTGGATGGACGTCTCTTTC  
TTGCTCTGCCCCAAGACACACCCTAGCCCCCTCCTTATTCTCAAAGGGGGAGCTGGGGAG  
10 CCTCCCCCTACCCTGGGGCCTCCCTGCCCCCTCCCCGCCCTGCCTGGCCGTCACCACTCCC  
CAGAGGGGCACAGGGCTCTGCTGTGCCTCAGAGCAAAAGTCCCAGAGCCAGCAGAGCAGG  
TGACGACCTGCAAGCCACAGTGGCTGCCCTGTGCGTGTGCGAGGTGGGGGACCCTGGGC  
AGGAAGCTGGCTGAGCCCCAAGACCCCGGGGGCCATGGGCGGGGATCTGGTGCTTGGCC  
TGGGGGCCTTGAGACGCCGAAAGCGCTTGCTGGAGCAGGAGAAGTCTCTGGCCGGCTGG  
15 GCACTGGTGCTGGCAGGAAGTGGCATTGGACTCATGGTGCTGCATGCAGAGATGCTGTGG  
TTCGGGGGGTGCTCGTGGGCGCTCTACCTGTTCTGGTTAAATGCACGATCAGCATTTC  
CCTTCTTACTCCTCTGCCTCATCGTGGCCTTTCATGCCAAAGAGGTCCAGCTGTTTCATGAC  
CGACAACGGGCTGCGGGACTGGCGCGTGGCGCTGACCGGGCGGCAGGCGGCGCAGATCG  
TGCTGGAGCTGGTGGTGTGTGGGCTGCACCCGGCGCCCGTGCGGGGCCCCGCCGTGCGTGC  
20 AGGATTTAGGGGCGCCGCTGACCTCCCCGCAGCCCTGGCCGGGATTCCTGGGCCAAGGGG  
AAGCGCTGCTGTCCCTGGCCATGCTGCTGCGTCTCTACCTGGTGCCCCGCGCCGTGCTCCT  
GCGCAGCGGCGTCTGCTCAACGCTTCTACCGCAGCATCGGCGCTCTCAATCAAGTCCG  
CTTCCGCCACTGGTTCGTGGCCAAGCTTTACATGAACACGCACCCTGGCCGCCTGCTGCTC  
GGCCTCACGCTTGGCCTCTGGCTGACCACCGCCTGGGTGCTGTCCGTGGCCGAGAGGCAG  
25 GCTGTTAATGCCACTGGGCACCTTTTCAGACACACTTTGGCTGATCCCCATCACATTCTGA  
CCATCGGCTATGGTGACGTGGTGCCGGGCACCATGTGGGGCAAGATCGTCTGCCTGTGCA  
CTGGAGTCATGGGTGTCTGCTGCACAGCCCTGCTGGTGGCCGTGGTGGCCCCGGAAGCTGG  
AGTTTAACAAGGCAGAGAAGCACGTGCACAACTTCATGATGGATATCCAGTATACCAAAG  
AGATGAAGGAGTCCGCTGCCCCAGTGCTACAAGAAGCCTGGATGTTCTACAAACATACTC  
30 GCAGGAAGGAGTCTCATGCTGCCCCGAGGCATCAGCGCAAGCTGCTGGCCGCCATCAACG  
CGTTCCGCCAGGTGCGGCTGAAACACCGGAAGCTCCGGGAACAAGTGAAGTCCATGGTG  
GACATCTCCAAGATGCACATGATCCTGTATGACCTGCAGCAGAATCTGAGCAGCTCACAC  
CGGGCCCTGGAGAAACAGATTGACACGCTGGCGGGGAAGCTGGATGCCCTGACTGAGCT  
GCTTAGCACTGCCCTGGGGCCGAGGCAGCTTCCAGAACCCAGCCAGCAGTCCAAGTAGCT  
35 GGACCCACGAGGAGGAACCAGGCTACTTTCCCCAGTACTGAGGTGGTGGACATCG

2/29

FIG. 1 (Continued)

5 TCTCTGCCACTCCTGACCCAGCCCTGAACAAAGCACCTCAAGTGCAAGGACCAAAGGGGG  
CCCTGGCTTGGAGTGGGTTGGCTTGCTGATGGCTGCTGGAGGGGACGCTGGCTAAAGTG  
GGTAGGCCTTGGCCACCTGAGGCCCCAGGTGGGAACATGGTCACCCCCACTCTGCATAC  
CCTCATCAAAAACACTCTCACTATGCTGCTATGGACGACCTCCAGCTCTCAGTTACAAGTG  
CAGGCGACTGGAGGCAGGACTCCTGGGTCCCTGGGAAAGAGGGTAATAGGGGCCCCGGAT  
10 XCAGGATTCTGGGAGGCTTCAGTTACCGCTGGCCGAGCTGAAGAACTGGGTATGAGGCTG  
GGGCGGGGCTGGAGGTGGCGCCCCCTGGTGGGACAACAAAGAGGACACCATTTTTCCAG  
AGCTGCAGAGAGCAXCTGGTGGGGAGGAAGAAGTGTAACCTACCAGCCTTTGCTCTTATT  
TTTGTAATAAATGTTAAAGCCAGAAAAAAAAAAAAAAAAAAAAAAAAAAAA 3'

15

SEQ ID NO:1

20

3/29

FIG.2

5 ATGGGCGGGGATCTGGTGCTTGGCCTGGGGGCCTTGAGACGCCGAAAGCGCTTGCTGGA  
GCAGGAGAAGTCTCTGGCCGGCTGGGCACTGGTGCTGGCAGGAAGTGGCATTGGACTCAT  
GGTGCTGCATGCAGAGATGCTGTGGTTTCGGGGGGTGCTCGTGGGCGCTCTACCTGTTCT  
GGTTAAATGCACGATCAGCATTTCCACCTTCTTACTCCTCTGCCTCATCGTGGCCTTTTCAT  
GCCAAAGAGGTCCAGCTGTTTCATGACCGACAACGGGCTGCGGGACTGGCGCGTGGCGCT  
10 GACCGGGCGGCAGGCGGCGCAGATCGTGCTGGAGCTGGTGGTGTGTGGGCTGCACCCGG  
CGCCCGTGCGGGGCCCCGCGTGCGTGCAAGATTTAGGGGCGCCGCTGACCTCCCCGCAGC  
CCTGGCCGGGATTCTGGGCCAAGGGGAAGCGCTGCTGTCCCTGGCCATGCTGCTGCGTC  
TCTACCTGGTGCCCCGCGCCGTGCTCCTGCGCAGCGGCGTCTGCTCAACGCTTCTTACCG  
CAGCATCGGCGCTCTCAATCAAGTCCGCTTCCGCCACTGGTTCGTGGCCAAGCTTTACATG  
15 AACACGCACCCTGGCCGCCTGCTGCTCGGCCTCACGCTTGGCCTCTGGCTGACCACCGCCT  
GGGTGCTGTCCGTGGCCGAGAGGCAGGCTGTTAATGCCACTGGGCACCTTTCAGACACAC  
TTTGGCTGATCCCCATCACATTCTGACCATCGGCTATGGTGACGTGGTGGCGGGCACCAT  
GTGGGGCAAGATCGTCTGCCTGTGCACTGGAGTCATGGGTGTCTGCTGCACAGCCCTGCT  
GGTGGCCGTGGTGGCCCGGAAGCTGGAGTTTAACAAGGCAGAGAAGCACGTGCACAAC  
20 TCATGATGGATATCCAGTATACCAAAGAGATGAAGGAGTCCGCTGCCCCGAGTGCTACAAG  
AAGCCTGGATGTTCTACAAACATACTCGCAGGAAGGAGTCTCATGCTGCCCCGAGGCATC  
AGCGCAAGCTGCTGGCCGCCATCAACGCGTTCCGCCAGGTGCGGCTGAAACACCGGAAG  
CTCCGGGAACAAGTGAACCTCCATGGTGGACATCTCCAAGATGCACATGATCCTGTATGAC  
CTGCAGCAGAATCTGAGCAGCTCACACCGGGCCCTGGAGAAACAGATTGACACGCTGGC  
25 GGGGAAGCTGGATGCCCTGACTGAGCTGCTTAGCACTGCCCTGGGGCCGAGGCAGCTTCC  
AGAACCCAGCCAGCAGTCCAAGTAG

30

SEQ ID NO:2

4/29

## FIG.3

5

MGGDLVLGLGALRRRKRLLEQEKSLAGWALVLAGTGIGLMVLHAEMLWFGGCSW  
ALYLFLVKCTISISTFLLLCLIVAFHAKEVQLFMTDNGLRDWRVALTGRQAAQIVLEL  
VVCGLHPAPVRGPPCVQDLGAPLTSPQPWPGFLGQGEALLSLAMLLRLYLVPRAVLL  
RSGVLLNASYRSIGALNQVRFRHWFVAKLYMNTHPGRLLLGLTLGLWLTAWVLSV  
10 AERQAVNATGHLSDTLWLIPITFLTIGYGDVVPGTMWGKIVCLCTGVMGVCCTALLV  
AVVARKLEFNKAEKHVHNFMMDIQYTKEMKESAARVLQEAWMFYKHTRRKESHA  
ARRHQRKLLAAINAFRQVRLKHRKLREQVNSMVDISKMHMILYDLQQNLSSSHRAL  
EKQIDTLAGKLDALTELLSTALGPRQLPEPSQQSK\*

15

SEQ ID NO:3

20

5/29

FIG.4

5  
MPGPRAACSEPNPCTQVVMNSHSYNGSVGRPLGSGPGALGRDPPDPEAGHPPQPPHS  
PGLQVVVAKSEPARPSPGSPRGQPQDQDDDEDDEEDEAGRQRASGKPSNVGHR LGH  
RRALFEKRKRLSDYALIFGMFGIVVMVTETELSWG VYTKESLYSFALKCLISLSTAILL  
GLVVLYHAREIQLFMVDNGADDWRIAMTCERVFLISLELAVCAIHPVPGHYRFTWTA  
10 RLAF TYAPSVAEADV DVLLSIPMFLRLYLLGRV MLLH SKIFTDASSRSIGALNKITFNT  
RFVMKTLMTICPGTVLLVFSISSWIIAAWTVRVCERYHDKQEVT SNFLGAMWLISITFL  
SIGYGDMVPHTYCGKGVCLLTGIMGAGCTALVVAVVARKLELTKA EKHVHNFMMMD  
TQLTKRVKNAAANVLRETWLIYKHTRLVKKPDQARVRKHQRKFLQAIHQAQKLRSV  
KIEQGKLNDQANTLTDLAKTQTVMYDLVSELHAQHEELEARLATLESRLDALGASLQ  
15 ALPGLIAQAIRPPPPPLPPRPGPGPQDQAARSSPCRWTPVAPSDCG

20

SEQ ID NO:4

6/29

FIG.5

5 MSSCRYNGGVMRPLSNLSSSRRLHEMDSEAQPLQPPASVVGGGGGASSPSAAAAA  
SSSAPEIVVSKPEHNNSNNLALYGTGGGGSTGGGGGGGGGGGGSGHGSSSGTKSSKK  
KNQNIGYKLGHRRALFKRKRLSDYALIFGMFGIVVMVIETELSWGAYDKASLYSLAL  
KCLISLSTILLGLIIVYHAREIQLFMVDNGADDWRIAMTYERIFFICLEILVCAIHPIPGN  
YTFTWTARLAFSYAPSTTTADVDIILSIPMFLRLYLIARVMLLHSLFTDASSRSIGALN  
10 KINFNTRFVMKTLMTICPGTVLLVFSISLWIIAAWTVRACERYHDQQDVTSNFLGAM  
WLISITFLSIGYGDMVPNTYCGKGVCLLTGIMGAGCTALVVAVVARKLELTKAEKHV  
HNFMMDTQLTKRVKNAAANVLRETWLIYKNTKLVKKIDHAKVRKHQRKFLQAIHQL  
RSVKMEQRKLNDQANTLVDLAKTQNIMYDMISDLNERSEDFEKRIVTLETKLETIGS  
IHALPGLISQTIQQQQRDFIETQMENYDKHVTYNAERSRSSSRRRRSSSTAPPTSSESS

15

SEQ ID NO:5

20



7/29

FIG.6

5  
MSSCKYSGGVMKPLSRLSASRRNLIEAEPEGQPLQLFSPSNPPEIISSREDNHAHQTL  
HHPNATHNHQHAGTTAGSTTFPKANKRKNQNIGYKLGHRRALFEKRKRLSDYALIFG  
MFGIVVMVIETELSWGLYSKDSMFSLALKCLISLSTIILLGLIIAYHTREVQLFVIDNGA  
DDWRIAMTYERILYISLEMLVCAIHPIPGYKFFWTARLAFSYTPSRAEADVDIILSIPM  
10 FLRLYLIARVMLLHSLFTDASSRSIGALNKINTRFVMKTLMTICPGTVLLMFSISLWII  
AAWTVRVCERYHDQQDVTSNFLGAMWLISITFLSIGYGDMVPHTYCGKGVCLLTGI  
MGAGCTALVVAVVARKLELTKAEKHVHNFMMDTQLTKRIKNAAANVLRETWLIYK  
HTKLLKKIDHAKVRKHQRKFLQAIHQLRGVKMEQRKLSAQANTLVDLSKMQNVMY  
DLITELNDRSEDLEKQIGSLESKLEHLTASFNSLPLLIADTLRQQQQQLLTAFVEARGIS  
15 VAVGTSHAPPSDSPIGISSTSFPEFLIF

20

SEQ ID NO:6

8/29

FIG.7

M S G C R Y X G G V M X P L S X L S X S R R N L X E - X S E Majority																														
10										20										30										
1	M	G	G	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
1	M	P	G	P	R	A	A	C	S	E	P	N	P	C	T	Q	V	V	M	N	S	H	S	Y	N	-	-	G	S	V
1	M	S	S	C	R	Y	N	G	G	V	M	R	P	L	S	N	L	S	S	S	R	R	N	L	H	E	M	D	S	E
1	M	S	S	C	K	Y	S	G	G	V	M	K	P	L	S	R	L	S	A	S	R	R	N	L	I	E	A	E	P	E
G Q P L Q - - X X V L G X G X - - - - - X X X P S S P Majority																														
40										50										60										
4	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
29	G	R	P	L	G	S	G	P	G	A	L	G	R	D	P	P	D	P	E	A	G	H	P	P	Q	P	P	H	S	P
31	A	Q	P	L	Q	P	P	A	S	V	V	G	G	G	A	S	S	P	S	A	A	A	A	A	S	S	S	A	A	A
31	G	Q	P	L	Q	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
P - E I V V S K X E X N X X X X - - L - G X X X G X X X X X Majority																														
70										80										90										
12	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
59	G	L	Q	V	V	V	A	K	S	E	P	A	R	P	S	P	-	-	-	-	G	S	P	R	G	Q	P	Q	D	Q
61	P	-	E	I	V	V	S	K	P	E	H	N	N	S	N	N	L	A	L	Y	G	T	G	G	G	S	T	G	G	G
43	P	-	E	I	I	I	S	S	R	E	D	N	H	A	H	Q	T	L	L	H	H	P	N	A	T	H	N	H	Q	H
X G X X X G - - - - - G - - S X S G - K - - K - K N Q N I G Majority																														
100										110										120										
12	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
85	D	D	D	E	D	D	E	E	D	E	A	G	R	Q	R	A	S	G	-	-	-	-	-	-	-	K	P	S	N	V
90	G	G	G	G	G	G	G	G	S	G	H	G	S	S	S	G	T	K	S	S	K	K	K	N	Q	N	I	G	I	G
72	A	G	T	T	A	G	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
Y K L G H R R A L F E K R K R L S D Y A L I F G M F G I V V Majority																														
130										140										150										
12	-	-	L	R	R	R	K	R	L	L	E	Q	E	K	S	L	A	G	W	A	L	V	L	A	G	T	G	I	G	L
109	H	R	L	G	H	R	R	A	L	F	E	K	R	K	R	L	S	D	Y	A	L	I	F	G	M	F	G	I	V	V
120	Y	K	L	G	H	R	R	A	L	F	E	K	R	K	R	L	S	D	Y	A	L	I	F	G	M	F	G	I	V	V
94	Y	K	L	G	H	R	R	A	L	F	E	K	R	K	R	L	S	D	Y	A	L	I	F	G	M	F	G	I	V	V

HKCA4.PRO = SEQ ID NO:3  
 HSK1.PRO = SEQ ID NO:4  
 RSK2.PRO = SEQ ID NO:5  
 RSK3.PRO = SEQ ID NO:6

9/29

FIG.7 CONT.D

M V I E T E L S W G X Y X K X S L Y S L A L K C L I S L S T																												Majority			
160														170														180			
40	M	V	L	H	A	E	M	L	W	F	G	G	C	S	W	A	L	Y	L	F	L	V	K	C	T	I	S	I	S	T	HKCA4.PRO
139	M	V	T	E	T	E	L	S	W	G	V	Y	T	K	E	S	L	Y	S	F	A	L	K	C	L	I	S	L	S	T	HSK1.PRO
150	M	V	I	E	T	E	L	S	W	G	A	Y	D	K	A	S	L	Y	S	L	A	L	K	C	L	I	S	L	S	T	RSK2.PRO
124	M	V	I	E	T	E	L	S	W	G	L	Y	S	K	D	S	M	F	S	L	A	L	K	C	L	I	S	L	S	T	RSK3.PRO
I I L L G L I V A Y H A R E V Q L F M V D N G A D D W R I A																												Majority			
190														200														210			
70	F	L	L	L	C	L	I	V	A	F	H	A	K	E	V	Q	L	F	M	T	D	N	G	L	R	D	W	R	V	A	HKCA4.PRO
169	A	I	L	L	G	L	V	V	L	Y	H	A	R	E	I	Q	L	F	M	V	D	N	G	A	D	D	W	R	I	A	HSK1.PRO
180	I	I	L	L	G	L	I	I	V	Y	H	A	R	E	I	Q	L	F	M	V	D	N	G	A	D	D	W	R	I	A	RSK2.PRO
154	I	I	L	L	G	L	I	I	A	Y	H	T	R	E	V	Q	L	F	V	I	D	N	G	A	D	D	W	R	I	A	RSK3.PRO
M T Y E R I F X I S L E L L V C A I H P I P G X - - - Y X F																												Majority			
220														230														240			
100	L	T	G	R	Q	A	A	Q	I	V	L	E	L	V	V	C	G	L	H	P	A	P	V	R	G	P	P	C	V	O	HKCA4.PRO
199	M	T	C	E	R	V	F	L	I	S	L	E	L	A	V	C	A	I	H	P	V	P	G	H	-	-	-	Y	R	F	HSK1.PRO
210	M	T	Y	E	R	I	F	F	I	C	L	E	I	L	V	C	A	I	H	P	I	P	G	N	-	-	-	Y	T	F	RSK2.PRO
184	M	T	Y	E	R	I	L	Y	I	S	L	E	M	L	V	C	A	I	H	P	I	P	G	E	-	-	-	Y	K	F	RSK3.PRO
T W T A R L A F S Y A P S X A E A D V D I L L S I P M F L R																												Majority			
250														260														270			
130	D	L	G	A	P	L	T	S	P	Q	P	W	P	G	F	L	G	Q	G	E	A	L	L	S	L	A	M	L	L	R	HKCA4.PRO
226	T	W	T	A	R	L	A	F	T	Y	A	P	S	V	A	E	A	D	V	D	V	L	L	S	I	P	M	F	L	R	HSK1.PRO
237	T	W	T	A	R	L	A	F	S	Y	A	P	S	T	T	T	A	D	V	D	I	I	L	S	I	P	M	F	L	R	RSK2.PRO
211	F	W	T	A	R	L	A	F	S	Y	T	P	S	R	A	E	A	D	V	D	I	I	L	S	I	P	M	F	L	R	RSK3.PRO
L Y L I A R V M L L H S K L F T D A S S R S I G A L N K I N																												Majority			
280														290														300			
160	L	Y	L	V	P	R	A	V	L	L	R	S	G	V	L	L	N	A	S	Y	R	S	I	G	A	L	N	Q	V	R	HKCA4.PRO
256	L	Y	L	L	G	R	V	M	L	L	H	S	K	I	F	T	D	A	S	S	R	S	I	G	A	L	N	K	I	T	HSK1.PRO
267	L	Y	L	I	A	R	V	M	L	L	H	S	K	L	F	T	D	A	S	S	R	S	I	G	A	L	N	K	I	N	RSK2.PRO
241	L	Y	L	I	A	R	V	M	L	L	H	S	K	L	F	T	D	A	S	S	R	S	I	G	A	L	N	K	I	N	RSK3.PRO

HKCA4.PRO = SEQ ID NO:3  
 HSK1.PRO = SEQ ID NO:4  
 RSK2.PRO = SEQ ID NO:5  
 RSK3.PRO = SEQ ID NO:6

10/29

FIG.7 CONT.D

F N T R F V M K T L M T I C P G T V L L V F S I S L W I I A Majority																																					
310														320														330									
190	F	R	H	W	F	V	A	K	L	Y	M	N	T	H	P	G	R	L	L	L	G	I	T	L	G	L	W	L	T	T	HKCA4.PRO						
286	F	N	T	R	F	V	M	K	T	L	M	T	I	C	P	G	T	V	L	L	V	F	S	I	S	S	W	I	I	A	HSK1.PRO						
297	F	N	T	R	F	V	M	K	T	L	M	T	I	C	P	G	T	V	L	L	V	F	S	I	S	L	W	I	I	A	RSK2.PRO						
271	F	N	T	R	F	V	M	K	T	L	M	T	I	C	P	G	T	V	L	L	M	F	S	I	S	L	W	I	I	A	RSK3.PRO						
A W T V R V C E R Y H D Q Q D V T S N F L G A M W L I S I T Majority																																					
340														350														360									
220	A	W	V	L	S	V	A	E	R	-	-	Q	A	V	N	A	T	G	H	L	S	D	T	L	W	L	I	P	I	T	HKCA4.PRO						
116	A	W	T	V	R	V	C	E	R	Y	H	D	K	Q	E	V	T	S	N	F	L	G	A	M	W	L	I	S	I	T	HSK1.PRO						
127	A	W	T	V	R	A	C	E	R	Y	H	D	Q	Q	D	V	T	S	N	F	L	G	A	M	W	L	I	S	I	T	RSK2.PRO						
101	A	W	T	V	R	V	C	E	R	Y	H	D	Q	Q	D	V	T	S	N	F	L	G	A	M	W	L	I	S	I	T	RSK3.PRO						
F L S I G Y G D M V P H T Y C G K G V C L L T G I M G A G C Majority																																					
370														380														390									
48	F	L	T	I	G	Y	G	D	V	V	P	G	T	M	W	G	K	I	V	C	L	C	T	G	V	M	G	V	C	C	HKCA4.PRO						
46	F	L	S	I	G	Y	G	D	M	V	P	H	T	Y	C	G	K	G	V	C	L	L	T	G	I	M	G	A	G	C	HSK1.PRO						
57	F	L	S	I	G	Y	G	D	M	V	P	N	T	Y	C	G	K	G	V	C	L	L	T	G	I	M	G	A	G	C	RSK2.PRO						
31	F	L	S	I	G	Y	G	D	M	V	P	H	T	Y	C	G	K	G	V	C	L	L	T	G	I	M	G	A	G	C	RSK3.PRO						
T A L V V A V V A R K L E L T K A E K H V H N F M M D T Q L Majority																																					
400														410														420									
78	T	A	L	L	V	A	V	V	A	R	K	L	E	F	N	K	A	E	K	H	V	H	N	F	M	M	D	I	Q	Y	HKCA4.PRO						
76	T	A	L	V	V	A	V	V	A	R	K	L	E	L	T	K	A	E	K	H	V	H	N	F	M	M	D	T	Q	L	HSK1.PRO						
87	T	A	L	V	V	A	V	V	A	R	K	L	E	L	T	K	A	E	K	H	V	H	N	F	M	M	D	T	Q	L	RSK2.PRO						
61	T	A	L	V	V	A	V	V	A	R	K	L	E	L	T	K	A	E	K	H	V	H	N	F	M	M	D	T	Q	L	RSK3.PRO						
T K R V K N A A A N V L R E T W L I Y K H T K L V K K I D H Majority																																					
430														440														450									
08	T	K	E	M	K	E	S	A	A	R	V	L	Q	E	A	W	M	F	Y	K	H	T	R	-	-	R	K	E	S	H	HKCA4.PRO						
06	T	K	R	V	K	N	A	A	A	N	V	L	R	E	T	W	L	I	Y	K	H	T	R	L	V	K	K	P	D	Q	HSK1.PRO						
17	T	K	R	V	K	N	A	A	A	N	V	L	R	E	T	W	L	I	Y	K	N	T	K	L	V	K	K	I	D	H	RSK2.PRO						
91	T	K	R	I	K	N	A	A	A	N	V	L	R	E	T	W	L	I	Y	K	H	T	K	L	L	K	K	I	D	H	RSK3.PRO						

HKCA4.PRO = SEQ ID NO:3  
 HSK1.PRO = SEQ ID NO:4  
 RSK2.PRO = SEQ ID NO:5  
 RSK3.PRO = SEQ ID NO:6

11/29

FIG.7 CONT.D


AKVRKHQRKFLQAIHQ - - - LRSVKMEQRKL																				Majority											
460										470										480											
336	A	-	A	R	R	H	Q	R	K	L	L	A	A	I	N	A	-	-	-	F	R	Q	V	R	L	K	H	R	K	L	HKCA4.PRO
436	A	R	V	R	K	H	Q	R	K	F	L	Q	A	I	H	Q	A	Q	K	L	R	S	V	K	I	E	Q	G	K	L	HSK1.PRO
447	A	K	V	R	K	H	Q	R	K	F	L	Q	A	I	H	Q	-	-	-	L	R	S	V	K	M	E	Q	R	K	L	RSK2.PRO
421	A	K	V	R	K	H	Q	R	K	F	L	Q	A	I	H	Q	-	-	-	L	R	G	V	K	M	E	Q	R	K	L	RSK3.PRO
NDQANTLVDLAKTQNVMYDLISELNXRSED																				Majority											
490										500										510											
362	R	E	Q	V	N	S	M	V	D	I	S	K	M	H	M	I	L	Y	D	L	Q	Q	N	L	S	S	S	H	R	A	HKCA4.PRO
466	N	D	Q	A	N	T	L	T	D	L	A	K	T	Q	T	V	M	Y	D	L	V	S	E	L	H	A	Q	H	E	E	HSK1.PRO
474	N	D	Q	A	N	T	L	V	D	L	A	K	T	Q	N	I	M	Y	D	M	I	S	D	L	N	E	R	S	E	D	RSK2.PRO
448	S	D	Q	A	N	T	L	V	D	L	S	K	M	Q	N	V	M	Y	D	L	I	T	E	L	N	D	R	S	E	D	RSK3.PRO
LEKQIXTLESKLDALTASXXALPGLIAQAL																				Majority											
520										530										540											
392	L	E	K	Q	I	D	T	L	A	G	K	L	D	A	L	T	E	-	-	-	-	-	-	L	L	S	T	A	L	HKCA4.PRO	
496	L	E	A	R	L	A	T	L	E	S	R	L	D	A	L	G	A	S	L	Q	A	L	P	G	L	I	A	Q	A	I	HSK1.PRO
504	F	E	K	R	I	V	T	L	E	T	K	L	E	T	L	I	G	S	I	H	A	L	P	G	L	I	S	Q	T	I	RSK2.PRO
478	L	E	K	Q	I	G	S	L	E	S	K	L	E	H	L	T	A	S	F	N	S	L	P	L	L	I	A	D	T	L	RSK3.PRO
RQQQXXLXX - - - E - - - - X V X X G X X X X X S S																				Majority											
550										560										570											
415	G	P	R	Q	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	HKCA4.PRO
526	R	P	P	P	P	P	L	P	P	-	-	-	-	-	-	-	-	R	P	G	P	G	P	Q	D	Q	A	A	R	S	HSK1.PRO
534	R	Q	Q	Q	R	D	F	I	E	T	Q	M	E	N	Y	D	K	H	V	T	Y	N	A	E	R	S	R	S	S	S	RSK2.PRO
508	R	Q	Q	Q	Q	Q	L	L	T	A	F	V	E	A	R	G	I	S	V	A	V	G	T	S	H	A	P	P	S	D	RSK3.PRO
SPXRXSSTAPPEXSXS -																				Majority											
580																															
419	-	-	-	-	-	-	-	-	-	L	P	E	P	S	Q	Q	S	K	-	-	-	-	-	-	-	-	-	-	-	-	HKCA4.PRO
548	S	P	C	R	W	T	P	V	A	P	S	-	-	-	D	C	G	-	-	-	-	-	-	-	-	-	-	-	-	-	HSK1.PRO
564	R	R	R	R	S	S	S	T	A	P	P	T	S	S	E	S	S	-	-	-	-	-	-	-	-	-	-	-	-	-	RSK2.PRO
538	S	P	I	G	I	S	S	T	S	F	P	E	F	L	-	I	F	-	-	-	-	-	-	-	-	-	-	-	-	-	RSK3.PRO

Decoration 'Decoration #1': Box residues that match the Consensus exactly.

HKCA4.PRO = SEQ ID NO:3  
 HSK1.PRO = SEQ ID NO:4  
 RSK2.PRO = SEQ ID NO:5  
 RSK3.PRO = SEQ ID NO:6

12/29

FIG.8




P = Pore Region

13/29

FIG.9

Northern Blot analysis  
<sup>32</sup>P labeled, overnight exposure



M=Molecular weight markers  
kb=kilobases

BEST AVAILABLE COPY

14/29

FIG.10

5  
CAGACCTACTGAATCCTGGAGGGTAGGGCCAGCAATCCACTTTAGCAAGCCTTCC  
AGGCAGGCTGTTAATGCCACTGGGCACCTTTCAGACACACTTTGGCTGATCCCCA  
TCACATTCTGACCATCGGCTATGGTGACGTGGTGCCGGGCACCATGTGGGGCAA  
GATCGTCTGCCTGTGCACTGGAGTCATGGGTGTCTGCTGCACAGCCCTGCTGGTG  
10 GCCGTGGTGGCCCGGAAGCTGGAGTTTAACAAGGCAGAGAAGCACGTGCACAAC  
TTCATGATGGATATCCAGTATACCAAAGAGATGAAGGAGTCCGCTGCCCCGAGTGC  
TACAAGAAGCCTGGATGTTCTACAAACATACTCGCAGGAAGGAGTCTCATGCTGC  
CCGCAGGCATCAGCGCAAGCTGCTGGCCGCCATCAACGCGTTCCGCCAGGTGCGG  
CTGAAACACCGGAAGCTCCGGGAACAAGTGAAGTCCATGGTGGACATCTCCAAG  
15 ATGCACATGATCCTGTATGACCTGCAGCAGAATCTGAGCAGCTCACACCGGGCCC  
TGGAGAAACAGATTGACACGCTGGCGGGGAAGCTGGATGCCCTGACTGAGCTGC  
TTAGCACTGCCCTGGGGCCGAGGCAGCTTCCAGAACCCAGCCAGCAGTCCAAGTA  
GCTGGACCCACGAGGAGGAACCAGGCTACTTTCCCCAGTACTGAGGTGGTGGAC  
ATCGTCTCTGCCACTCCTGACCCAGCCCTGAACAAAGCACCTCAAGTGCAAGGAC  
20 CAAAGGGGGCCCTGGCTTGGAGTGGGTTGGCTTGCTGATGGCTGCTGGAGGGGAC  
GCTGGCTAAAGTGGGTAGGCCTTGGCCACCTGAGGCCCCAGGTGGGAACATGGT  
CACCCCCACTCTGCATACCCTCATCAAAAACACTCTCACTATGCTGCTATGGACG  
ACCTCCAGCTCTCAGTTACAAGTGCAGGCGACTGGAGGCAGGACTCCTGGGTCCC  
TGGGAAAGAGGGTACTAGGGGGCCCGATCCAGGATTCTGGGAGGCTTCAGTTAC  
25 CGCTGGCCGAGCTGAAGAACTGGGTATGAGGCTGGGGCGGGGCTGGAGGTGGCG  
CCCCCTGGTGGGACAACAAAGAGGACACCATTTTTCCAGAGCTGCAGAGAGCAA  
CTGGTGGGGAGGAAGAAGTGTAATCACCAGCCTTTGCTTTTATTTTTGTAATAA  
ATGTTAAAGCCAGXAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

30

SEQ ID NO:8



15/29

FIG.11

5 CAAAGAGGTCCAGCTGTTTCATGACCGACAACGGGCTGCGGGACTGGCGCGTGGC  
GCTGACCGGGCGGCAGGCGGCGCAGATCGTGCTGGAGCTGGTGGTGTGTGGGCT  
GCACCCGGCGCCCGTGCGGGGCCCCGCCGTGCGTGCAGGATTTAGGGGCGCCGCT  
GACCTCCCCGCAGCCCTGGCCGGGATTCTGGGCCAAGGGGAAGCGCTGCTGTCC  
CTGGCCATGCTGCTGCGTCTCTACCTGGTGCCCCGCGCCGTGCTCCTGCGCAGCG  
10 GCGTCCTGCTCAACGCTTCCTACCGCAGCATCGGCGCTCTCAATCAAGTCCGCTTC  
CGCCACTGGTTTCGTGGCCAAGCTTTACATGAACACGCACCCTGGCCGCCTGCTGC  
TCGGCCTCACGCTTGGCCTCTGGCTGACCACCGCCTGGGTGCTGTCCGTGGCCGA  
GAGGCAGGCTGTTAATGCCACTGGGCACCTTTCAGACACACTTTGGCTGATCCCC  
ATCACATTCCTGACCATCGGCTATGGTGACGTGGTGCCGGGCACCATGTGGGGCA  
15 AGATCGTCTGCCTGTGCACTGGAGTCATGGGTGTCTGCTGCACAGCCCTGCTGGT  
GGCCGTGGTGGCCCGGAAGCTGGAGTTTAACAAGGCAGAGAAGCACGTGCACAA  
CTTCATGATGGATATCCAGTATACCAAAGAGATGAAGGAGTCCGCTGCCCCGAGTG  
CTACAAGAAGCCTGGATGTTCTACAAACATACTCGCAGGAAGGAGTCTCATGCTG  
CCCGCAGGCATCAGCGCAAGCTGCTGGCCGCCATCAACGCGTTCCGCCAGGTGCG  
20 GCTGAAACACCGGAAGCTCCGGGAACAAGTGAAGTCCATGGTGGACATCTCCAA  
GATGCACATGATCCTGTATGACCTGCAGCAGAATCTGAGCAGCTCACACCGGGGCC  
CTGGAGAAACAGATTGACACGCTGGCGGGGAAGCTGGATGCCCTGACTGAGCTG  
CTTAGCACTGCCCTGGGGCCGAGGCAGCTTCCAGAAGCCAGCCAGCAGTCCAAGT  
AGCTGGACCCACGAGGAGGAACCAGGCTACTTTCCCCAGTACTGAGGTGGTGGA  
25 CATCGTCTCTGCCACTCCTGACCCAGCCCTGAACAAAGCACCTCAAGTGCAAGGA  
CCAAAGGGGGCCCTGGCTTGGAGTGGGTTGGCTTGCTGATGGCTGCTGGAGGGGA  
CGCTGGCTAAAGTGGGTAGGCCTTGGCCACCTGAGGCCCCAGGTGGGAACATG  
GTCACCCCCACTCTGCATACCCTCATCAAAAACACTCTCACTATGCTGCTATGGAC  
GACCTCCAGCTCTCAGTTACAAGTGCAGGCGACTGGAGGCAGGACTCCTGGGTCC  
30 CTGGGAAAGAGGGTAATAGGGGCCCCGGATXCAGGATTCTGGGAGGCTTCAGTTA  
CCGCTGGCCGAGCTGAAGAACTGGGTATGAGGCTGGGGCGGGGCTGGAGGTGGC  
GCCCCCTGGTGGGACAACAAAGAGGACACCATTTTTCCAGAGCTGCAGAGAGCA

16/29

**FIG. 11 CONT.D**

5 XCTGGTGGGGAGGAAGAAGTGTAACCTCACCAGCCTTTGCTCTTATTTTGTAAATA  
AATGTTAAAGCCAGAAAAAAAAAAAAAAAAAAAAAAAAA

10

SEQ ID NO:9

17/29

FIG. 12

5  
5'  
GTCCTTCGGTGTCTGGGTGTGGTGAGTAGAGGTGTGTGTACAAAGTACAGACCATTGTGTGTGACAAAGCCCAT  
CGTGTGTCTGTGTGTGTCTTTATCCACGTGGATGGACGTCTCTTTCTTGCTCTGCCCCAAGACACACCCTAGCCC  
CTCCTTATTCTCAAAAGGGGGAGCTGGGGAGCCTCCCCCTACCCCTGGGGCCTCCCCCTGCCCTCCCCGCCCTGCC  
10 TGGCCGTCACCACTCCCCAGAGGGCACAGGGCTCTGCTGTGCCTCAGAGCAAAAGTCCCAGAGCCAGCAGAGCAG  
GCTGACGACCTGCAAGCCACAGTGGCTGCCCTGTGCGTGTGCGAGGTGGGGGACCCTGGGCAGGAAGCTGGCTG  
AGCCCCAAGACCCCGGGGGCCATGGGCGGGGATCTGGTGCCTTGGCCCTGGGGGCCCTTGAGACGCCGAAAGCGCTTG  
CTGGAGCAGGAGAAGTCTCTGGCCGGCTGGGCACTGGTGTGCGCAGGAAGTGGCATTGGACTCATGGTGTGCTGCAT  
GCAGAGATGTGTGGTTTCGGGGGGTGTCTCGTGGGCGCTCTACCTGTTTCTGGTTAAATGCACGATCAGCATTTTCC  
15 ACCTTCTTACTCCTCTGCCTCATCGTGGCCTTTTCATGCCAAAGAGGTCCAGCTGTTTCATGACCGACAACGGGCTG  
CGGACTGGCGCTGGCGCTGACCGGGCGGCAGGCGGCGCAGATCGTGCTGGAGCTGGTGGTGTGTGGGCTGCAC  
CCGGCGCCCCGTGGGGGCCCCCGTGCCTGCAGGATTTAGGGGCGCCGCTGACCTCCCCGAGCCCTGGCCGGGA  
TTCCTGGGCCAAGGGGAAGCGCTGCTGTCCCTGGCCATGCTGCTGCGTCTCTACCTGGTGCCCCGCGCCGTGCTC  
CTGCGCAGCGCGCTCTGCTCAACGCTTCCTACCGCAGCATCGGCGCTCTCAATCAAGTCCGCTTCCGCCACTGG  
20 TTCGTGGCCAAGCTTTACATGAACACGCACCCTGGCCGCCTGCTGCTCGGCCTCACGCTTGGCCTCTGGCTGACC  
ACCGCCTGGGTGCTGTCCGTGGCCGAGAGGCAGGCTGTTAATGCCACTGGGCACCTTTTCAGACACACTTTGGCTG  
ATCCCCATCACATTCTGACCATCGGCTATGGTGACGTGGTGCCGGGCACCATGTGGGGCAAGATCGTCTGCCTG  
TGCACTGGAGTCATGGGTGTCTGCTGCACAGCCCTGCTGGTGGCCGTGGTGGCCCGGAAGCTGGAGTTTAAACAAG  
GCAGAGAAGCACGTGCACAACCTTCATGATGGATATCCAGTATACCAAAGAGATGAAGGAGTCCGCTGCCCGAGTG  
25 CTACAAGAAGCCTGGATGTTCTACAAACATACTCGCAGGAAGGAGTCTCATGCTGCCCGCAGGCATCAGCGCAAG  
CTGCTGGCCGCCATCAACGCGTTCCGCCAGGTGCGGCTGAAACACCGGAAGCTCCGGGAACAAGTGAACCTCCATG  
GTGGACATCTCCAAGATGCACATGATCCTGTATGACCTGCAGCAGAATCTGAGCAGCTCACACCGGGCCCTGGAG  
AAACAGATTGACACGCTGGCGGGGAAGCTGGATGCCCTGACTGAGCTGCTTAGCACTGCCCTGGGGCCGAGGCAG  
CTTCCAGAACCCAGCCAGCAGTCCAAGTAGCTGGACCCACGAGGAGGAACCAGGCTACTTTCCCCAGTACTGAGG  
30 TGGTGGACATCGTCTCTGCCACTCCTGACCCAGCCCTGAACAAAGCACCTCAAGTGCAAGGACCAAGGGGGCCC  
TGGCTTGGAGTGGGTGGCTTGGCTGATGGCTGCTGGAGGGGACGCTGGCTAAAGTGGGTAGGCCTTGGCCACCT  
GAGGCCCCAGGTGGGAACATGGTCACCCCCACTCTGCATACCCCTCATCAAAAACACTCTCACTATGCTGCTATGG  
ACGACCTCCAGCTCTCAGTTACAAGTGCAGGCGACTGGAGGCAGGACTCCTGGGTCCCTGGGAAAGAGGGTACTA  
GGGGCCCCGATCCAGGATTCTGGGAGGCTTCAGTTACCGCTGGCCGAGCTGAAGAAGTGGGTATGAGGCTGGGGC  
35 GGGGCTGGAGGTGGCGCCCCCTGGTGGGACAACAAAGAGGACACCATTTTTCCAGAGCTGCAGAGAGCACCTGGT  
GGGAGGAAGAAGTGTAACCTACCAGCCTCTGCTCTTATCTTTGTAATAAATGTTAAAGCCAG 3'

SEQ ID NO: 7

18/29

FIG. 13





FIG.14




20/29

FIG.15


A

-396 gtccttcgggtgtctctgggtgtggtgagtagaggtgtgtgtcacaaagtagacaccattgtgtgtgacaaagcc  
 -324 catcgtgtgtctgtgtgtgtctttatccaagtgatggagcgtctctttctgtctgtgcccacagacacacc  
 -252 tagccccctcttattctcaaaaggggagctggggagcctccccctaccctggggcctccccctgccccccc  
 -180 cgccccctgctggcgtcaccactccccagagggcagagggctctgctgtgacctcagagcaaaagtcccagag  
 -108 ccagcagagcaggtgacgacctgcaagccagtggtgacctgtgctgctgaggtgggggacctgtg  
 -36 gcaggaagctggtgagcccccaagacccgggggcatggcggggatctgggtgctggcctggggggccttg  
 1 M G G D L V L G L G A L  
 37 AGACGCCGAAAGCGCTTGCTGGAGCAGGAGAAGTCTCTGGCCGGCTGGGCACCTGGTGGTGGCAGGAAGTGGC  
 13 R R R K R L L E Q E K S L A G W A L V L A G T G  
 109 ATTGGACTCATGGTGTGCTGCATGCAGAGATGCTGTGGTTCGGGGGGTGTCTGGGGGCTCTACCTGTTCCTG  
 37 L G L M V L H A E S M L W F G G C S W A L Y L F L  
 181 GTTAAATGCACGATCAGCATTTCCACCTTCTTACTCCTCTGCCTCATCGTGGCCTTTCATGCCAAAGAGGTC  
 61 V K C T I S I S T F L L L C L I V A F H A K E V  
 253 CAGCTGTTTCATGACGACACACGGGCTGCGGACTGCGCGTGGCGCTGACCGGGCGGCGAGCGGGCAGATC  
 85 Q L F M T D N G L R D W R V A L T G R Q A A Q I  
 325 GTGCTGGAGCTGGTGTGTGTGGCTGCACCGGGCGCCGTGCGGGGCCCGCGTGGGTGCGAGGATTTAGGG  
 109 N L E L V V C G L H P A P V R G P P C V Q D L G  
 397 GCGCGCTGACCTCCCCGAGCCTTGGCGGGATTCTGGGCCAAGGGGAAGCGCTGCTGTCCCTGGCCATG  
 133 A P L T S P Q P W P G F L G Q G E A L L S L A M  
 469 CTGCTGGCTCTCTACCTGGTGGCCCGCGCGTCTCTGCGCAGCGCGCTCTGCTCAACGCTTCTACCGC  
 157 L L R L Y L V P R A V L L R S G V L L N A S Y R  
 541 AGCATCGGCGCTCTCAATCAAGTCCGCTTCCGCCACTGGTTCGTGGCCAAAGCTTTACATGAACACGACCCCT  
 181 S I G A L N Q V R F R H W F V A K L Y M N T H P  
 613 GGCGCGCTGCTGCTCGGCCCTCACGCTTGGCCTCTGGCTGACCACCGCCTGGGTGCTGTCCGTGGCCGAGAGG  
 205 G R L L L G L T L G L W L T T A W V L S V A E R  
 685 CAGGCTGTTAATGCCACTGGGCACTTTAGACACACTTTGGCTGATCCCCATCACATTECTGACCATCGGC  
 229 Q A V A A T G H L S D T L W L I P I T F L T I G  
 757 TATGGTGACGTGGTGGCGGACCATGTGGGGCAAGATCGTCTGCTGTGCACTGGAGTCATGGGTGCTGTC  
 253 Y G D V V P G T M W G K I V C L C T G V M G V C  
 829 TGCACAGCCCTGCTGGTGGCCGTGGTGGCCCGGAAGCTGGAGTTTAAACAGGAGAGAACGACGTGCACAA  
 277 C T A L L V A V V A R K L E F N K A E K H V H N  
 901 TTCATGATGGATATCCAGTATACCAAGAGATGAAGGAGTCCGCTGCCGAGTGCTACAAGAAGCCTGGATG  
 301 F M M D I Q Y T K E M K E S A A R V L Q E A W M  
 973 TTCTACAACATACTCGCAGGAAGGAGTCTCATGTGCCCCGAGGCATCAGCGCAAGCTGTGGCCGCGCATC  
 325 F Y K H R R K E S H A A R R R H Q R K L L A A I  
 1045 AACCGCTTCCGCCAGGTGCGGCTGAAACACCGGAAGCTCCGGGAACAAGTGAATCCATGGTGGACATCTCC  
 349 N A F R Q V R L K H R K L R E Q V N S M V D I S  
 1117 AAGATGCACATGATCCTGTATGACCTGCAGCAGAATCTGAGCAGCTCACACCGGGCCCTGGAGAAACAGATT  
 373 K M H M I L Y D L Q Q N L S S S H R A L E K Q I  
 1189 GACACGCTGGCGGGGAAGCTGGATGCCCTGACTGAGCTGCTTAGCACTGCCCTGGGGCCGAGGCAGCTTCCA  
 397 D T L A G K L D A L T E L L S T A L G P R Q L P  
 1261 GAACCCAGCCAGCAGTCCAAGTAGctggaccacagaggaggaaccaggctactttccccagtactgaggtgg  
 421 E P S Q Q S K \*  
 1333 tggacatcgtctctgccaactcctgacccagccctgaacaaagcacctcaagtgaaggaccaaagggggccc  
 1405 tggcttggagtgggttggcttctgatggctgctggaggggacgctggctaaagtgggtaggccttggccca  
 1477 cctgagggccccaggtgggaacatggtcaccacccactctgcataccctcatcaaaaacactctcactatgctg  
 1549 ctatggagcagctccagctctcagttacaagtgcagggcagctggagggcaggaactcctgggtccctgggaaag  
 1621 aggttactagggggcccgatccaggattctgggaggtctcagttaccgctggccgagctgaagaactgggtga  
 1693 tgaggtctggggcggggctggaggtggcgccccctggtgggacaacaaaggagacaccattttccagagctg  
 1765 cagagagcacctggtggggaggaagaagtgaactcaccagcctctgctcttatcttgaataaagttaa  
 1837 agccag


B



C



21/29



5


10

FIG. 16

BEST AVAILABLE COPY

22/29

## hKCa4 Expression



5

Text in paranthesis represent channel numbers/cell.


10

15

FIG. 17




23/29



5

FIG. 18




5

Western Blot analysis of KCa4 antibodies

FIG. 19

BEST AVAILABLE COPY

Kca channel in anergic T cells (3926 cells)



Summary

Human T Cells	Resting	"Anergic"	Activated
Current (pA)	120 (n=7)	73 (n=11)	414 (n=11)
SE	62	33	82
P Value	<0.05	<0.0005	

5


p value for resting and anergic groups were compared to the activated group.

Down-regulation of KCa4 currents in T cell-lines induced to anergy

10

FIG.20

26/29



5

Deletion analysis of KCa4 C-tail to narrow the CAM interaction site

10

FIG.21

27/29

**mKCa4 cDNA sequence**

5

CAAGCTGGCTGCCATTTGCGTACTGAGAGGTGGGGGACCTGGGCAGGAAGCTGGCTGAGCCCC  
AAGACCTCAGGGGGCCATGGGCGGGGAGCTGGTGACTGGCCTGGGGGCCCTGAGACGGAGAAAG  
CGCCTGCTGGAGCAGGAGAAGAGGGTGGCCGGCTGGGCGTTGGTGCTGGCGGGAAC TGGCATC  
GGACTCATGGTTCTGCACGCTGAGATGTTGTGGTTCTTGGGCTGCAAGTGGGTGCTGTACCTG  
10 CTCCTGGTTAAGTGTTTGATCACCCCTGTCCACTGCCTTCCTCCTTTGTCTTATTGTGGTCTTC  
CATGCCAAGGAGGTCCAGCTGTTTCATGACTGACAACGGGCTCCGGGACTGGCGCGTGGCGCTG  
ACCCGGCGGCAGGTGGCGCAGATCCTGCTGGAGCTGTTGGTGTGCGGGGTGCACCCGGTGCCC  
CTACGGAGCCCGCACTGCGCCCTGGCGGGGAGGCCACCGACGCGCAGCCCTGGCCGGGTTTC  
CTGGGCGAAGGCGAGGCGTTGCTGTCCCTGGCCATGCTCCTGCGTCTCTACCTGGTGCCCCGC  
15 GCGGTGCTGCTGCGCAGCGGGGTCCTGCTCAACGCGTCCTACCGCAGCATCGGGGCGCTCAAC  
CAAGTCCGCTTCCGCCACTGGTTTCGTGGCCAAGCTGTACATGAACACGCACCCGGGTGCGCTG  
CTGCTGGGCCTCACGCTGGGTCTCTGGCTCACCACAGCTTGGGTGCTGTCTGTGGCTGAGAGG  
CAGGCTGTCAATGCCACGGGGCACCTCACAGACACACTGTGGCTGATTCCGATCACATTCCTG  
ACCATTTGGCTATGGGGACGTGGTACCTGGCACCATGTGGGGCAAGATTGTCTGCCTGTGCACC  
20 GGAGTCATGGGGGTCTGCTGCACAGCTCTCCTGGTGGCTGTGGTGGCTCGGAAGCTGGAGTTC  
AACAAGGCGGAGAAACACGTGCACAACTTCATGATGGACATCCATTATGCCAAAGAGATGAAG  
GAGTCAGCGGCGCGGCTGCTGCAGGAAGCCTGGATGTACTACAAGCACACTCGAAGGAAGGAC  
TCCCGGGCTGCCCCGAGACATCAGCGCAAGATGCTGGCCGCCATCCACACGTTCCGCCAGGTA  
CGGCTGAAACACCGGAAGCTCCGGGAACAAGTGAATTCCATGGTGGACATCTCCAAGATGCAC  
25 ATGATCCTGTGCGACCTGCAGCTGGGTCTCAGCTCCTCGCACCGTGCCCTGGAGAAGAGAATC  
GACGGTCTGGCAGGAAAGCTGGATGCCCTGACAGAGCTGCTCGGCACTGCTCTGCAGCAACAG  
CAGCTACCAGAACCCAGTCAGGAGGCCACATAGCTCCACATGAACTCACAGAAGAACA

30

**FIG.22**

28/29

**mKCa4 amino acid sequence**

5 MGGEELVTGLGALRRRKRLLEQEKRVAGWALVLAGTGIGLMVLHAEMLWFLGCKWVLYLLLVKC  
LITLSTAFLLCLIVVFHAKVQLFMTDNGLRDWRVALTRRQVAQILLELLVCGVHPVPLRSPH  
CALAGEATDAQPWPGFLGEGEALLSLAMLLRLYLVPRAVLLRSGVLLNASYRSIGALNQVRFR  
HWFVAKLYMNTHPGRLLLGLTLGLWLTTAWVLSVAERQAVNATGHLTDTLWLIPITFLTIGYG  
DVVPGTMWGKIVCLCTGVMGVCCTALLVAVVARKLEFNKAEKHVHNFMMDIHYAKEMKESAAR  
10 LLQEAWMYKHTRRKDSRAARRHQKMLAAIHTFRQVRLKHRKLREQVNSMVDISKMHMILCD  
LQLGLSSSHRALEKRIDGLAGKLDALTELLGTALQQQQLPEPSQEAT\*

15

**FIG.23**

20

25

30

29/29

Alignment Report of hvsmKCa4.MEG, using Clustal method with PAM250 residue weight table.  
 Wednesday, June 17, 1998 4:23 PM

Page 1

		MGGDLVLGLGALRRRKRLLLEQEKSVACWALVLAGTGIGLMVLHAEMLWFG																																								Majority
		10										20										30										40										50
1		MGGDLVLGLGALRRRKRLLLEQEKSLAGWALVLAGTGIGLMVLHAEMLWFG																																								Hkca4.pro
1		MGC[ELV]T[GLGALRRRKRLLLEQEK]RV[AGWALVLAGTGIGLMVLHAEMLW]FL																																								mKCa4.pro
		GCSWALYLLLVKCLISLSTALLLCLIVAFHAKREVQLFMTDNGLRDWRVAL																																								Majority
		60										70										80										90										100
51		GCSWALYLLLVKCTISISTFLLCLIVAFHAKREVQLFMTDNGLRDWRVAL																																								Hkca4.pro
51		GC[KW]V[L]Y[L]L[V]KCL[LI]T[ST]A[F]LLCLIV[V]FHAKEVQLFMTDNGLRDWRVAL																																								mKCa4.pro
		TGROAAQIVLELVVCGVHPAPVRGPHCVQALGGEATSAQFWPGFLGEGEA																																								Majority
		110										120										130										140										150
101		TGROAAQIVLELVVCGVHPAPVRGPHCVQALGGEATSAQFWPGFLGEGEA																																								Hkca4.pro
101		TRROV[AOI]L[LEL]LVCGV[HPV]P[PLR]S[PHC]--AL[AGEAT]DA[QFWPGFLGEGEA]																																								mKCa4.pro
		LLSLAMLLRLLYLVPRAVLLRSGVLLNASYRSIGALNOVRFRRHWFVAKLYM																																								Majority
		160										170										180										190										200
151		LLSLAMLLRLLYLVPRAVLLRSGVLLNASYRSIGALNOVRFRRHWFVAKLYM																																								Hkca4.pro
149		LLSLAMLLRLLYLVPRAVLLRSGVLLNASYRSIGALNOVRFRRHWFVAKLYM																																								mKCa4.pro
		NTHPGRLLLLGLTLGLWLTTANVLSVAERQAVNATGHLSDTLWLIPITFLT																																								Majority
		210										220										230										240										250
201		NTHPGRLLLLGLTLGLWLTTANVLSVAERQAVNATGHLSDTLWLIPITFLT																																								Hkca4.pro
199		NTHPGRLLLLGLTLGLWLTTANVLSVAERQAVNATGHL[DTLWLIPITFLT]																																								mKCa4.pro
		IGYGDVVPGTMMNSKIVCLCTGVHGVCTALLVAVVARKLEFNKAEKRVHN																																								Majority
		260										270										280										290										300
251		IGYGDVVPGTMMNSKIVCLCTGVHGVCTALLVAVVARKLEFNKAEKRVHN																																								Hkca4.pro
249		IGYGDVVPGTMMNSKIVCLCTGVHGVCTALLVAVVARKLEFNKAEKRVHN																																								mKCa4.pro
		FMNDIQYAKEMKESAAARVLOSAMHFKYKHTRRKDSHAARRHQRKLLAAINA																																								Majority
		310										320										330										340										350
301		FMNDIQYAKEMKESAAARVLOSAMHFKYKHTRRKDSHAARRHQRKLLAAINA																																								Hkca4.pro
299		FMNDI[HYAKEMKESAAAR]L[LOEAM]FYKHTRRK[DSR]AARRHQRK[HLAAI]HT																																								mKCa4.pro
		FROVRLKHKRLREOVNSMVDISKNNHILYDLOGLSSSHRALEKQIDGLA																																								Majority
		360										370										380										390										400
351		FROVRLKHKRLREOVNSMVDISKNNHILYDLOGLSSSHRALEKQIDGLA																																								Hkca4.pro
349		FRQVRLKHKRLREOVNSMVDISKNNHIL[C]D[LQ]LGLSSSHRALEK[R]I[DGLA]																																								mKCa4.pro
		GKLDALTELLGTALGQQQLPEPSQEAT-																																								Majority
		410										420																														
401		GKLDALTELLGTALGQQQLPEPSQ[QSK].																																								Hkca4.pro
399		GKLDALTELLGTALGQQQLPEPSQ[EAT].																																								mKCa4.pro

Sequence pair distances of hvsmKCa4.MEG, using Clustal method with PAM250 residue weight table.  
 Wednesday, June 17, 1998 4:23 PM

Page 1

Percent Similarity					
	1	2		Hkca4.pro	mKCa4.pro
1		86.6	1		
2	13.7		2		
	1	2			

## ALIGNMENT OF MOUSE AND HUMAN KCa4

FIG.24

**PCT**WORLD INTELLECTUAL PROPERTY ORGANIZATION  
International Bureau

## INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

<b>(51) International Patent Classification <sup>6</sup> :</b> C12N 15/12, C07K 14/705, 16/28, C12N 5/10, A61K 38/17, G01N 33/53, 33/68		<b>A3</b>	<b>(11) International Publication Number:</b> WO 99/03882 <b>(43) International Publication Date:</b> 28 January 1999 (28.01.99)
<b>(21) International Application Number:</b> PCT/GB98/02058 <b>(22) International Filing Date:</b> 13 July 1998 (13.07.98) <b>(30) Priority Data:</b> 9714760.7 15 July 1997 (15.07.97) GB 9721366.4 9 October 1997 (09.10.97) GB <b>(71) Applicant:</b> ZENECA LIMITED [GB/GB]; 15 Stanhope Gate, London W1Y 6LN (GB). <b>(72) Inventors:</b> AIYAR, Jayashree; 1800 Concord Pike, P.O. Box 15437, Wilmington, DE 19850-8437 (US). LOGSDON, Naomi, Jean; 1800 Concord Pike, P.O. Box 15437, Wilm- ington, DE 19850-8437 (US). <b>(74) Agent:</b> PHILLIPS, Neil, Godfrey, Alasdair; Zeneca Pharmaceu- ticals, Intellectual Property Dept., Mereside, Alderley Park, Macclesfield, Cheshire SK10 4TG (GB).		<b>(81) Designated States:</b> AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB, GE, GH, GM, HU, ID, IL, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, UZ, VN, YU, ZW, ARIPO patent (GH, GM, KE, LS, MW, SD, SZ, UG, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG).  <b>Published</b> <i>With international search report.</i> <b>(88) Date of publication of the international search report:</b> 8 April 1999 (08.04.99)	
<b>(54) Title:</b> HUMAN LEUKOCYTE CALCIUM ACTIVATED POTASSIUM CHANNEL POLYPEPTIDE			
<b>(57) Abstract</b>  A novel human leukocyte calcium activated potassium channel polypeptide is described which is expressed at a high level in activated T-cells. A full length cDNA which encodes the novel calcium activated potassium channel polypeptide is disclosed as well as the interior structural region and the amino acid residue sequence of the native biological molecule. Methods are provided to identify compounds that modulate the biological activity of a human leukocyte calcium-activated potassium channel.			



**FOR THE PURPOSES OF INFORMATION ONLY**

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AL	Albania	ES	Spain	LS	Lesotho	SI	Slovenia
AM	Armenia	FI	Finland	LT	Lithuania	SK	Slovakia
AT	Austria	FR	France	LU	Luxembourg	SN	Senegal
AU	Australia	GA	Gabon	LV	Latvia	SZ	Swaziland
AZ	Azerbaijan	GB	United Kingdom	MC	Monaco	TD	Chad
BA	Bosnia and Herzegovina	GE	Georgia	MD	Republic of Moldova	TG	Togo
BB	Barbados	GH	Ghana	MG	Madagascar	TJ	Tajikistan
BE	Belgium	GN	Guinea	MK	The former Yugoslav Republic of Macedonia	TM	Turkmenistan
BF	Burkina Faso	GR	Greece	ML	Mali	TR	Turkey
BG	Bulgaria	HU	Hungary	MN	Mongolia	TT	Trinidad and Tobago
BJ	Benin	IE	Ireland	MR	Mauritania	UA	Ukraine
BR	Brazil	IL	Israel	MW	Malawi	UG	Uganda
BY	Belarus	IS	Iceland	MX	Mexico	US	United States of America
CA	Canada	IT	Italy	NE	Niger	UZ	Uzbekistan
CF	Central African Republic	JP	Japan	NL	Netherlands	VN	Viet Nam
CG	Congo	KE	Kenya	NO	Norway	YU	Yugoslavia
CH	Switzerland	KG	Kyrgyzstan	NZ	New Zealand	ZW	Zimbabwe
CI	Côte d'Ivoire	KP	Democratic People's Republic of Korea	PL	Poland		
CM	Cameroon	KR	Republic of Korea	PT	Portugal		
CN	China	KZ	Kazakhstan	RO	Romania		
CU	Cuba	LC	Saint Lucia	RU	Russian Federation		
CZ	Czech Republic	LI	Liechtenstein	SD	Sudan		
DE	Germany	LK	Sri Lanka	SE	Sweden		
DK	Denmark	LR	Liberia	SG	Singapore		
EE	Estonia						

# INTERNATIONAL SEARCH REPORT

International Application No

PCT/GB 98/02058

## A. CLASSIFICATION OF SUBJECT MATTER

IPC 6 C12N15/12 C07K14/705 C07K16/28 C12N5/10 A61K38/17  
G01N33/53 G01N33/68

According to International Patent Classification (IPC) or to both national classification and IPC

## B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

IPC 6 C07K C12N

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practical, search terms used)

## C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category *	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
A	KOHLER M ET AL: "Small-conductance, calcium-activated potassium channels from mammalian brain 'see comments!'" SCIENCE, SEP 20 1996, 273 (5282) P1709-14, XP002091031 UNITED STATES cited in the application see abstract; figure 2 ---	1-13, 16, 17, 23-26
P, X	WO 98 11139 A (UNIV OREGON HEALTH SCIENCES ;ADELMAN JOHN P (US); BOND CHRIS T (US) 19 March 1998 see claim 2; figure SEQ.ID.32 --- -/--	1-13, 16, 17

☒ Further documents are listed in the continuation of box C.

☒ Patent family members are listed in annex.

### \* Special categories of cited documents:

- "A" document defining the general state of the art which is not considered to be of particular relevance
- "E" earlier document but published on or after the international filing date
- "L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)
- "O" document referring to an oral disclosure, use, exhibition or other means
- "P" document published prior to the international filing date but later than the priority date claimed

- "T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention
- "X" document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone
- "Y" document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art.
- "&" document member of the same patent family

Date of the actual completion of the international search

26 January 1999

Date of mailing of the international search report

04/02/1999

Name and mailing address of the ISA

European Patent Office, P.B. 5818 Patentlaan 2  
NL - 2280 HV Rijswijk  
Tel. (+31-70) 340-2040, Tx. 31 651 epo nl,  
Fax: (+31-70) 340-3016

Authorized officer

Gurdjian, D

# INTERNATIONAL SEARCH REPORT

International Application No

PCT/GB 98/02058

## C.(Continuation) DOCUMENTS CONSIDERED TO BE RELEVANT

Category *	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
P,X	<p>JOINER WJ ET AL: "hSK4, a member of a novel subfamily of calcium-activated potassium channels."</p> <p>PROC NATL ACAD SCI U S A, SEP 30 1997, 94 (20) P11013-8, XP002091032</p> <p>UNITED STATES</p> <p>see the whole document</p> <p>---</p>	1-13
P,X	<p>LOGSDON NJ ET AL: "A novel gene, hKCa4, encodes the calcium-activated potassium channel in human T lymphocytes."</p> <p>J BIOL CHEM, DEC 26 1997, 272 (52) P32723-6, XP002091033</p> <p>UNITED STATES</p> <p>see the whole document</p> <p>-----</p>	1-13

# INTERNATIONAL SEARCH REPORT

International application No.

PCT/GB 98/ 02058

## Box I Observations where certain claims were found unsearchable (Continuation of item 1 of first sheet)

This International Search Report has not been established in respect of certain claims under Article 17(2)(a) for the following reasons:

1. ☒ Claims Nos.: 14, 15, 18-22 and 27, 28  
because they relate to subject matter not required to be searched by this Authority, namely:  
see FURTHER INFORMATION sheet PCT/ISA/210
2. ☐ Claims Nos.:  
because they relate to parts of the International Application that do not comply with the prescribed requirements to such an extent that no meaningful International Search can be carried out, specifically:
3. ☐ Claims Nos.:  
because they are dependent claims and are not drafted in accordance with the second and third sentences of Rule 6.4(a).

## Box II Observations where unity of invention is lacking (Continuation of item 2 of first sheet)

This International Searching Authority found multiple inventions in this international application, as follows:

1. ☐ As all required additional search fees were timely paid by the applicant, this International Search Report covers all searchable claims.
2. ☐ As all searchable claims could be searched without effort justifying an additional fee, this Authority did not invite payment of any additional fee.
3. ☐ As only some of the required additional search fees were timely paid by the applicant, this International Search Report covers only those claims for which fees were paid, specifically claims Nos.:
4. ☐ No required additional search fees were timely paid by the applicant. Consequently, this International Search Report is restricted to the invention first mentioned in the claims; it is covered by claims Nos.:

Remark on Protest

- ☐ The additional search fees were accompanied by the applicant's protest.
- ☐ No protest accompanied the payment of additional search fees.

FURTHER INFORMATION CONTINUED FROM PCT/ISA/ 210

Although claims 14,18-20 ,22 and claim 15 , as far as it concerns an in vivo method , are directed to a method of treatment of the human/animal body, the search has been carried out and based on the alleged effects of the compound/composition.

-----  
Claims Nos.: 14,15,18-22 and 27,28

Claims 21,22 and 27,28 relating to modulators of the polypeptide of claim 6,23 could not be searched as its subject-matter was insufficiently disclosed .

# INTERNATIONAL SEARCH REPORT

Information on patent family members

International Application No

PCT/GB 98/02058

Patent document cited in search report	Publication date	Patent family member(s)	Publication date
WO 9811139 A	19-03-1998	AU 4266097 A	02-04-1998

**PCT**WORLD INTELLECTUAL PROPERTY ORGANIZATION  
International Bureau

## INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

<b>(51) International Patent Classification <sup>6</sup> :</b> C12N 15/12, C07K 14/705, 16/28, C12N 5/10, A61K 38/17, G01N 33/53, 33/68	<b>A3</b>	<b>(11) International Publication Number:</b> WO 99/03882 <b>(43) International Publication Date:</b> 28 January 1999 (28.01.99)
<b>(21) International Application Number:</b> PCT/GB98/02058 <b>(22) International Filing Date:</b> 13 July 1998 (13.07.98)  <b>(30) Priority Data:</b> 9714760.7      15 July 1997 (15.07.97)      GB 9721366.4      9 October 1997 (09.10.97)      GB  <b>(71) Applicant:</b> ZENECA LIMITED [GB/GB]; 15 Stanhope Gate, London W1Y 6LN (GB).  <b>(72) Inventors:</b> AIYAR, Jayashree; 1800 Concord Pike, P.O. Box 15437, Wilmington, DE 19850-8437 (US). LOGSDON, Naomi, Jean; 1800 Concord Pike, P.O. Box 15437, Wilmington, DE 19850-8437 (US).  <b>(74) Agent:</b> PHILLIPS, Neil, Godfrey, Alasdair; Zeneca Pharmaceuticals, Intellectual Property Dept., Mereside, Alderley Park, Macclesfield, Cheshire SK10 4TG (GB).		<b>(81) Designated States:</b> AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB, GE, GH, GM, HU, ID, IL, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, UZ, VN, YU, ZW, ARIPO patent (GH, GM, KE, LS, MW, SD, SZ, UG, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG).  <b>Published</b> <i>With international search report.</i>  <b>(88) Date of publication of the international search report:</b> 8 April 1999 (08.04.99)
<b>(54) Title:</b> HUMAN LEUKOCYTE CALCIUM ACTIVATED POTASSIUM CHANNEL POLYPEPTIDE  <b>(57) Abstract</b>  A novel human leukocyte calcium activated potassium channel polypeptide is described which is expressed at a high level in activated T-cells. A full length cDNA which encodes the novel calcium activated potassium channel polypeptide is disclosed as well as the interior structural region and the amino acid residue sequence of the native biological molecule. Methods are provided to identify compounds that modulate the biological activity of a human leukocyte calcium-activated potassium channel.		

**FOR THE PURPOSES OF INFORMATION ONLY**

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AL	Albania	ES	Spain	LS	Lesotho	SI	Slovenia
AM	Armenia	FI	Finland	LT	Lithuania	SK	Slovakia
AT	Austria	FR	France	LU	Luxembourg	SN	Senegal
AU	Australia	GA	Gabon	LV	Latvia	SZ	Swaziland
AZ	Azerbaijan	GB	United Kingdom	MC	Monaco	TD	Chad
BA	Bosnia and Herzegovina	GE	Georgia	MD	Republic of Moldova	TG	Togo
BB	Barbados	GH	Ghana	MG	Madagascar	TJ	Tajikistan
BE	Belgium	GN	Guinea	MK	The former Yugoslav Republic of Macedonia	TM	Turkmenistan
BF	Burkina Faso	GR	Greece			TR	Turkey
BG	Bulgaria	HU	Hungary	ML	Mali	TT	Trinidad and Tobago
BJ	Benin	IE	Ireland	MN	Mongolia	UA	Ukraine
BR	Brazil	IL	Israel	MR	Mauritania	UG	Uganda
BY	Belarus	IS	Iceland	MW	Malawi	US	United States of America
CA	Canada	IT	Italy	MX	Mexico	UZ	Uzbekistan
CF	Central African Republic	JP	Japan	NE	Niger	VN	Viet Nam
CG	Congo	KE	Kenya	NL	Netherlands	YU	Yugoslavia
CH	Switzerland	KG	Kyrgyzstan	NO	Norway	ZW	Zimbabwe
CI	Côte d'Ivoire	KP	Democratic People's Republic of Korea	NZ	New Zealand		
CM	Cameroon	KR	Republic of Korea	PL	Poland		
CN	China	KZ	Kazakhstan	PT	Portugal		
CU	Cuba	LC	Saint Lucia	RO	Romania		
CZ	Czech Republic	LI	Liechtenstein	RU	Russian Federation		
DE	Germany	LK	Sri Lanka	SD	Sudan		
DK	Denmark	LR	Liberia	SE	Sweden		
EE	Estonia			SG	Singapore		



# INTERNATIONAL SEARCH REPORT

International Application No

PCT/GB 98/02058

## A. CLASSIFICATION OF SUBJECT MATTER

IPC 6 C12N15/12 C07K14/705 C07K16/28 C12N5/10 A61K38/17  
G01N33/53 G01N33/68

According to International Patent Classification (IPC) or to both national classification and IPC

## B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

IPC 6 C07K C12N

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practical, search terms used)

## C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category *	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
A	KOHLER M ET AL: "Small-conductance, calcium-activated potassium channels from mammalian brain 'see comments!'" SCIENCE, SEP 20 1996, 273 (5282) P1709-14, XP002091031 UNITED STATES cited in the application see abstract; figure 2 ---	1-13, 16, 17, 23-26
P, X	WO 98 11139 A (UNIV OREGON HEALTH SCIENCES ;ADELMAN JOHN P (US); BOND CHRIS T (US) 19 March 1998 see claim 2; figure SEQ.ID.32 ---	1-13, 16, 17
	-/--	

☒ Further documents are listed in the continuation of box C.

☒ Patent family members are listed in annex.

\* Special categories of cited documents:

"A" document defining the general state of the art which is not considered to be of particular relevance

"E" earlier document but published on or after the international filing date

"L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)

"O" document referring to an oral disclosure, use, exhibition or other means

"P" document published prior to the international filing date but later than the priority date claimed

"T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention

"X" document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone

"Y" document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art

"&" document member of the same patent family

Date of the actual completion of the international search

26 January 1999

Date of mailing of the international search report

04/02/1999

Name and mailing address of the ISA

European Patent Office, P.B. 5818 Patentlaan 2  
NL - 2280 HV Rijswijk  
Tel. (+31-70) 340-2040, Tx. 31 651 epo nl,  
Fax: (+31-70) 340-3016

Authorized officer

Gurdjian, D

# INTERNATIONAL SEARCH REPORT

International Application No

PCT/GB 98/02058

## C.(Continuation) DOCUMENTS CONSIDERED TO BE RELEVANT

Category *	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
P,X	<p>JOINER WJ ET AL: "hSK4, a member of a novel subfamily of calcium-activated potassium channels."</p> <p>PROC NATL ACAD SCI U S A, SEP 30 1997, 94 (20) P11013-8, XP002091032</p> <p>UNITED STATES</p> <p>see the whole document</p> <p>---</p>	1-13
P,X	<p>LOGSDON NJ ET AL: "A novel gene, hKCa4, encodes the calcium-activated potassium channel in human T lymphocytes."</p> <p>J BIOL CHEM, DEC 26 1997, 272 (52) P32723-6, XP002091033</p> <p>UNITED STATES</p> <p>see the whole document</p> <p>-----</p>	1-13

# INTERNATIONAL SEARCH REPORT

International application No.

PCT/GB 98/ 02058

## Box I Observations where certain claims were found unsearchable (Continuation of item 1 of first sheet)

This International Search Report has not been established in respect of certain claims under Article 17(2)(a) for the following reasons:

1. ☒ Claims Nos.: 14, 15, 18-22 and 27, 28  
because they relate to subject matter not required to be searched by this Authority, namely:  
see FURTHER INFORMATION sheet PCT/ISA/210
2. ☐ Claims Nos.:  
because they relate to parts of the International Application that do not comply with the prescribed requirements to such an extent that no meaningful International Search can be carried out, specifically:
3. ☐ Claims Nos.:  
because they are dependent claims and are not drafted in accordance with the second and third sentences of Rule 6.4(a).

## Box II Observations where unity of invention is lacking (Continuation of item 2 of first sheet)

This International Searching Authority found multiple inventions in this international application, as follows:

1. ☐ As all required additional search fees were timely paid by the applicant, this International Search Report covers all searchable claims.
2. ☐ As all searchable claims could be searched without effort justifying an additional fee, this Authority did not invite payment of any additional fee.
3. ☐ As only some of the required additional search fees were timely paid by the applicant, this International Search Report covers only those claims for which fees were paid, specifically claims Nos.:
4. ☐ No required additional search fees were timely paid by the applicant. Consequently, this International Search Report is restricted to the invention first mentioned in the claims; it is covered by claims Nos.:

Remark on Protest

☐ The additional search fees were accompanied by the applicant's protest.

☐ No protest accompanied the payment of additional search fees.

FURTHER INFORMATION CONTINUED FROM PCT/ISA/ 210

Although claims 14,18-20 ,22 and claim 15 , as far as it concerns an in vivo method , are directed to a method of treatment of the human/animal body, the search has been carried out and based on the alleged effects of the compound/composition.

-----  
Claims Nos.: 14,15,18-22 and 27,28

Claims 21,22 and 27,28 relating to modulators of the polypeptide of claim 6,23 could not be searched as its subject-matter was insufficiently disclosed .

# INTERNATIONAL SEARCH REPORT

Information on patent family members

International Application No

PCT/GB 98/02058

Patent document cited in search report	Publication date	Patent family member(s)	Publication date
WO 9811139 A	19-03-1998	AU 4266097 A	02-04-1998